

Apr 08, 2024

Immunofluorescence staining protocol with Antigen Retrieval

DOI

dx.doi.org/10.17504/protocols.io.bp2l62nk1gqe/v1

madalynn.erb Erb¹

¹Van Andel Research Institute

Team Lee



Jane Balster

ASAP - Team Lee

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l62nk1gqe/v1>

Protocol Citation: madalynn.erb Erb 2024. Immunofluorescence staining protocol with Antigen Retrieval. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bp2l62nk1gqe/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: February 09, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 97924

Keywords: ASAPCRN, immunofluorescence staining, protocol with antigen retrieval, antigen retrieval, protocol details the immunofluorescence, immunofluorescence, antigen, staining protocol

Abstract

This protocol details the immunofluorescence staining protocol with antigen retrieval.

Materials


⊗ ProLong[®]; Diamond Antifade Mountant **Thermo Fisher Catalog #P36961**

⊗ Universal HIER antigen retrieval reagent (10X) **Abcam Catalog #ab208572**

Troubleshooting




Day 1

- 1 Staining protocol for $\pm 35\ \mu\text{m}$ free floating mouse brain sections.
- 2 Wash tissue sections 3 times (5 minutes each wash) in PBS to remove cryoprotectant solution.
- 2.1 Wash tissue sections for ⌚ 00:05:00 in PBS to remove cryoprotectant solution (1/3). 5m
- 2.2 Wash tissue sections for ⌚ 00:05:00 in PBS to remove cryoprotectant solution (2/3). 5m
- 2.3 Wash tissue sections for ⌚ 00:05:00 in PBS to remove cryoprotectant solution (3/3). 5m
- 3 Mount tissue onto positively charged slides.
- 4 Allow slides to air dry then place in oven set at $37\ ^\circ\text{C}$ - $50\ ^\circ\text{C}$ ⌚ Overnight . 5m

Note

Can use slide warmer or incubator for this step.

Day 2

- 5 If you would like to use a PAP pen to draw a box around your tissue sections, do that now and wait for slides to dry.
- 6 Wash slides in PBS for ⌚ 00:05:00 . 5m
- 7 Preheat steamer with a Coplin jar containing the antigen retrieval buffer until temperature reaches $95\ ^\circ\text{C}$ - $100\ ^\circ\text{C}$ (about ⌚ 00:10:00 10m


**Note**




This protocol uses Universal HIER antigen reagent as the buffer. This comes at 10X reagent, therefore, you need to dilute 1:10 with distilled water for 1x concentration of the desired volume.

- 8 Immerse slides into the preheated buffer in the steamer for 00:30:00 . 30m
- 9 After 30 minutes turn off the steamer but allow to cool to Room temperature (about 20- 00:30:00) before progressing to the next step. 30m
- 10 Wash slides in PBS 0.1% Triton X-100 (PBST) 2 × 2 min.
- 10.1 Wash slides in PBS 0.1% Triton X-100 (PBST) for 00:02:00 (1/2). 2m
- 10.2 Wash slides in PBS 0.1% Triton X-100 (PBST) for 00:02:00 (2/2). 2m
- 11 Wash slides in PBS 3 × 5 min.
- 11.1 Wash slides in PBS for 00:05:00 (1/3). 5m
- 11.2 Wash slides in PBS for 00:05:00 (2/3). 5m
- 11.3 Wash slides in PBS for 00:05:00 (3/3). 5m
- 12 Block for 02:00:00 at Room temperature in PBS (0.4% BSA, 10% Goat normal serum, 0.3% triton X-100). 2h



Note

- Move slides to black slide box for this step and subsequent steps.
- Use 1-2ml solution for washes and  700 μL for antibody incubations.

- 13 Incubate  Overnight with primary antibody in PBS (2% Goat normal serum, 0.3% triton X-100)  Overnight at  Room temperature .

20m



Day 3

9h

- 14 Wash slides in PBST 3 \times 10 min.



- 14.1 Wash slides in PBST for  00:10:00 (1/3).



10m

- 14.2 Wash slides in PBST for  00:10:00 (2/3).

10m

- 14.3 Wash slides in PBST for  00:10:00 (3/3).

10m

- 15 Incubate slides with secondary antibody in PBS for  04:00:00 at  Room temperature or at  4 $^{\circ}\text{C}$  Overnight .

8h



- 16 Wash slides in PBST 3 \times 10 min.



- 16.1 Wash slides in PBST for  00:10:00 (1/3).

10m

- 16.2 Wash slides in PBST for  00:10:00 (2/3).

10m

- 16.3 Wash slides in PBST for  00:10:00 (3/3).

10m



- 17 Incubate slides with Hoescht stain (1:5000) for ⌚ 00:30:00 at 🌡 Room temperature 30m
.
- 18 Wash slides in PBST 3 × 10 min. 🧴
- 18.1 Wash slides in PBST for ⌚ 00:10:00 (1/3). 10m
- 18.2 Wash slides in PBST for ⌚ 00:10:00 (2/3). 10m
- 18.3 Wash slides in PBST for ⌚ 00:10:00 (3/3). 10m
- 19 Mount slides using ProLong™ Diamond Antifade Mountant.
- 20 Dry slides on a flat surface at 🌡 Room temperature in the dark.
- Note

Put them under aluminum foil or a box on the bench.
- 21 Store slides at 🌡 4 °C the next day.