

Feb 23, 2026

Immunofluorescence under confinement

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvjn48ngk5/v1>

Ankita Jha¹

¹National Institutes of Health



Ankita Jha

National Institutes of Health

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvjn48ngk5/v1>

Protocol Citation: Ankita Jha 2026. Immunofluorescence under confinement. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.n2bvjn48ngk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 26, 2024



Last Modified: February 23, 2026

Protocol Integer ID: 108415

Keywords: immunofluorescence under confinement cell, agarose confinement, immunofluorescence, confinement cell, desired antibody

Abstract

Cells under agarose confinement were fixed and labelled with desired antibodies.



Immunofluorescence under confinement

5h 20m

- 1 Samples under agarose pad confinement were fixed with 4% paraformaldehyde (PFA; Electron Microscopy Sciences, 15710) in cytoskeleton buffer (CB; 10 mM MES, 3 mM MgCl_2 , 138 mM KCl, 2 mM EGTA) for 20 min at 37 °C, after which the agarose pad was gently removed. 20m
- 2 Cells were permeabilized with 0.5% Triton X-100 in CB for 5 min at room temperature (RT) and quenched with 10 mM glycine in CB for 5 min. 10m
- 3 Samples were washed with Tris-buffered saline (TBS; 20 mM Tris, pH 7.6, 137 mM NaCl) (2 × 5 min, then 2 × 10 min) and blocked for 1 h at RT in blocking buffer (2% IgG- and protease-free BSA and 0.1% Tween-20 in TBS). 1h 30m
- 4 Cells were incubated for 2h at RT or overnight at 4 °C with primary antibodies against CD44 (BJ18; BioLegend), Ezrin (3C12; Invitrogen), EGFR (H11; Invitrogen), or phospho-ERM (Thr567/564/558; Cell Signaling Technology) diluted in TBS. 2h
- 5 Following three TBS washes, samples were incubated for 1 h at RT with fluorophore-conjugated secondary antibodies (1:500; Jackson ImmunoResearch) and Alexa Fluor 647-phalloidin (1:200; Thermo Fisher Scientific) diluted in blocking buffer. 1h
- 6 Samples were washed with TBS (2 × 10 min) and mounted in Dako mounting medium (S3023) on glass-bottom dishes. 20m