

Apr 07, 2023

Immunofluorescent staining for neuronal marker MAP2

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlkwj46l5r/v1>

Qing Wang^{1,2}

¹Massachusetts General Hospital; ²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network



Jackson Schumacher

Massachusetts General Hospital

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlkwj46l5r/v1>

Protocol Citation: Qing Wang 2023. Immunofluorescent staining for neuronal marker MAP2. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.j8nlkwj46l5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 07, 2023

Last Modified: December 12, 2024

Protocol Integer ID: 80149

Keywords: ASAPCRN, neuronal marker map2, immunofluorescent staining

Funders Acknowledgements:

From cancer associations to altered immunity in the pathogenesis of Parkinson's disease

Grant ID: ASAP-000312

Abstract

This is the protocol for immunofluorescent staining for neuronal marker MAP2.

- 1 1. Treat differentiated SH-SY5Y cells with 40 ug/mL eumelanin or pheomelanin or PBS for 24 hours.
- 2 2. Process cells using Cytofix/Cytoperm™ fixation/permeabilization solution (BD554714, Thermo Fisher Scientific) and block with 5% normal goat serum.
- 3 3. Add primary antibody against microtubule-associated protein 2 (MAP2)(30 µg/mL, OSM00036G, Thermo Fisher Scientific) and incubate at 4°C overnight.
- 4 4. Add secondary antibody (1:1000, Alexa 594-conjugated, A11012, Thermo Fisher Scientific) and incubate for 2 hours at room temperature.
- 5 5. Stain nuclei with DAPI.
- 6 6. For MAP2-positive cell counting, three images are captured randomly from each well using FluoView FV300 confocal microscope under a 60x objective lens. MAP2 and DAPI channels are merged, and MAP2-positive cells in each random visual field of 0.045000 mm² are counted using ImageJ.

Protocol references

Cai, W., Srivastava, P., Feng, D., Lin, Y., Vanderburg, C. R., Xu, Y., ... & Chen, X. (2022). Melanocortin 1 receptor activation protects against alpha-synuclein pathologies in models of Parkinson's disease. *Molecular neurodegeneration*, 17(1), 1-21.

Royo, N. C., Vandenberghe, L. H., Ma, J. Y., Hauspurg, A., Yu, L., Maronski, M., ... & Watson, D. J. (2008). Specific AAV serotypes stably transduce primary hippocampal and cortical cultures with high efficiency and low toxicity. *Brain research*, 1190, 15-22.

Ye, Q., Srivastava, P., Al-Kuwari, N., & Chen, X. (2022). Oncogenic BRAFV600E induces microglial proliferation through extracellular signal-regulated kinase and neuronal death through c-Jun N-terminal kinase. *Neural Regeneration Research*.