

Aug 29, 2024

Immunofluorescent staining of dopaminergic neurons in brains of *D. melanogaster*

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2lyweyevx9/v1>

Natalie Kaempf^{1,2,3}, Uli Pech^{1,2}, Patrik Verstreken^{1,2}

¹VIB-KU Leuven Center for Brain & Disease Research, 3000 Leuven, Belgium;

²KU Leuven, Department of Neurosciences, Leuven Brain Institute, 3000 Leuven, Belgium;

³Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815



Natalie Kaempf

VIB KU Leuven

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2lyweyevx9/v1>

Protocol Citation: Natalie Kaempf, Uli Pech, Patrik Verstreken 2024. Immunofluorescent staining of dopaminergic neurons in brains of *D. melanogaster*. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.eq2lyweyevx9/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 17, 2024

Last Modified: August 29, 2024

Protocol Integer ID: 103595

Keywords: ASAPCRN, Immunohistochemistry, dopaminergic neurons, Drosophila brains, dopaminergic neurons in drosophila brain, dopaminergic neuron, dopaminergic neurons in brain, drosophila brain, brain

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000430

EMBO long-term postdoctoral fellowship

Grant ID: ALTF_299-2019

Research project, FWO Vlaanderen

Grant ID: G0A5219N

Research project, FWO Vlaanderen

Grant ID: G0B8119N

Methusalem project

Grant ID: METH/21/05 (3M210778)

Research project, KU Leuven Parkinson Fonds

Grant ID: EQZ-PARFON-O2010

Opening the Future grant, Leuven Universiteitsfonds (LUF)

Grant ID: EQZ-OPTFUP-O2010

Research project, FWO Vlaanderen

Grant ID: G031324N

Abstract

This protocol is used to visualize for dopaminergic neurons in Drosophila brains, but can be adapted for other targets as well.

Materials

Rabbit polyclonal anti-TH, Millipore, Cat# AB152, RRID: AB_390204

Mouse monoclonal anti-DLG, DSHB, DSHB Cat# 4F3 anti-discs large; RRID: AB_528203

Goat anti-Rabbit IgG Alexa FluorTM 488, Life Technologies, Cat# A-11034 RRID: AB_2576217

Goat anti-Mouse IgG Alexa FluorTM 555, Life Technologies, Cat# A-21424 RRID: AB_141780

RapiClear 1.47, Sunjin Lab, Cat# RC147001

Triton X-100 Solution, Sigma Aldrich, Cat# 93443-100ML

Formaldehyde solution (37 wt. % in H₂O) F1635-500ml, Sigma-Aldrich



Dissection of fly brains, fixation, block and primary antibody incubation

1h 30m

- 1 Prepare fresh 3.7% paraformaldehyde in 1x PBS, 0.2% Triton X-100 (PBX),  500 μ L per genotype and store  On ice .
- 2 Dissect fly brains in ice-cold PBS under stereomicroscope, collect with a PBX coated pipette and transfer in the 3.7% paraformaldehyde solution.
- 2.1 Store the dissected fly brains in paraformaldehyde solution on ice, but do not keep fly brains longer than  01:00:00  On ice 1h
- 3 Incubate the dissected fly brains in paraformaldehyde solution on rotation wheel for  00:30:00 at  Room temperature 30m
- 4 Take the paraformaldehyde solution off and wash a first time with PBX 0.2% ( 500 μ L) and take it off right away
- 5 Wash 3 times with PBX 0.2% ( 500 μ L) for 15 min on a rotator at  Room temperature
- 6 Block with 10% NGS in PBX ( 500 μ L) for  01:00:00 on rotator at  Room temperature 1h
- 7 Add primary antibody (in  500 μ L of 10% NGS in PBX) put it on rotation wheel at  4 $^{\circ}$ C  48:00:00 2d
- 7.1 Use rabbit a-TH, Sigma, 1:200, and mouse a-DLG, DSHB, 1:100

secondary antibody incubation

2d



- 8 Take the primary antibody solution off and wash a first time with PBX 0.2% ( 500 μ L) and take it off right away
- 9 Wash 3 times with PBX 0.2% ( 500 μ L) for 15 min on a rotator at  Room temperature
- 10 Add secondary antibody (in  500 μ L of 10% NGS in PBX), cover with tinfoil and put it on rotation wheel at  4 $^{\circ}$ C  Overnight 2d
- 10.1 Use Goat anti-Rabbit IgG Alexa Fluor™ 488, Life Technologies, and Goat anti-Mouse IgG Alexa Fluor™ 555, Life Technologies, at 1:500 dilution.

Mounting

15m

- 11 Take the secondary antibody solution off and wash a first time with PBX 0.2% ( 500 μ L) and take it off right away
- 12 Wash 3 times with PBX 0.2% ( 500 μ L) for 15 min on a rotator at  Room temperature
- 13 prepare slides:
 - 13.1 clean slides with EtOH
 - 13.2 place 2 book binder rings on top of each other
 - 13.3 fill book binder rings with 0.075% PLL solution, leave minimum  00:15:00 , take the solution off and leave a thin layer to dry 15m
- 14 Mount fly brains by pipetting the brains next to the book binders, fill the book binder chamber with PBX and transfer the brains individually with forceps to the book binder chamber, orient the brains correctly with anterior facing up (antenna lobe up) and softly push down until they stick
- 15 remove PBX



- 16 add 8 μ l of mounting medium (RapiClear 1.47) by gently pipetting on top of brains
- 17 add glass coverslip, seal with nail polish and let cure overnight at  4 °C in dark