

Jul 28, 2023

immunofluorescent staining with anti-GFP and anti-CD63 antibodies

DOI

dx.doi.org/10.17504/protocols.io.14egn2pypg5d/v1

rosanne.wouters¹, Peter Vangheluwe¹

¹KU Leuven



rosanne.wouters

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn2pypg5d/v1>

Protocol Citation: rosanne.wouters, Peter Vangheluwe 2023. immunofluorescent staining with anti-GFP and anti-CD63 antibodies. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.14egn2pypg5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 01, 2023



Last Modified: May 31, 2024

Protocol Integer ID: 82751

Keywords: ASAPCRN, fixed hela cell, hela cell, confocal imaging, cell, gfp

Funders Acknowledgements:

ASAP


Abstract

This protocol was used for immunofluorescent staining in fixed HeLa cells with anti-GFP and anti-CD63 antibodies, followed by confocal imaging.

Protocol materials

 FluorSave™ Reagent Merck MilliporeSigma (Sigma-Aldrich) Catalog #345789

Troubleshooting

- 1 plate cells on glass coverslips and grown until they reach 60% confluency
- 2 fix cells with 4% paraformaldehyde for 20 min at room temperature
- 3 permeabilize cells with 0.1% Triton X-100 in PBS for 5 min
- 4 block for 1 h with blocking buffer (PBS with 0.5% Tween20, 0.1% BSA, 0.2% FBS)
- 5 incubate coverslips with primary antibodies for 2 h at room temperature
(anti-CD63, exbio, 11-343-C100, mouse; anti-GFP, abcam, ab13970, chicken)
- 6 wash coverslips 3 times with PBS-T
- 7 incubated coverslips 30 min with secondary antibodies
(goat-anti-mouse-AlexaFluor647, goat-anti-chicken-AlexaFluor488)
- 8 wash coverslips 3 times with PBS-T
- 9 incubate coverslips with DAPI
- 10 wash coverslips 3 times with PBS-T
- 11 mount coverslips using
 FluorSave™ Reagent **Merck MilliporeSigma (Sigma-Aldrich) Catalog #345789**
- 12 images were acquired using an LSM780 confocal microscope (Zeiss) with a 10x or 40x objective



13 colocalization analysis was performed with Fiji plugin Jacop