

Apr 23, 2020

# 🌐 Immunohistochemical labelling of spinal cord neurons involved in bladder activity

📁 In 1 collection

DOI

[dx.doi.org/10.17504/protocols.io.bakkicuw](https://dx.doi.org/10.17504/protocols.io.bakkicuw)

Janet R Keast<sup>1</sup>, Peregrine B Osborne<sup>1</sup>, Nicole Wiedmann<sup>1</sup>

<sup>1</sup>University of Melbourne


SPARC

Tech. support email: [info@neuinfo.org](mailto:info@neuinfo.org)



Janet R Keast

University of Melbourne



OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.bakkicuw](https://dx.doi.org/10.17504/protocols.io.bakkicuw)

**Protocol Citation:** Janet R Keast, Peregrine B Osborne, Nicole Wiedmann 2020. Immunohistochemical labelling of spinal cord neurons involved in bladder activity. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bakkicuw>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** December 17, 2019

**Last Modified:** November 27, 2023

**Protocol Integer ID:** 31084

**Keywords:** activity-mapping, neuroanatomy, immunohistochemistry



## Abstract

This protocol is used for immunohistochemical visualisation of immediate early gene expression (c-Fos or Egr-1) in cryosections of rat lumbosacral spinal cord. Free-floating sections are processed in a double labelling protocol to distinguish immediate early gene expression in different neurochemical classes of spinal cord neurons:

- ChAT [choline acetyltransferase]: preganglionic neurons
- TH [tyrosine hydroxyls]: dopaminergic neurons
- Pax2: inhibitory interneurons

## Materials

### MATERIALS

- ⊗ Horse serum **Merck MilliporeSigma (Sigma-Aldrich) Catalog #12449C**
- ⊗ OCT (Optimal Cutting Temperature compound) **Sakura Finetek Catalog #4583**
- ⊗ Cy3 Donkey anti-goat IgG **Jackson ImmunoResearch Laboratories, Inc. Catalog #705-165-147**
- ⊗ AF488 Donkey anti-rabbit IgG **Jackson ImmunoResearch Laboratories, Inc. Catalog #711-545-152**
- ⊗ NeuroTrace®; 640/660 Deep-Red Fluorescent Nissl Stain - Solution in DMSO **Thermo Fisher Catalog #N21483**
- ⊗ Mouse anti-cFos antibody **Santa Cruz Biotechnology Catalog #sc166940**
- ⊗ Rabbit anti-Pax2 antibody **Invitrogen - Thermo Fisher Catalog #71-6000**
- ⊗ Rabbit anti-Egr-1 (588) antibody **Santa Cruz Biotechnology Catalog #sc-110**
- ⊗ Rabbit anti-TH antibody **Merck Millipore (EMD Millipore) Catalog #AB152**
- ⊗ Cy3 Donkey anti-mouse IgG **Jackson ImmunoResearch Laboratories, Inc. Catalog #715-165-150**
- ⊗ Goat anti-ChAT antibody **Merck Millipore (EMD Millipore) Catalog #AB144P**

### Solutions:

- PBS: phosphate-buffered saline, 0.1 M, pH 7.2
- PBS containing 0.1% sodium azide
- PBS containing 30% sucrose (w/v)
- Blocking solution: PBS containing 10% normal horse serum and 0.5% triton X-100
- PBS containing 0.1% sodium azide, 2% normal horse serum and 0.5% triton X-100

### Primary Antibodies:

	Abbreviation	Gene name	Synonym	RRID	Host Species	Dilution
	cFos	Fos	cFos	AB_10609634	Mouse	1:100
	ChAT	Chat	Choline acetyltransferase	AB_2079751	Goat	1:500



Egr-1/ Zif268	Egr1/ Zif268	Early growth h receptor 1	AB_2097174	Rabbit	1:5000
Pax2	Pax2	Paired box gene 2	AB_2533990	Rabbit	1:1000
TH	Th	Tyrosine hydroxylase	AB_390204	Rabbit	1:2000

## Secondary Antibodies:

Tag- antibody	Host species	Dilution
Cy3 anti-mouse	Donkey	1:2000
AF488 anti-rabbit	Donkey	1:1000
Cy3 anti-goat	Donkey	1:1000

## Preparation of cryosections

- 1 Cryoprotect fixed tissue (L5-S2 spinal cord) in phosphate-buffered saline (PBS; 0.1 M, pH7.2) containing 30% sucrose. This should be performed at 4C, 24-72h prior to cutting.
- 2 Embed tissue in cryomold using OCT, freeze in cryostat and cut sections (40  $\mu$ m), collecting sections progressively across sets of 4 wells to collect 160  $\mu$ m spaced series.

## Immunostaining

- 3 Wash sections in PBS (3  $\times$  10 min)
- 4 Incubate sections in blocking solution at room temperature for 2 h
- 5 Incubate sections in appropriate dilutions of primary antibodies (or combinations of primary antibodies) for 48-72h. Antibodies are diluted in PBS containing 0.1% sodium azide, 2% horse serum, and 0.5% triton-X.
- 6 Wash sections in PBS (3  $\times$  10 min)
- 7 Incubate sections in appropriate dilutions of secondary antibodies (or combinations of secondary antibodies) 4 h in the dark. Antibodies are diluted in PBS containing 2% horse serum, and 0.5% triton-X.

### Note

A useful counterstain to visualise spinal regions can be included here, by adding NeuroTrace (fluorescent Nissl stain; 1:100) to the secondary antibody incubation.

- 8 Wash sections in PBS (3  $\times$  10 min)
- 9 Mount sections onto glass slides and coverslip in preferred anti-fade mountant.



## Microscope

- 10 Labeled neurons are counted and classified according to their immunoreactivity, including only nucleated neuronal profiles in the analysis.

### Note

For digital analysis, tile-scanning of complete spinal cord sections is recommended, ensuring that the order of sections (rostral to caudal) is noted.