



Nov 17, 2017

Immunohistochemical staining

DOI

dx.doi.org/10.17504/protocols.io.kujcwun

Xiao Bing Tang¹, Li Huan¹, Jin Zhang¹, Wei Lin Wang¹, Zheng Wei Yuan², Yu Zuo Bai¹

¹Department of Pediatric Surgery, Shengjing Hospital, China Medical University;

²The Key Laboratory of Health Ministry for Congenital Malformation



Xiaobing Tang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kujcwun>

Protocol Citation: Xiao Bing Tang, Li Huan, Jin Zhang, Wei Lin Wang, Zheng Wei Yuan, Yu Zuo Bai 2017. Immunohistochemical staining. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kujcwun>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: November 17, 2017

Last Modified: March 23, 2018

Protocol Integer ID: 8811

Keywords: immunohistochemical staining, staining

Troubleshooting

- 1 Slides are deparaffinized by immersion in xylene 2 times for 2 hr.
- 2 Slides are rehydrated by immersion in 100% ethanol (10 min.), in 95% ethanol (10 min), in 85% ethanol (10min). Slides are washed with PBS (3 times, 3 min).
- 3 Antigen retrieval is performed by heating the slides in 10 mM citrate buffer (pH 6.0) at 98°C for 10 min in a microwave oven. Slides are washed with PBS (3 times, 3 min).
- 4 The endogenous peroxidase activity was blocked by incubation in 3% H₂O₂ for 20 min. Slides are washed with PBS (3 times, 3 min).
- 5 Antibody incubations were performed in phosphate-buffered saline (PBS) supplemented with 10% goat serum for 20 min at room temperature.
- 6 Slides are incubated 16hr at 4°C in a humidifying box with primary antibody. Negative controls are made with PBS alone.
- 7 Slides are washed with PBS (3 times, 3 min).
- 8 Slides are incubated with biotinylated secondary antibody for 30 min in a humidifying box at room temperature.
- 9 Slides are washed with PBS (3 times, 3 min).
- 10 Slides are incubated with horseradish peroxidase (HRP)-conjugated secondary antibody (Santa Cruz Bio-technology) for 20 min in a humidifying box at room temperature.
- 11 Slides are washed with PBS (3 times, 3 min).
- 12 Signals were visualized by using 3 '3Pdiaminobenzidine (DAB; Sigma, UK) for 1 minute and terminated by incubated in distilled water.
- 13 Slides are washed with distilled water and counterstained with hematoxylin for 10 seconds.
- 14 Slides are washed with running water for 6 min.



- 15 Slides are dehydrated by immersion in 70% ethanol (1 min), in 95% ethanol (5 min), in 100% ethanol (2 times, 5 min) and in xylene (2 times, 20 min).
- 16 Cover slides are mounted with gum.