




Feb 05, 2024

Immunohistochemical staining, vibratome sections

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2lypnkqlx9/v1>



Joris Van Asselberghs¹, Veerle Baekelandt¹

¹KU Leuven



Joris Van Asselberghs

KU Leuven

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2lypnkqlx9/v1>

Protocol Citation: Joris Van Asselberghs, Veerle Baekelandt 2024. Immunohistochemical staining, vibratome sections. protocols.io <https://dx.doi.org/10.17504/protocols.io.eq2lypnkqlx9/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 07, 2020

Last Modified: May 31, 2024

Protocol Integer ID: 44295

Keywords: vibratome, sections, brain, neurological, IHC, staining, Immunohistochemical, ASAPCRN, general procedures for immunohistochemical staining, immunohistochemical staining, vibratome section, using vibratome section, vibratome sections this protocol

Abstract

This protocol outlines general procedures for Immunohistochemical staining, using vibratome sections.

Attachments



Immunohystochemical_



139KB



Materials

Reagents Needed

- 3% H₂O₂ in 1 x PBS
- 1 x PBS 0.1% Triton X-100
- 1 x PBS 0.1% Triton X-100 + 10% serum (accordingly to the secondary Abs)
- Primary Antibodies (TH Ab152 1/5000)
- Biotinylated Secondary Antibodies (SAR 1/300)
- Streptavidine-HRP complex (1/1000)
- DAB + H₂O₂ (1.4 µl H₂O₂ for 5 ml filtered DAB solution)
 - DAB solution: 10 mg DAB (=1 tablet) for 25 ml 0.05 TRIS (TBS) pH 7.6; dissolve and filter through 0.22 µm filter, add H₂O₂ just before
- PBS or PBS + 0.1% Na Azide
- ½ PBS + ½ AD
- Ethanol:
 - 70% ethanol
 - 90% ethanol
 - 100% ethanol
- HistoClear II
- DPX mountant

Consumables Needed

- Coverslips
- Slides

Equipment Needed

- Wobbler

Safety warnings

 Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.

Before start

This protocol uses floating sections in 24-well plate in `[M] 1 X PBS` .



Quench Endogenous Peroxidase Activity

20m

- 1 Incubate samples in 1 x PBS 3% H₂O₂, for 00:10:00 at Room temperature on the wobbler (to quench endogenous peroxidase activity) using 500 µL / well. 10m 
- 2 Rinse with 0.1 % (v/v) Triton X-100 in 1x PBS . 
- 3 Rinse with 0.1 % (v/v) Triton X-100 in 1x PBS for 00:05:00 on wobbler. 5m 
- 4 Repeat: rinse with 0.1 % (v/v) Triton X-100 in 1x PBS for 00:05:00 on wobbler. 5m 

Primary Antibody Incubation

40m

- 5 Add 250 µL Primary Abs in 1 x PBS 0.1% Triton X-100 + 10% serum (accordingly to the secondary Abs) Overnight at Room temperature on wobbler. 40m  

Note

Dilution: 1st AB: TH Ab152 1/5000

- 6 Rinse with 0.1 % (v/v) Triton X-100 in 1x PBS . 
- 7 Rinse with 0.1 % (v/v) Triton X-100 in 1x PBS for 00:05:00 on wobbler. 5m 
- 8 Repeat: rinse with 0.1 % (v/v) Triton X-100 in 1x PBS for 00:05:00 on wobbler. 5m 

Secondary Antibody Incubation

20m



- 9 Add 250 μ L biotinylated Secondary Abs in 1 x PBS 0.1% Triton X-100 for 00:30:00 at Room temperature on wobbler.

30m

**Note**Dilution: 2nd AB: SAR 1/300

- 10 Rinse with 0.1 % (v/v) Triton X-100 in 1x PBS .



- 11 Rinse with 0.1 % (v/v) Triton X-100 in 1x PBS for 00:05:00 on wobbler.

5m



- 12 Repeat: rinse with 0.1 % (v/v) Triton X-100 in 1x PBS for 00:05:00 on wobbler.

5m

**Streptavidine-HRP complex Incubation**

30m

- 13 Add 250 μ L Streptavidine-HRP complex 1/1000 (in 1 x PBS 0.1% Triton X-100) for 00:30:00 at Room temperature on wobbler.

30m

**Note**

Dilution: STRP-HRP 1/1000

- 14 Rinse with 0.1 % (v/v) Triton X-100 in 1x PBS .



- 15 Rinse with 0.1 % (v/v) Triton X-100 in 1x PBS for 00:05:00 on wobbler.

5m



- 16 Repeat: rinse with 0.1 % (v/v) Triton X-100 in 1x PBS for 00:05:00 on wobbler.

5m





DAB + H₂O₂

17

Safety information

Work on DAB Bench!

Prepare DAB solution: 10 mg DAB (=1 tablet) for 25 mL 0.05 TRIS (TBS) pH 7.6 ; dissolve and filter through 0.22 µm filter, add H₂O₂ just before use! Or Vector SG (TH). (Add 1.4 µL H₂O₂ for 5 mL filtered DAB solution).

18

Add 250 µL DAB + H₂O₂ **(ON DAB BENCH!)**. Allow reaction to proceed for a few minutes at Room temperature **without wobbling**.

19

Rinse with 0.1 % (v/v) Triton X-100 in 1x PBS .



20

Replace TBS Triton X-100 with PBS or 0.1 % (v/v) Sodium Azide in PBS in case sections are not to be mounted immediately. Store at 4 °C .

½ PBS + ½ AD Rinse

30m

21

Briefly rinse sections in ½ PBS + ½ AD and mount on gelatin coated microscopy slides. Allow to dry for 00:30:00 in flow or couple of hours on bench.

30m

Dehydration

25m

22

Dehydrate samples in 70 % (v/v) ethanol for 00:05:00 .

5m

23

Dehydrate samples in 90 % (v/v) ethanol for 00:05:00 .

5m

24

Dehydrate samples in 100 % (v/v) ethanol for 00:05:00 .

5m

25

Repeat: dehydrate samples in 100 % (v/v) ethanol for 00:05:00 .

5m



26 Replace ethanol with HistoClear II for  00:05:00 .

5m

Mounting

30m

27 Mount coverslips on top of the slides with DPX and allow to dry  Overnight in the flow.

30m



28 Press out bubbles next morning.