



Apr 03, 2025

Immunohistochemistry & Microscopy-Falasconi & Kanodia et al 2025

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g7m7k9gwz/v1>

Antonio Falasconi^{1,2,3}, Harsh Kanodia^{1,2,3}, Silvia Arber^{1,2,3}

¹Biozentrum, University of Basel; ²Friedrich Miescher Institute for Biomedical Research;

³Aligning Science Across Parkinson's (ASAP)



Christine Weber-Schmidt

Allen Institute

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7m7k9gwz/v1>

Protocol Citation: Antonio Falasconi, Harsh Kanodia, Silvia Arber 2025. Immunohistochemistry & Microscopy-Falasconi & Kanodia et al 2025. [protocols.io https://dx.doi.org/10.17504/protocols.io.q26g7m7k9gwz/v1](https://dx.doi.org/10.17504/protocols.io.q26g7m7k9gwz/v1)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 19, 2025

Last Modified: April 03, 2025

Protocol Integer ID: 120855

Keywords: ASAPCRN, kanodia et al 2025 immunohistochemistry, immunohistochemistry, falasconi, microscopy method, microscopy, kanodia et al

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-020551

European Union's Horizon 2020 research and innovation programme

Grant ID: InterAct, grant agreement No 101018151

Swiss National Science Foundation

Grant ID: N/A

Abstract

Immunohistochemistry and microscopy methods used in Falasconi & Kanodia et al 2025.

Guidelines

Anesthesia:

Animals are anesthetized with a ketamine–xylazine solution.

Materials

Antibodies:

	A	B	C	D
	Chicken anti-GFP	Invitrogen	RRID: AB_2534023	Cat# A10262
	Rabbit anti-RFP	Rockland	RRID: AB_2209751	Cat# 600-401-379
	chicken anti-TH	Neuromics	RRID: AB_2201403	Cat# CH23006
	goat anti-ChAT	Millipore	RRID: AB_2079751	Cat# AB144P
	Donkey anti-rabbit Cy3	Jackson Immuno Research	RRID: AB_2307443	Cat#711-165-152
	Donkey anti-goat Cy5	Invitrogen	RRID: AB_2535864	Cat# A-21447
	Donkey anti-chicken 488	Jackson Immuno Research	RRID: AB_2340375	Cat#703-545-155
	Donkey anti-chicken Cy5	Jackson Immuno Research	RRID: AB_2340379	Cat#703-605-155
	Donkey anti-goat 488	Invitrogen	RRID: AB_2534102	Cat# A-11055

Safety warnings

! Wear appropriate PPE as required by your institution.

Ethics statement

Prior ethics approval should be obtained before performing these experiments. Approval was obtained from the Cantonal Veterinary Office Basel-Stadt and performed in compliance with the Swiss Veterinary Law guidelines.

Before start

Read through the entire protocol before starting.

Transcardial Perfusion and Tissue Fixation

- 1 Prepare PBS and 4% PFA in PBS for perfusion.
- 2 Set up the perfusion pump and fill the tubing with PBS.
- 3 Pin the mouse to a Styrofoam plate.
- 4 Make a midline incision to expose the thoracic cavity, cutting through the skin and rib cage.
- 5 Insert the needle into the left ventricular chamber of the heart.
- 6 Switch on the perfusion pump and open the right atrium using sharp forceps or small scissors immediately.
- 7 Perfuse the mouse with PBS to flush out the blood.
- 8 Switch the perfusion to 4% PFA and perfuse until limbs are stiff.
- 9 Dissect the brain and spinal cord then place it in 10 mL of 4% PFA for post-fixation overnight at 4°C.
- 10 Transfer the brain to a 30% sucrose solution in PBS at 4°C until the tissue sinks (24–48 hours).
- 11 Embed the brain in OCT compound and freeze at -80°C.

Sectioning

- 12 Mount brain onto a cyrostat mount using OCT to be able to section in a coronal orientation.



- 13 Cut 80 μm sections and place sections in 24 well plates with PBS

Staining (IHC)

- 14 Incubate floating sections for 1 hour in blocking solution (1% BSA, 0.2% Triton X-100, PBS) at room temperature.
- 15 Add primary antibodies to the blocking solution and incubate for 1-3 days at 4°C.
- 16 Floating sections are washed x3 with PBS for 10 minutes each wash.
- 17 Add fluorophore-coupled secondary antibodies in blocking solution and incubate for 1 day at 4°C.
- 18 Floating sections are washed x3 with PBS for 10 minutes each wash.
- 19 Mount sections using fluorescence-compatible mounting medium and carefully place coverslips.

Imaging

- 20 For low-resolution overview imaging: scan slides with an Axioscan light microscope (Zeiss).
- 21 For higher resolution imaging use an Axio Imager M2 microscope (Zeiss) with a Yokogawa CSU W1 Dual camera T2 spinning disk confocal scanning unit.