



Jan 08, 2026

Immunohistochemistry and immunofluorescence protocols

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov1ko5kgr2/v1>

Junyu Xiang¹

¹Department of Gastroenterology, Chongqing Key Laboratory of Digestive Malignancies, Daping Hospital, Chongqing 400042, China



Junyu Xiang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1ko5kgr2/v1>

Protocol Citation: Junyu Xiang 2026. Immunohistochemistry and immunofluorescence protocols. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.ewov1ko5kgr2/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: January 08, 2026

Last Modified: January 08, 2026

Protocol Integer ID: 238193

Keywords: immunofluorescence protocol, immunohistochemistry, immunofluorescence, phosphorylation level, indicated protein, phosphorylation, ihc, such as nprl2

Abstract

A detailed description for immunohistochemistry (IHC) and immunofluorescence (IF) experiments. The expected results can clearly demonstrate the indicated proteins and their phosphorylation levels, such as NPRL2 and pS235/236-S6.

Immunohistochemistry protocol

- 1 Formalin-fixed paraffin-embedded specimens or organoids are cut into 5-10 μm thick sections.
- 2 The sections are deparaffinized, rehydrated, and blocked with 3% H_2O_2 for 15 min.
- 3 To remove melanin in immunohistochemistry of melanoma tumor tissues, sections are emersed in 0.5% permanganate potassium solution at 55°C for 3 minutes, then treated with 1% oxalic acid solution for 3 cycles of 30 seconds followed by water rinse for 1 minute.
- 4 Tissue antigens are retrieved by heating the sections in a pressure cooker for 150 sec.
- 5 After cooling to room temperature, each section is incubated with primary antibody at 4°C overnight, and then with goat secondary antibody against rabbit and mouse immunoglobulins (DAKO).
- 6 For visualization, 1-5% 3,3'-diaminobenzidine tetrahydrochloride is used as a chromogen.
- 7 After washing with water, the sections are counterstained with hematoxylin.
- 8 Negative control slides with the primary antibody replaced with phosphate-buffered saline (PBS) are included.
- 9 All immunostained sections are independently evaluated by two pathologists blinded to the patients' clinical status, with interobserver discrepancy being less than 10%. In the case of discrepancy, a consensus interpretation is reached with a third pathologist.
- 10 The quantification method is based on cell number counting, or multiplicative index of the average staining intensity (score 0-3) and extent of staining (score 0, 0%; 1, 0-25%; 2, 25-50%; 3, 50-75%; 4, 75-100%) in the cores (such as pS6 IHC score), yielding a staining index ranging from 0 (no staining) to 12 (extremely extensive, extremely strong staining).

Immunofluorescence protocol



- 11 Tissue sections are permeabilized with PBS-T (containing 0.25% Triton X-100) for 30 minutes and blocked with 5% goat serum in PBS for 60 minutes.
- 12 Primary antibodies are diluted in TBS-T (containing 0.1% Tween-20) (dilution 1: 25~200) and incubated at 4°C overnight, followed by incubation with appropriate secondary antibodies
- 13 Coverslips are mounted on slides using anti-fade mounting medium with DAPI.
- 14 Images are acquired with a ZEISS 880 Microscope.
- 15 For multiple-immunofluorescence for samples array, Opal 7-color Manual IHC Kit (NEL801001KT, PerkinElmer) is used and performed according to manufactures instructions.
- 16 **Fluorescence signals** were detected by Vectra Polaris Automated Quantitative Pathology Imaging System (PerkinElmer) and the data was analyzed by StrataQuest (ver. 7.0.1.176, TissueGnostics).
- 17 Pearson's correlation in cells were quantified by the Coloc2 plugins in ImageJ (ver. 1.53q).

