


May 29, 2023

🌐 Immunohistochemistry for CD8+ staining in Breast Cancer Tissue

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l69odrlqe/v1>



Freda Halim¹

¹Faculty of Medicine, Pelita Harapan University, Tangerang



Freda Halim

Faculty of Medicine, Pelita Harapan University, Tangerang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l69odrlqe/v1>

Protocol Citation: Freda Halim 2023. Immunohistochemistry for CD8+ staining in Breast Cancer Tissue. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bp2l69odrlqe/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: May 28, 2023

Last Modified: May 29, 2023

Protocol Integer ID: 82572

Keywords: using primary antibody cd8, primary antibody cd8, study of cd8, slides for cd8, immunohistochemistry, cd8, breast cancer tissue, using primary antibody, primary antibody, antibody, cell count, negative bc patient, staining

Abstract

These are protocols used for study of CD8+ Cell Count in Luminal B Her-2 negative BC patients. We used using paraffin sections of 66 samples and stained the slides for CD8+ antibody, using primary antibody CD8 (clone C8/144B, dilution 1:100; Dako



Deparaffinization and rehydration

23m

- 1 Incubate slides in Xylenes for 3 minutes 9m
- 2 Rehydrate slides in 100% Ethanol, 96% Ethanol, 70% Ethanol each for 3 minutes 9m
- 3 Rinse with running tap water and aquadest 5m

Blockage of Endogenous Peroxidase

20m

- 4 Incubate slides in 3% H₂O₂ for 15 minutes 15m
- 5 Rinse with running tap water and aquadest 5m

Antigen Retrieval

1h

- 6 Antigen Retrieval with Tris EDTA (pH9) with pressure cooker, in 95 degree Celcius temperature 20m
- 7 Open the lid and cool down in room temperature 15m
- 8 Rinse slides with running water and aquadest for 5 minutes 5m
- 9 Rinse in PBS (Phosphate Buffer Saline) in pH 7.40-7.60 5m
- 10 Excell Block 10m
- 11 Rinse in PBS in pH 7.40-7.60 5m



Primary Antibody

1h 5m

12 Wipe excess liquid around the tissue, apply primary antibody :CD8 (clone C8/144B, dilution 1:100; Dako

13 Incubate for 60 minutes and then rinse with PBS for 5 minutes

1h 5m

Secondary Antibody

40m

14 Apply Excell Link as secondary antibody for 15 minutes, then rinse with PBS for 5 minutes

20m

15 Apply Excell HRP as secondary antibody

20m

Signal Detection/Histochemistry

27m

16 Apply DAB (Diamino-benzidine) 80-100 μ L for 10 minutes , Rinse with running tap water and aquadest for 5 minutes

15m

17 Apply Hematoxyline for 1 minutes, Rinse with running tap water and aquadest for 5 minutes

6m

18 Apply Tatcha's bluing solution and rinse with running tap water and aquadest for 5 minutes

6m

Dehydration and Clearing

20m

19 Clear excess water from slides and Dehydrate Slides in 70%, 96%, 100% for 5 minutes each

15m

20 Incubate slides in Xylenes for 5 minutes

5m

21 Mount the slides



Protocol references

Asano Y, Kashiwagi S, Goto W, Kurata K, Noda S, Takashima T, et al. Tumour-infiltrating CD8 to FOXP3 lymphocyte ratio in predicting treatment responses to neoadjuvant chemotherapy of aggressive breast cancer. *Br J Surg.* 2016;103(7):845–54.