




Mar 20, 2023

Immunohistochemistry for p53 staining in Breast Cancer Tissue

DOI

dx.doi.org/10.17504/protocols.io.5qpvard7dv4o/v1



Freda Halim^{1,2}

¹Doctoral Degree, Faculty of Medicine, Universitas Padjajaran;

²Surgery Department, Faculty of Medicine, Pelita Harapan University



Freda Halim

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvard7dv4o/v1>

Protocol Citation: Freda Halim 2023. Immunohistochemistry for p53 staining in Breast Cancer Tissue. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5qpvard7dv4o/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: March 13, 2023

Last Modified: March 20, 2023

Protocol Integer ID: 78640

Keywords: immunohistochemistry for p53, study of p53 expression difference, p53 expression difference, p53 expression, slides for p53 expression, breast cancer tissue, p53, immunohistochemistry, using primary antibody clone, primary antibody clone, cancer, primary endocrine therapy resistance, antibody

Abstract

These are protocols used for study of p53 expression differences in Luminal B Her-2 negative patients. The aim of the study is to show p53 expression differences in Luminal B Her-2 negative patients with and without primary endocrine therapy resistance. We used using paraffin sections of 67 samples and stained the slides for p53 expression, using primary antibody Clone DO-7 (DAKO).

Troubleshooting



Deparaffinization and Rehydration

- 1 Incubate slides in Xylenes for 3 minutes 3m
- 2 Incubate slides in Xylenes for 3 minutes 3m
- 3 Incubate slides in Xylenes for 3 minutes 3m
- 4 Rehydrate slides in 100% Ethanol for 3 minutes 3m
- 5 Rehydrate slides in 96% Ethanol for 3 minutes 3m
- 6 Rehydrate slides in 70% Ethanol for 3 minutes 3m
- 7 Rinse with running tap water and aquadest for 5 minutes 5m

Blockage of Endogenous Peroxidase

- 8 Incubate slides in 3% H₂O₂ for 15 minutes 15m
- 9 Rinse slides with running tap water and aquadest for 5 minutes 5m

Antigen Retrieval

40m

- 10 Antigen Retrieval
- 10.1 Antigen Retrieval with Tris EDTA (pH9) with pressure cooker, in 95⁰ Celcius temperature 20m



10.2 Open the lid and cool down in room temperature

15m

10.3 Rinse slides with running water and aquadest for 5 minutes

5m

10.4 Rinse in PBS (Phosphate Buffer Saline) in pH 7.40-7.60

5m

10.5 Excell Block

10m

10.6 Rinse in PBS in pH 7.40-7.60

5m

Primary Antibody

40m

11 Wipe excess liquid around the tissue

12 apply primary antibody (clone DO7, Dako) 120μL

13 Incubate for 60 minutes

1h

14 Rinse with PBS

5m

Secondary Antibody

40m

15 Apply Excell Link as secondary antibody

15m

16 Rinse with PBS

5m



17 Apply Excell HRP as secondary antibody

20m

Signal Detection/ Histochemistry

27m

18 Apply DAB (Diamino-benzidine) 80-100 μ L for 10 minutes , Rinse with running tap water and aquadest for 5 minutes

15m

19 Apply Hematoxylline for 1 minutes, Rinse with running tap water and aquadest for 5 minutes

6m

20 Apply Tatcha's bluing solution and rinse with running tap water and aquadest for 5 minutes

6m

Dehydration and Clearing

20m

21 Clear excess water from the slides

22 Dehydrate Slides in 70%, 96%, 100% for 5 minutes each

15m

23 Incubate slides in Xylenes for 5 minutes

5m

24 Mount the slides

Protocol references

Kikuchi, S., Nishimura, R., Osako, T., et al. Definition of p53 Overexpression and its Association with the Clinicopathological Features in Luminal/HER2-negative Breast Cancer. *ANTICANCER RESEARCH* 33: 3891-3898 (2013)