


Mar 25, 2024

# 🌐 Immunohistochemistry (IHC) Staining Mouse Brain Sections

DOI

[dx.doi.org/10.17504/protocols.io.5qpvo3b7bv4o/v1](https://dx.doi.org/10.17504/protocols.io.5qpvo3b7bv4o/v1)



Naveen Ouellette<sup>1</sup>, Daphne Toglia<sup>1,2</sup>, Holly Myers<sup>1</sup>

<sup>1</sup>Allen Institute for Neural Dynamics; <sup>2</sup>Boston University

Allen Institute for Neural...



Holly Myers

Allen Institute

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvo3b7bv4o/v1>

**Protocol Citation:** Naveen Ouellette, Daphne Toglia, Holly Myers 2024. Immunohistochemistry (IHC) Staining Mouse Brain Sections. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5qpvo3b7bv4o/v1>



**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** July 26, 2023

**Last Modified:** June 18, 2024

**Protocol Integer ID:** 85540

**Keywords:** Histology, Immunohistochemistry, IHC, Immunolabeling, Antibody, mouse brain sections the immunohistochemistry, staining mouse brain section, micron mouse brain tissue slice, mouse brain sections protocol detail, immunohistochemistry, secondary antibody stain, ihc, mouse brain, secondary antibody, staining step, staining duration

## Abstract

The Immunohistochemistry (IHC) Staining for Mouse Brain Sections protocol details the blocking, primary, and secondary antibody staining of 50-100 micron mouse brain tissue slices fixed in 4% PFA. The protocol includes suggested staining duration based on slice thickness for each staining step and tables of the most frequently requested blocking serums and primary and secondary antibody stains and dilutions.


## Materials

 Sample sections to be stained:

- 50-100 microns thick
- Fixed in paraformaldehyde
- Free floating, stored in PBS or PBS + 0.01% Azide
- Protected from light

### Reagents:

 10xPBS **Ambion Catalog #AM9624**


 Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8787-50ML**


 Sodium Azide 5% **Ricca Chemical Company Catalog #71448-16**

 Milli-Q water

### Optional reagents - determined by stain requested:

- Note these reagents and stained tissue must be protected from light

 DAPI (4',6-diamidino-2-phenylindole, dihydrochloride) **Thermo Fisher Catalog #62247**

 Primary Antibody

 Secondary Antibody

 Normal Goat Serum **Vector Laboratories Catalog #S-1000-20**

 Normal Donkey Serum **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D9663-10ML**

 Block Aid **Invitrogen - Thermo Fisher Catalog #B10710**

 Urea powder **Merck MilliporeSigma (Sigma-Aldrich) Catalog #U5378-100G**



	Materials	Product number
	Serological pipet filler	Thermofisher, 9501
	48 well plate	Costar, 3548
	Cell culture inserts	Netwell, 734-1589
	Serological pipets	Fisher Scientific, 13-678-11E
	Nutating mixer	Fisherbrand, 88-861-043
	Stir plate	Merck MilliporeSigma, Z693510



	Materials	Product number
	Manual single channel pipettes: P20, P200, P1000	Rainin, 30456871
	Manual single channel pipettes: P2	Rainin, 17014393
	Manual dingle channel pipettes: P5000	Rainin, 17011790
	Stir bar	Grainger, 21R590




## Recipes - blocking and DAPI solutions should be made depending on experimental need

### 1L 1xPBS:

Combine the following reagents into a container with a stir bar. A graduated cylinder may be used to measure MilliQ water and a graduated cylinder or serological pipet may be used to measure 10xPBS. Mix well on a stir plate at high speed (300 RPM or higher) for  00:02:00 or until solution is mixed. Store at  Room temperature for 1 month.

	Reagent	Volume
	Milli-Q water	900 mL
	10xPBS	100 mL

### 1L 1xPBS & Sodium Azide 0.01% :



Combine the following reagents into a container with a stir bar, using a graduated cylinder to measure the 1xPBS and a P5000 pipette to measure the 5% sodium azide. Mix well on a stir plate at high speed (300 RPM or higher) for  00:02:00 or until solution is mixed. Store at  Room temperature or  4 °C for up to 1 year.

	Reagent	Volume
	1xPBS	998 mL





	Reagent	Volume
	0.01% Sodium Azide	2 mL

**1L 1xPBS & Triton X-100 0.2%:**

Measure 1xPBS into a container with a stir bar using a graduated cylinder and slowly pipet Triton X-100 in using a P5000 pipette. Mix well on a stir plate at high speed (300 RPM or higher) for  00:15:00 to ensure Triton X-100 goes into solution. Store at  Room temperature for 1 month.


	Reagent	Volume
	1xPBS	998 mL
	Triton X-100	2 mL

**1L 1xPBS & Triton X-100 0.06%:**

Measure 1xPBS into a container with a stir bar using a graduated cylinder and slowly pipet Triton X-100 in using a P5000 pipette. Mix well on a stir plate at high speed (300 RPM or higher) for  00:15:00 to ensure Triton X-100 goes into solution. Store at  Room temperature for 1 month.



	Reagent	Volume
	1xPBS	994 mL
	Triton X-100	6 mL

**2mL 5mg/mL DAPI solution:**

Add Milli-Q water to DAPI powder in 10mg vial using P5000 pipette. Vortex until powder completely mixes into solution. Store at  4 °C and vortex before use.

	Reagent	Volume
	Milli-Q water	2 mL
	DAPI powder	10 mg



### 200 mL 5% Normal Goat Serum & Triton X-100 0.06% & 4M Urea (NGSTU)

Measure urea on scale and add into a container for mixing NGSTU. Add 10x PBS into the container with a serological pipet, then add normal goat serum, Triton and water to container using manual pipette. Add stir bar to container and mix on stir plate at high speed (300 RPM or higher) for  00:15:00 to ensure all reagents fully homogenize. This working solution can be stored at  4 °C and used for up to 1 week. Beyond this time, there would be concern for possible microbial contamination that would affect performance.

\*Amount of blocking serum made should be determined by number of sections to be stained.

	Reagent	Amount
	Normal Goat Serum	10 mL
	Urea	48 g
	10X PBS	20 mL
	Triton X-100	120 uL
	Milli-Q Water	Fill to 200 mL



### 200 mL 5% Normal Donkey Serum & Triton X-100 0.2% in PBS (NDST)

Measure 10X PBS into container with a serological pipet, then add Normal Donkey Serum, Triton, and water to container using manual pipette. Add stir bar to container and mix on stir plate at high speed (300 RPM or higher) for  00:15:00 to ensure all reagents fully homogenize. This working solution can be stored at  4 °C and used for up to 1 week. Beyond this time, there would be concern for possible microbial contamination that would affect performance.

\*Amount of blocking serum made should be determined by number of sections to be stained.

	Reagent	Amount
	Normal Donkey Serum	10 mL
	10X PBS	20 mL
	Triton X-100	400 uL
	Milli-Q Water	Fill to 200 mL

### 200 mL 5% Normal Goat Serum & Triton X-100 0.06% (NGST):

Measure 10X PBS into container with a serological pipet, then add normal goat serum, Triton, and water to container using manual pipette. Add stir bar to container and mix on stir plate at high speed (300 RPM or higher) for  00:15:00 to ensure all reagents fully homogenize. This working solution can be stored at  4 °C and used for up to 1 week. Beyond this time, there would be concern for possible microbial contamination that would affect performance.

\*Amount of blocking serum made should be determined by number of sections to be stained.

	Reagent	Amount
	Normal Goat Serum	10 mL
	10X PBS	20 mL
	Triton X-100	400 uL
	Milli-Q Water	Fill to 200 mL

## Troubleshooting



## Safety warnings

⚠ DAPI is a mutagen and should be handled with care. Wear PPE and dispose into hazardous waste stream. Please consult your immediate supervisor or the EH&S manager/representative if you have questions or concerns.

Sodium Azide is toxic and carcinogenic. It should be handled and prepared with care. Do not breathe dust, do not use metal utensils. Wear gloves when handling this chemical.

Paraformaldehyde is carcinogenic. Wear gloves at all times when handling specimen tissue slices fixed in paraformaldehyde, as the tissue may contain trace amounts of the chemical.

Personal Protective Equipment (PPE) should be used at all times while operating this protocol. If you are unsure what PPE you should be using, see your immediate supervisor.

## Before start

This protocol details staining of mouse brain tissue that has already been sectioned. Reference the following protocol for instructions on sectioning a whole mouse brain fixed in 4% PFA: **[Sectioning Mouse Brain with Sliding Microtome](#)**

For all steps, protect tissue from light in order to preserve fluorescence by wrapping well plate in foil between steps and during washes and stains.





## Setup

### 1 Well plate setup:

#### 1.1 Select either a plastic well plate or Netwell cell culture inserts for tissue section staining.

##### Note

A plastic well plate may be helpful for conducting multiple different stains simultaneously, while a Netwell kit may be efficient for applying the same stain to a large number of tissue sections.

#### 1.2 Fill wells with 1xPBS and place tissue sections to be stained into wells.

If there are a large amount of tissue sections, the sections may be double or triple stacked in each well in the well plate (or more if using Netwell insert), but ideally sections in the same well should be separated by at least a few hundred microns (ex: every 6th section or more) so it is easier to see anatomical differences and determine section order when mounting.


#### 1.3 If tissues will be stained in different conditions (ex: staining sections with different primary or secondary antibodies) it is best to draw a map of these conditions on the plastic well plate lid and a corresponding map in lab notebook.

## Staining

**3d 0h 10m**

### 2 Wash mouse brain sections three times in 1xPBS for 00:05:00 at

**15m**

 Room temperature on shaker using either a plastic well plate or Netwell kit with mesh inserts.

##### Note


In this step and all subsequent steps, protect tissue from light to prevent bleaching by wrapping well plate in foil while it is on the shaker.

### 3 Blocking:

### 3.1 Select blocking solution. See examples of frequently requested blocking solutions below:


	Blocking solution	Notes	Example
	BlockAid	Can be paired with any secondary antibody	Example: BlockAid + chicken anti-GFP (primary antibody) + (goat OR donkey) anti-chicken (secondary antibody)
	Normal goat serum with triton (NGST)	Must match secondary antibody	Example: NGST (blocking solution) + chicken anti-GFP (primary antibody) + goat anti-chicken (secondary antibody)
	Normal donkey serum with triton (NDST)	Must match secondary antibody	Example: NDST (blocking solution) + chicken anti-GFP (primary antibody) + donkey anti-chicken (secondary antibody)
	Normal goat serum with triton and urea (NGSTU)	Must match secondary antibody, urea helps with background	Example: NDSTU (blocking solution) + chicken anti-GFP (primary antibody) + donkey anti-chicken (secondary antibody)

### 3.2 Choose a blocking duration based on section thickness. See table below. Note that suggested blocking durations are minimums and sections may be blocked up to

 Overnight

	Section thickness	Blocking duration
	50 microns	1 hour
	100 microns	2 hours

### 3.3 Wash sections in chosen blocking solution for chosen blocking duration on shaker at


 Room temperature

## 4 Primary staining:

### 4.1 Choose a primary antibody and dilution based on tissue labeling. See table below for commonly requested primary antibodies and dilutions. If using NDST, NGST, or NGSTU as blocking serum, dilute the primary antibody in tube with blocking serum. If using Block



Aid as blocking serum, dilute the primary antibody in tube with 1xPBS & Triton X-100 0.2%. Mix primary antibody dilution mixture in tube by gently swirling.

Frequently requested primary antibodies	Suggested dilution
Immunostar mouse anti-Tyrosine hydroxylase (22941)	1:1000
NovusBio rabbit anti-Tph2 antibody (NB100-74555)	1:1000
Rockland rabbit anti-RFP antibody pre-adsorbed (600-401-379)	1:1000 or 1:800
Invitrogen mouse anti-HA (26183)	1:500
AVES chicken anti-GFP (GFP-1020)	1:800
Immunostar goat anti-5HT (26183)	1:800
FujiFilm Wako rabbit anti-IBA1 (019-19741)	1:1000
Sigma mouse anti-GFAP (G6171)	1:1000

- 4.2 Chose a staining duration based on section thickness. Note that the durations below are minimums and sections may be stained up to  72:00:00

3d

Section thickness	Primary antibody staining duration minimums
50 microns	24+ hours
100 microns	48+ hours

- 4.3 Remove blocking solution and incubate tissue in chosen primary antibody solution on shaker at  4 °C or  Room temperature for chosen staining duration.

#### Note

**If no secondary antibody stain is required, skip to step 7.**



- 5 Wash sections in new wells containing 1xPBS & Triton X-100 0.2% five or more times for

00:05:00 on shaker at Room temperature

25m

- 6 Secondary staining:

- 6.1 Choose a secondary antibody and dilution based on the primary antibody used. See table below for commonly requested secondary antibodies and dilutions. If using NDST, NGST, or NGSTU as blocking serum, dilute the secondary antibody in tube with blocking serum. If using Block Aid as blocking serum, dilute the secondary antibody in tube with 1xPBS & Triton X-100 0.2%. Mix secondary antibody dilution mixture in tube by gently swirling.

Frequently requested secondary antibodies	Suggested dilution
DAPI (5 mg/mL stock concentration, then diluted to working concentration) combine with secondary antibody	1:1000 or 1:5000
Invitrogen goat anti-mouse 488 cross-adsorbed (A11001)	1:500
Invitrogen goat-anti-rabbit 488 (A11012)	1:500
Invitrogen goat anti-mouse 647 (A-21236)	1:500
Invitrogen goat anti-rabbit 405 (A-31556)	1:500
Invitrogen goat anti-rabbit 647 (A21244)	1:500
Jackson Immuno goat anti-chicken 488 (703-545-155)	1:500
Invitrogen goat anti-rabbit 594 (A-11012)	1:500
Invitrogen donkey anti-goat 647 (A-21447)	1:500

- 6.2

#### Note

**If additional DAPI stain is required in combination with secondary antibody, continue with step 6.2. If not, skip to step 6.3.**



Vortex  DAPI solution and dilute it to preferred concentration in secondary antibody solution. 1:5000 is a common dilution for DAPI.

#### Safety information

DAPI is a mutagen, so wear PPE (gloves, lab coat, safety goggles) and dispose in hazardous waste stream.

- 6.3 Choose a staining duration based on section thickness. Suggestion below are minimums and secondary staining duration may run up to  :

2d

	Section thickness	Secondary staining duration minimums
	50 microns	1+ hour
	100 microns	2+ hours

- 6.4 Incubate sections in chosen secondary antibody solution for chosen staining duration on shaker at

- 7 Wash sections in 1xPBS & Triton X-100 0.2% three times for  on shaker at

25m

- 8 Wash sections in 1xPBS two times for  on shaker at

30m

- 9 Store samples in 1xPBS for short term use or 1xPBS & Azide 0.01% for long term use (  or up to several weeks) at  until ready for mounting.



### Safety information

Sodium Azide is toxic and carcinogenic. It should be handled and prepared with care. Do not breathe dust, do not use metal utensils. Wear gloves when handling this chemical.

### Note

Proceed to **Mounting and Coverslipping Mouse Brain Sections** protocol.