



Mar 16, 2024

Immunohistochemistry (neural organoids)

DOI

dx.doi.org/10.17504/protocols.io.n92ldp22nl5b/v1

anita.adami¹

¹Laboratory of Molecular Neurogenetics, Department of Experimental Medical Science, Wallenberg Neuroscience Center and Lund Stem Cell Center, BMC A11, Lund University, 221 84 Lund, Sweden.



Anita Adami

Lund University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldp22nl5b/v1>

Protocol Citation: anita.adami 2024. Immunohistochemistry (neural organoids). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n92ldp22nl5b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 03, 2023

Last Modified: June 01, 2024

Protocol Integer ID: 76335

Keywords: ASAPCRN, neural organoid, immunohistochemistry

Funders Acknowledgements:

Aligning Science Across Parkinson's through the Michael J. Fox Foundation for Parkinson's Research

Grant ID: ASAP-000520

Swedish Research Council

Grant ID: 2018-02694

Swedish Brain Foundation

Grant ID: FO2019-0098

Cancerfonden

Grant ID: 190326

Barncancerfonden

Grant ID: PR2017-0053

NIHR Cambridge Biomedical Research Centre

Grant ID: NIHR203312

Swedish Society for Medical Research

Grant ID: S19-0100

National Institutes of Health

Grant ID: HG002385

Swedish Research Council

Grant ID: 2021-03494

Swedish Research Council

Grant ID: 2020-01660

Abstract

This protocols describes how to perform immunohistochemistry on neural organoids

Troubleshooting



Fixing and mounting

4h

- 1 The organoids were fixed in 4% paraformaldehyde for 02:00:00 at Room temperature .
- 2 They were subsequently washed three times with KPBS and left in a 1:1 OCT:30% sucrose solution and OCT (HistoLab) mixture Overnight .
- 3 The organoids were then transferred to a cryomold containing OCT and frozen on dry ice and stored at -80°C.

Sectioning

- 4 Prior to the staining protocol, the organoids were cryosectioned at a thickness of 20 µM.

Staining

2h 10m

- 5 For immunohistochemistry, sections were washed in PBS1X for 00:10:00 and then blocked and permeabilized in **Blocking solution** (0.3% Triton X-100 and 4% normal donkey serum in PBS1X) for at least 01:00:00 .
- 6 The sections were then incubated with primary antibodies in **blocking solution** at 4 °C Overnight .

The primary antibodies used for neural organoids characterisation were rabbit anti-PAX6 (1:300; BioLegend, Cat# 3700, RRID:[AB_2242334](#)) and rat anti-ZO1 NB110-68140, RRID:[AB_1111431](#)(1:300; Novus)).
- 7 After incubation with the primary antibodies, the sections were incubated for 1 h with the appropriate secondary antibodies (Alexa Fluor 488, 594, 647 used at 1:400; Molecular Probes).

Mounting

- 8 Finally, the sections were mounted on gelatin-coated slides and coverslipped with PVA-DABCO containing DAPI (1:1000).