

Mar 06, 2025

Immunohistochemistry Protocol (brain post-mortem tissue)

Forked from a private protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.n92ldroyog5b/v1>

**Immunohistochemistry
Protocol**

Roger Barker¹

¹University of Cambridge



Laura Castilla-Vallmanya

Lund University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldroyog5b/v1>

Protocol Citation: Roger Barker 2025. Immunohistochemistry Protocol (brain post-mortem tissue). protocols.io
<https://dx.doi.org/10.17504/protocols.io.n92ldroyog5b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 06, 2025

Last Modified: March 06, 2025

Protocol Integer ID: 123911

Keywords: immunohistochemistry, antibody, sequential rehydration, counterstaining, ASAPCRN, immunohistochemistry protocol, immunohistochemistry steps for antibody, immunohistochemistry step, antibody, tissue

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-00520

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-025170

Abstract

This protocol details the immunohistochemistry steps for antibody blocking over a period of 2 days.

Attachments



[Immunohistochemistry..](#)

↓

107KB



Materials

REAGENTS

- Xylene
- 70% EtOH
- 90% EtOH
- 100% EtOH
- dH₂O
- Sodium citrate buffer: 2.94g Tri-sodium citrate (dehydrate), 1000ml Distilled Water (adjust to pH 6 with 1N HCl) then add 0.5ml Tween 20
- TBPS: 1x PBS + 0.05% Tween 20
- 3% H₂O₂ diluted in TPBS
- Blocking solution (e.g. 5% Normal horse serum in TPBS)
- Primary and Secondary Antibody
- ABC solution (DAB substrate): 5ml TPBS with 2 drops of reagent A and 2 drops of reagent B
- DAB stain: 5 ml dH₂O with: 2 drops of buffer (pH 7.5), 4 drops of ABC and 2 drops of hydrogen peroxide
- Hematoxylin
- 1% Acid alcohol: 350ml industrial methylated spirit (denatured alcohol), 5ml HCl, 145ml dH₂O
- Mounting medium: DPX

CONSUMABLES

- Slides
- Tip box for hematoxylin

EQUIPMENT

- Microwave

Safety warnings

 For safety warnings and hazards, see the Safety Data Sheet (SDS).



Day 1

5m

1 Sequential rehydration

Xylene, 100% EtOH, 90% EtOH, 70% EtOH, dH₂O – 5 min in each:

1.1 Rehydrate samples in Xylene for 00:05:00 .

5m

1.2 Rehydrate samples in 100% EtOH for 00:05:00 .

5m

1.3 Rehydrate samples in 90% EtOH for 00:05:00 .

5m

1.4 Rehydrate samples in 70% EtOH for 00:05:00 .

5m

1.5 Rehydrate samples in dH₂O for 00:05:00 .

5m

2 Antigen Retrieval

20m

Microwave slides at high power for 00:20:00 in

[M] 0.05 % (v/v) Tween 20 in Sodium Citrate Buffer , pH 6 .

Note

Sodium citrate recipe:

🧪 2.94 g Tri-sodium citrate (dehydrate)

🧪 1000 mL Distilled water

Mix to dissolve, adjust to pH 6 with 1N HCl and then add 🧪 0.5 mL Tween 20 and mix well.

*Store at RT for 3 months or at 4°C for longer storage.

3 Wash



3 × 5 min with TBPS ([M] 0.05 % (v/v) Tween 20 in 1xPBS):



3.1 Wash 1 of 3: ⌚ 00:05:00 with TBPS (1x PBS + 0.05% Tween 20)

5m

3.2 Wash 2 of 3: ⌚ 00:05:00 with TBPS (1x PBS + 0.05% Tween 20)

5m

3.3 Wash 3 of 3: ⌚ 00:05:00 with TBPS (1x PBS + 0.05% Tween 20)

5m

4 Endogenous peroxidase blocking

15m

Incubate slides in [M] 3 % (v/v) hydrogen peroxide (H₂O₂) in TPBS at

🌡 Room temperature for ⌚ 00:15:00 .



5 Wash

3 × 5 min with TBPS:



5.1 Wash 1 of 3: ⌚ 00:05:00 with TBPS

5m

5.2 Wash 2 of 3: ⌚ 00:05:00 with TBPS

5m

5.3 Wash 3 of 3: ⌚ 00:05:00 with TBPS

5m

6 Blocking for 1 hour

Blocking Solution (e.g. 5% Normal horse serum in TPBS).

7 Primary antibody – Overnight incubation

15m

Dilute in TPBS in appropriate concentration and incubate ⌚ Overnight at 🌡 4 °C .



Day 2

4h 3m 31s

8 Wash

3 × 5 min with TPBS:



8.1 Wash 1 of 3: ⌚ 00:05:00 with TBPS

5m



- 8.2 Wash 2 of 3: ⌚ 00:05:00 with TBPS 5m
- 8.3 Wash 3 of 3: ⌚ 00:05:00 with TBPS 5m
- 9 **Secondary antibody** 2h
Dilute antibody in blocking solution and incubate for ⌚ 02:00:00 at
🌡 Room temperature
- 10 **Wash**
3 × 5 min with 1xTPBS:
- 10.1 Wash 1 of 3: ⌚ 00:05:00 with 1xTBPS 5m
- 10.2 Wash 2 of 3: ⌚ 00:05:00 with 1xTBPS 5m
- 10.3 Wash 3 of 3: ⌚ 00:05:00 with 1xTBPS 5m
- 11 **ABC solution (DAB substrate)** 1h
Note
*Prepare ⌚ 00:30:00 before use!
Make 🧪 5 mL TPBS with 2 drops of reagent A and 2 drops of reagent B. Leave ABC on tissue for ⌚ 01:00:00
- 12 **Wash**
3 × 5 min with TPBS:
- 12.1 Wash 1 of 3: ⌚ 00:05:00 with TBPS 5m



12.2 Wash 2 of 3:  00:05:00 with TBPS

5m

12.3 Wash 3 of 3:  00:05:00 with TBPS

5m

13 DAB Stain

Make  5 mL dH₂O with: 2 drops of buffer ( 7.5), 4 drops of ABC and 2 drops of hydrogen peroxide. Incubate with DAB at  Room temperature until colour change.



14 Wash

2 × 5 min with dH₂O:



14.1 Wash 1 of 2:  00:05:00 with dH₂O

5m

14.2 Wash 2 of 2:  00:05:00 with dH₂O

5m

15 Counterstaining with Hematoxylin

Place slides in tip box of hematoxylin for  00:00:30 .

30s

15.1 1 x  00:05:00 wash with dH₂O.

5m



15.2 Put slides in acid alcohol,  00:00:01 quickly in and out. Keep the slides in running water for  00:03:00 for blueing.

3m 1s

Note

1% Acid alcohol:

 350 mL Industrial methylated spirit (denatured alcohol)

 5 mL HCl

 154 mL dH₂O

16 Dehydration process



70% EtOH, 90% EtOH, 100% EtOH, Xylene:

16.1 Dehydrate samples in 70% EtOH.

16.2 Dehydrate samples in 90% EtOH.

16.3 Dehydrate samples in 100% EtOH.

16.4 Dehydrate samples in Xylene.

17 **Mounting Slides**

Mount slides using medium DPX.