


Mar 18, 2024

Immunohistochemistry protocol optimized at iMM-JLA

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlkon2wv5r/v1>



Ana M Biscaia Santos¹, Ana Cristóvão Pinto¹, Ana R Pires¹, Joana Antunes¹

¹Comparative Pathology Unit, Instituto de Medicina Molecular - João Lobo Antunes, Lisboa, Portugal

Ana M Biscaia Santos: Biomedical Scientist Specialist

Ana Cristóvão Pinto: Biomedical Scientist

Ana R Pires: Biomedical Scientist Specialist

Joana Antunes: Biomedical Scientist - corresponding author: joana.antunes@medicina.ulisboa.pt

iMM Comparative Pathol...



Joana G Antunes

Gulbenkian Institute for Molecular Medicine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlkon2wv5r/v1>



Protocol Citation: Ana M Biscaia Santos, Ana Cristóvão Pinto, Ana R Pires, Joana Antunes 2024. Immunohistochemistry protocol optimized at IMM-JLA. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.j8nlkon2wv5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 31, 2023

Last Modified: March 18, 2024

Protocol Integer ID: 90197

Keywords: immunohistochemistry, histology, immunostaining, FFPE, frozen, Envision HRP, immunohistochemistry protocol, jla the immunohistochemistry, immunohistochemistry, specific antigen in the protein, antibody, specific antigen, primary antibody, secondary antibody, concentration of the antibody, ph of the antigen retrieval solution, antigen retrieval solution, antigen retrieval, location of specific protein, type of antigen retrieval, location of the protein, specific protein, hydrogen peroxide in the substract, protein, heat induced epitope retrieval, horseradish peroxidase, hydrogen peroxide, protein of interest, jla, dab solution, ihc

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.



Abstract

The immunohistochemistry (IHC) protocol is used to showcase the location of specific proteins in relation to the architecture of the tissue.

In this method, we use a primary antibody which binds to a specific antigen in the protein of interest. This antibody has a specific host and can be ready to use or undiluted. The secondary antibody is an anti-host IgG which will bind to the primary antibody ¹. It's also conjugated with an Horseradish Peroxidase (HRP).

The HRP will react with the Hydrogen Peroxide in the substrate of the DAB solution, oxidizing the DAB, creating a brown deposit on the location of the protein of interest.

Depending on location of the protein, we can have staining in the nuclei, the cytoplasm or the cellular membrane.

Possible variables are:

the concentration of the antibody;

type of antigen retrieval;

the pH of the antigen retrieval solution (for Heat Induced Epitope Retrieval (HIER));

the time of incubation;

the tissue used ².

You can use this protocol for Formalin-Fixed Paraffin-Embedded samples (FFPE) and Fixed Frozen samples (Gelatin and OCT).

Materials

Coplin Jars
Incubation chamber
Slide Dipper
Slide dishes with Lids
Wash bottles with nozzle
Magnets
Magnetic mixer
Tweezers
Micropipettes p20, p200 and p1000 and respective tips
Hydrophobic PAP pen
Eppendorfs 1.5mL
Plastic pipettes
Equipment Pretreatment ModuleTM Deparaffinization and Heat-Induced Epitope Retrieval by Thermo Scientific

Reagents:

Ethanol 70%
Ethanol 96%
Ethanol 100%
Xylene
Distilled Water

Protocol materials

- ✕ EnvisionTM Flex Wash Buffer **Dako Catalog #Ref. K8007**
- ✕ EnVisionTM Flex conjugated w/ HRP **Dako Catalog #Ref. K4061**
- ✕ Protein Block Serum free **Dako Catalog #Ref. X0909**
- ✕ Hydrogen Peroxide Solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #7722-84-1**
- ✕ EMSURE Methanol **Supelco Catalog #1.06009.2500**
- ✕ Entellan ® **Merck MilliporeSigma (Sigma-Aldrich) Catalog #1.00869**
- ✕ EnvisionTM Dako DAB **Dako Catalog #Ref. K3468**
- ✕ Harris hematoxylin for histology **Bio-optica Catalog #05-06004E**
- ✕ Antibody Diluent **Dako Catalog #Ref. 8006**
- ✕ EnVision FLEX Target Retrieval Solution Low pH **Dako Catalog #K8005**
- ✕ EnVision FLEX Target Retrieval Solution High pH **Dako Catalog #K8004**



Safety warnings

! This protocol uses the following hazardous chemicals:

Xylene

Methanol

Hydrogen peroxide

DAB

Resinous Mounting Medium

Please comply with all safety measures and risk assessment protocols.

Latex gloves can be used for the entire protocol, except when manipulating the reagents previously stated, during which you should use nitrile gloves.

Type IV biological disposable bin should be readily available.

Before start

If FFPE - begin protocol after adhering sections to the slides for 1h at 60°C



Preparation of Wash Buffer

- 1 Measure  20 mL of 20X
 Envision™ Flex Wash Buffer **Dako Catalog #Ref. K8007**
- 2 Make up with distilled water to  1 L of total volume

Sample Preparation after sectioning - FFPE

- 3 Deparaffinize slides with Xylene for  00:10:00 10m
- 4 Change Xylene and deparaffinize for another  00:10:00 10m
- 5 Place slides in Ethanol 100% for  00:05:00 5m
- 6 Place slides in Ethanol 96% for  00:05:00 5m
- 7 Place slides in Ethanol 70% for  00:05:00 5m
- 8 To complete the hydration, place slides in distilled Water for  00:05:00 5m

Sample Preparation after sectioning - Frozen

- 9 **If frozen samples:** 10m
Remove from storage and let dry.
If embedded in OCT, put in distilled water for  00:05:00 at  Room temperature .
If embedded in gelatin, put in PBS for  00:05:00 at  37 °C .

Antigen Retrieval



- 10 Perform antigen retrieval using PT Module by Thermo Scientific with a cycle of pre-heat to 65°C and heat to 95°C for 20 minutes.

Note

This equipment, when optimized, allows you to skip the step of deparaffinization in this section, by putting the slides directly through the antigen retrieval step.

This equipment allows two buffers with different pH, in which we use

⊗ EnVision FLEX Target Retrieval Solution High pH **Dako Catalog #K8004** and

⊗ EnVision FLEX Target Retrieval Solution Low pH **Dako Catalog #K8005** ,

depending on the optimization of the antibody.

Blocking Endogenous Peroxidase and Total Proteins

35m

- 11 Wash in 1X ⊗ Envision™ Flex Wash Buffer **Dako Catalog #Ref. K8007** 3 times for
⌚ 00:05:00 each

5m



- 12 Put slides in 🧪 200 mL of 100%

⊗ EMSURE Methanol **Supelco Catalog #1.06009.2500**

Note

Extemporaneous solution. 100% methanol can be used up to 3 times.

- 13 Put 🧪 6 mL of

⊗ Hydrogen Peroxide Solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #7722-84-1**

at 🌡 4 °C in the 100% ⊗ EMSURE Methanol **Supelco Catalog #1.06009.2500** with the slides

- 14 Cover the slide dish and incubate slides in the solution for ⌚ 00:30:00

30m

Note

Hydrogen Peroxide must be added each time when immersing the slides, to obtain a concentration of 3%.

- 15 Wash in 1X  EnvisionTM Flex Wash Buffer **Dako Catalog #Ref. K8007** 3 times for  00:05:00 each



- 16 Draw a limit on the edge of the tissue with the hydrophobic PAP Pen

- 17 Put  Protein Block Serum free **Dako Catalog #Ref. X0909** on the tissue

Incubation and Washing

1h 54m 10s

- 18 Incubate with Primary antibody in the optimized concentration (if needed, dilute in  Antibody Diluent **Dako Catalog #Ref. 8006**) for  01:00:00 at  Room temperature

1h



Note

Normally, around 100µL of diluted antibody per section, depending on the size of the tissue.

- 19 Wash in 1X  EnvisionTM Flex Wash Buffer **Dako Catalog #Ref. K8007** 3 times for  00:05:00 each



- 20 Incubate with secondary antibody  EnVisionTM Flex conjugated w/ HRP **Dako Catalog #Ref. K4061** for  00:30:00  Room temperature

30m





Note

Be sure to match the secondary antibody to the host of the primary antibody (ex: if using mouse anti-CD3 primary, use anti-mouse secondary)

- 21 Wash in 1X  EnvisionTM Flex Wash Buffer **Dako Catalog #Ref. K8007** 3 times for  00:05:00 each



- 22 Incubate with  EnvisionTM Dako DAB **Dako Catalog #Ref. K3468** for  00:02:00 to  00:05:00

7m



Safety information

After incubation, pour DAB excess on a paper towel and dispose on type IV waste, with tips and eppendorfs used, avoiding spills at any cost.
Clean the incubation chamber and surfaces with Bleach.

- 23 Wash in distilled water for  00:05:00

5m



Counterstain and mount

17m 10s

- 24 Counterstain with  Harris hematoxylin for histology **Bio-optica Catalog #05-06004E** for  00:00:10

10s

Note

filter Hematoxylin before use

- 25 Put slides in lukewarm running water for  00:05:00

5m



- 26 Dehydrate with 96% ethanol for 00:02:00 2m
- 27 Place in 100% ethanol for 00:02:00 2m
- 28 Clear in Xylene for 00:10:00 10m
- 29 Mount with appropriate mounting medium (i.e.
 Entellan[®] Merck MilliporeSigma (Sigma-Aldrich) Catalog #1.00869)

Note

Let the slides dry in a fume hood until complete polymerization of the mounting medium is achieved

Protocol references

- 1 - Magaki S, Hojat SA, Wei B, So A, Yong WH. An Introduction to the Performance of Immunohistochemistry. Methods Mol Biol. 2019;1897:289-298. doi: 10.1007/978-1-4939-8935-5_25. PMID: 30539453; PMCID: PMC6749998.
- 2 - Kim SW, Roh J, Park CS. Immunohistochemistry for Pathologists: Protocols, Pitfalls, and Tips. J Pathol Transl Med. 2016 Nov;50(6):411-418. doi: 10.4132/jptm.2016.08.08. Epub 2016 Oct 13. PMID: 27809448; PMCID: PMC5122731.