

Aug 01, 2023

Immunohistochemistry (using floating section)

DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvr36yovmk/v1>

Tae-Un Han¹, Ellen Sidransky¹

¹National Institute of Health



Tae-Un Han

NIH/NHGRI

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr36yovmk/v1>

Protocol Citation: Tae-Un Han, Ellen Sidransky 2023. Immunohistochemistry (using floating section). protocols.io <https://dx.doi.org/10.17504/protocols.io.6qpvr36yovmk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 28, 2023



Last Modified: May 31, 2024

Protocol Integer ID: 85658

Keywords: ASAPCRN, phospho s129 alpha synuclein protein expression, immunohistochemistry, immunohistochemistry protocol, floating section, mouse brain, pathology, slide, evaluation of pathology, protein

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000458

Abstract

- This is an immunohistochemistry protocol used for evaluation of pathology in mouse brain.
- We successfully detected Gfap, Iba1, CD68, TH and alpha-Syn and phospho S129 alpha synuclein protein expression using this protocol.
- In this protocol, floating sections are mounted on slides first before starting blocking step. In this way, large number of sections can be stained more efficiently.

Materials

30 um thickness of floating tissue sections in PBS with 0.1% sodium azide.

- Washing: Glass Coplin jar

- PBST: PBS 1X + 0.05% Tween20

- Normal Goat Serum Blocking Solution:

2% goat serum (blocking)

1%BSA (stabilizer)

0.1% cold fish skin gelatin (blocking)

0.1% Triton X-100 (penetration enhancer)

0.05% Tween 20 (detergent and surface tension reducer)

0.05% sodium azide (preservative)

0.01M PBS, pH 7.2

Mix well and store at 4 °C for up to 6 months.

- Primary Antibody Dilution Buffer

1%BSA (blocking & stabilizer)

0.1% cold fish skin gelatin (blocking)

0.5% Triton X-100 (penetration enhancer)

0.05% sodium azide (preservative)

0.01M PBS, pH 7.2-7.4

Mix well and store at 4 °C. up to 6 months.

- Secondary Antibody Dilution Buffer:

0.01M PBS, pH 7.2

0.05% Tween 20

- List of antibodies used for neuropathology evaluation with IHC

	A	B	C	D
	Antibodies	Cat#	Company	dilution
	Anti-Synuclein α Antibody	AB5038	Sigma- Aldrich	1:1000
	Recombinant Anti-Alpha-synuclein (phospho S129) antibody [EP1536Y]	ab51253	Abcam	1:1000



	A	B	C	D
	GFAP Antibody , Rabbit Polyclonal	16825-1- AP	Proteinte ch	1:4000
	Recombinant Anti-Iba1 antibody [EPR16588]	ab178846	Abcam	1:1000
	Anti-Tyrosine Hydroxylase	AB152	Sigma- Aldrich	1:1000
	Anti-CD68 antibody	ab125212	Abcam	1:1000
	Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) : 2nd antibody	ab150077	Abcam	1:1000



- 1 Mount floating sections onto slides and dry for 00:30:00 (<45 min), While section drying, PAP pen drawing line and dry 00:05:00 . 35m
- 2 PBS rehydrate 00:10:00 , Wash with PBS for 00:05:00 three times. 15m
- 3 Blocking sol 500 μ L , incubate 01:00:00 at Room temperature . 1h
- 4 Replace buffer to primary antibody sol, incubate Overnight at 4 °C . 10m
- 5 Wash in PBST three times for 00:05:00 , 00:10:00 and 00:05:00 . 20m
- 6 Secondary antibody sol, incubate for 01:30:00 at Room temperature . 1h 30m
- 7 Wash in PBS three times for 00:05:00 , 00:10:00 and 00:05:00 . 20m
- 8 Mounting with Vectorshield mounting media