



Apr 02, 2023

## Immunohistochemistry using paraffin embedded tissue

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2ly72ywlx9/v1>

Michael J Hurley<sup>1,2</sup>

<sup>1</sup>Department of Clinical and Movement Neurosciences, UCL Queen Square Institute of Neurology, London, UK;

<sup>2</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, Maryland, USA



Michael J Hurley

UCL

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.eq2ly72ywlx9/v1>

**Protocol Citation:** Michael J Hurley 2023. Immunohistochemistry using paraffin embedded tissue. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.eq2ly72ywlx9/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** February 27, 2023

**Last Modified:** May 31, 2024

**Protocol Integer ID:** 77671

**Keywords:** ASAPCRN, Immunohistochemistry, paraffin wax, immunofluorescence, synuclein by immunofluorescence, immunofluorescence, immunohistochemistry, antibody, peroxidase immunohistochemistry, synuclein, tissue this protocol, optimal dilution of the primary antibody, primary antibody, embedded tissue, other antibody, using paraffin, most antibody, many other antibody, only change necessary for most antibody, antigen, paraffin wax, best type of antigen, peroxidase, paraffin, sections of tissue, tissue

**Funders Acknowledgements:**

Aligning Science Across Parkinson's

Grant ID: ASAP-000420

## Abstract

This protocol describes how to stain paraffin wax embedded sections of tissue for alpha-synuclein by immunofluorescence or DAB/peroxidase immunohistochemistry and is based on Hurley et al., 2013 (<https://doi.org/10.1093/brain/awt134>).

The protocol works for many other antibodies and the only change necessary for most antibodies is determining the optimal dilution of the primary antibody and the best type of antigen-retrieval to use.

## Materials

Rabbit anti-alpha-synuclein ab212184 Abcam

Goat anti-rabbit IgG biotin111-065-144-JIR Stratech

Goat anti-rabbit IgG Alexa Fluor™ 488 A27034 Thermo Fisher

Streptavidin POD11089153001



## Wax embedding

- 1 Dissect tissue (fresh or trans-cardially perfused with formalin (10% neutral buffered formaldehyde) or 4% paraformaldehyde in PBS).

### Safety information

Formaldehyde is toxic by inhalation  
[HT501128 \(sigmaaldrich.com\)](#)

- 2 Place tissue in at least 5 volumes of fixative at 4 °C for > 48:00:00 if tissue was perfused or > 96:00:00 if fresh. The tissue can stay in fixative indefinitely. If the tissue was fresh and the solution is turbid with blood, change the fixative after 24 hours.
- 3 Embed formalin fixed tissue in wax using an automated tissue processor. It can be done manually if you do not have access to an automated tissue processor. Histo-Clear™ II can be used instead of xylene.

6d

### Safety information

Xylene is toxic by inhalation and flammable  
[Xylene - Safety Data Sheet \(chemicalbook.com\)](#)

Ethanol is toxic inhalation and flammable  
[Ethanol - Safety Data Sheet \(chemicalbook.com\)](#)

### 3.1 Gut tissue

13h 42m

- |    |               |          |                  |
|----|---------------|----------|------------------|
| 1. | 10 % formalin | 00:01:00 | Room temperature |
| 2. | 10 % formalin | 00:01:00 | Room temperature |
| 3. | 70 % ethanol  | 01:00:00 | 30 °C            |
| 4. | 90 % ethanol  | 01:00:00 | 30 °C            |
| 5. | 100 % ethanol | 01:00:00 | 30 °C            |
| 6. | 100 % ethanol | 01:00:00 | 30 °C            |
| 7. | 100 % ethanol | 01:00:00 | 30 °C            |



|     |               |          |       |
|-----|---------------|----------|-------|
| 8.  | 100 % ethanol | 01:20:00 | 30 °C |
| 9.  | xylene        | 01:00:00 | 30 °C |
| 10. | xylene        | 01:00:00 | 30 °C |
| 11. | xylene        | 01:20:00 | 30 °C |
| 12. | paraffin wax  | 01:20:00 | 62 °C |
| 13. | paraffin wax  | 01:20:00 | 62 °C |
| 14. | paraffin wax  | 01:20:00 | 62 °C |

### 3.2 Brain tissue (whole mouse)

1d 2h 2m

|     |               |          |                  |
|-----|---------------|----------|------------------|
| 1.  | 10 % formalin | 00:01:00 | Room temperature |
| 2.  | 10% formalin  | 00:01:00 | Room temperature |
| 3.  | 70% ethanol   | 02:00:00 | 30 °C            |
| 4.  | 90% ethanol   | 02:00:00 | 30 °C            |
| 5.  | 100% ethanol  | 02:00:00 | 30 °C            |
| 6.  | 100% ethanol  | 02:00:00 | 30 °C            |
| 7.  | 100% ethanol  | 02:00:00 | 30 °C            |
| 8.  | 100% ethanol  | 02:00:00 | 30 °C            |
| 9.  | xylene        | 02:00:00 | 30 °C            |
| 10. | xylene        | 02:00:00 | 30 °C            |
| 11. | xylene        | 04:00:00 | 30 °C            |
| 12. | paraffin wax  | 02:00:00 | 62 °C            |
| 13. | paraffin wax  | 02:00:00 | 62 °C            |
| 14. | paraffin wax  | 02:00:00 | 62 °C            |

4 Mount the tissue in a block of wax using an embedding station and metal or plastic moulds.

Store at Room temperature

5 Section the tissue ( 5 µm to 8 µm ) using a microtome.

Store sections at Room temperature



## Dewaxing and antigen retrieval

- 6 Select slides to be stained and place in a microscope slide rack.
- 7 Heat sections  60 °C  00:30:00 30m
- 8 Dewax sections in two changes of Histo-Clear II (National Diagnostics). 20m
1.  00:10:00
2.  00:10:00
- 9 Rehydrate sections through graded (v/v) ethanol solutions. 25m
1. 100 %  00:05:00
2. 100 %  00:05:00
3. 95 %  00:05:00
4. 80 %  00:05:00
5. 70 %  00:05:00
- 10 Wash with PBS  00:05:00 5m
- 11 Quench endogenous peroxidase with 0.3 % hydrogen peroxide (v/v)  00:05:00 5m
- 12 Wash with PBS  00:05:00 5m
- 13 For antigen retrieval heat  100 °C the sections in 15 mM sodium citrate buffer containing 0.1% Tween® 20, pH 6 for  00:25:00 in an 800-watt microwave set on 30 % power. 25m
- 14 Cool with running tap water  00:05:00 5m
- 15 Wash with PBS 2 x  00:05:00 5m



16 For **immunofluorescence** staining go to step 17

For **DAB/peroxidase** staining go to step 27

## Immunofluorescence staining

23h 25m

17 Draw around the sections with a wax pen.

18 Block non-specific binding sites with 10 % goat serum in PBS containing 0.005 % Triton™ X-100 and 0.05 % thimerosal ⌚ 01:30:00

Depending on the size and number of sections 150-250 microliters of solution should be sufficient to cover the sections.

19 Incubate with anti-alpha-synuclein antibody (1:500) (ab212184) in 10 % goat serum containing 0.005 % Triton™ X-100 and 0.05 % thimerosal (antibody buffer)

20h

⌚ 20:00:00

🌡 Room temperature

20 Wash with PBS 4 x ⌚ 00:05:00

5m

21 incubated sections with Alexa Fluor goat anti-rabbit IgG (1:200) (A27034) ⌚ 01:30:00

1h 30m

22 Wash with PBS 2 x ⌚ 00:05:00

5m

23 incubate sections with DAPI (1 µg/ml) ⌚ 00:05:00

5m

24 Wash with PBS 3 x ⌚ 00:05:00

5m

25 Dip slides in water to remove salts

26 Mount sections using Hydromount™ (National Diagnostics).



## DAB/peroxidase staining

43m 15s

- 27 Draw around the sections with a wax pen.
- 28 Block non-specific binding sites with 10 % goat serum in PBS containing 0.005 % Triton™ X-100 and 0.05 % thimerosal  01:30:00 1h 30m  
Depending on the size and number of sections 150-250 microliters of solution should be sufficient to cover the sections.
- 29 Incubate with anti-alpha-synuclein antibody (1:500) (ab212184) in 10 % goat serum containing 0.005 % Triton™ X-100 and 0.05 % thimerosal (antibody buffer) 20h  
 20:00:00  Room temperature
- 30 Wash PBS 4 x  00:05:00 5m
- 31 Incubate sections in biotinylated goat anti-rabbit IgG (1:250) secondary antibody for 2 hours.
- 32 Wash PBS 4 x  00:05:00 5m
- 33 Incubate sections with streptavidin-peroxidase conjugate (1:250) for 1 hour.
- 34 Wash PBS 4 x  00:05:00 5m
- 35 Visualise staining by incubation with PBS containing 0.05 % (w/v) 3,3'-Diaminobenzidine/0.015 % (v/v) H<sub>2</sub>O<sub>2</sub>/0.05 % (w/v) nickel ammonium sulphate for up to  00:03:00 , 3m
- 36 Wash sections for >  00:05:00 in running tap water. 5m
- 37 Counter stain with haematoxylin  00:00:05 to  00:00:10 if desired. 15s
- 38 Wash sections with running tap water until it runs clear to remove excess stain.



39 Dehydrate sections through graded (v/v) ethanol solutions

25m

70 %  00:05:00

80 %  00:05:00

95 %  00:05:00

100 %  00:05:00

100 %  00:05:00

40 Clear sections with Histo-Clear II (National Diagnostics) **2 x**  00:05:00

5m

41 Mount with Omnimount™ (National Diagnostics) or equivalent permanent resin.

## Protocol references

<https://doi.org/10.1093/brain/awt134>