


Jun 06, 2023

Version 1

## Immunolabeling of a Whole Mouse Brain V.1

DOI

[dx.doi.org/10.17504/protocols.io.ewov1okwylr2/v1](https://dx.doi.org/10.17504/protocols.io.ewov1okwylr2/v1)



Andrew Recknagel<sup>1</sup>, Kevin Cao<sup>2</sup>, Judith Baka<sup>2</sup>, Naveen Ouellette<sup>2</sup>, Jayaram Chandrashekar<sup>2</sup>, Molly Logsdon<sup>2</sup>

<sup>1</sup>Janelia Farm; <sup>2</sup>Allen Institute for Neural Dynamics

Allen Institute for Neural...



**Naveen Ouellette**

Allen Institute

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1okwylr2/v1>

**Protocol Citation:** Andrew Recknagel, Kevin Cao, Judith Baka, Naveen Ouellette, Jayaram Chandrashekar, Molly Logsdon 2023. Immunolabeling of a Whole Mouse Brain. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.ewov1okwylr2/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** March 24, 2023

**Last Modified:** June 18, 2024

**Protocol Integer ID:** 79369

**Keywords:** immunolabeling, whole brain, antibody, delipidated brain sample, fluorescent protein, viral labeling of target neuron population, encoded fluorescent protein, whole mouse brain, whole brain sample, target neuron population, antibody labeling, accessible to antibody labeling, mice, antibody, brain, immunolabeling, fluorophore, concentration of fluorophore, organic solvent delipidation step, delipidation, signal from endogenous fp


**Funders Acknowledgements:**

Allen Institute

## Abstract


Our mice have sparse viral labeling of target neuron populations with genetically encoded fluorescent proteins (FPs). Delipidation ensures that whole brain samples become optically clear and more accessible to antibody labeling. However, the signal from endogenous FPs is often dimmed or quenched during organic solvent delipidation steps. Tissue expansion also dilutes the concentration of fluorophores. To amplify signals of FPs for high resolution, multi-scale imaging, we immunolabel delipidated brain samples.

## Materials

 PBS - Phosphate-Buffered Saline (10X) pH 7.4 **Invitrogen - Thermo Fisher Catalog #AM9625**

 Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8787-50ML**

 Tween 20 **Catalog #P1379**

 5% Sodium Azide **Fisher Scientific Catalog #71448-16**

Equipment

VWR® Standard Orbital Shaker, Model 3500	NAME
Shaker	TYPE
VWR	BRAND
89032-092	SKU
<a href="https://us.vwr.com/store/product/4901055/vwr-standard-orbital-shaker-model-3500">https://us.vwr.com/store/product/4901055/vwr-standard-orbital-shaker-model-3500</a>	LINK



## Equipment

**WHEATON® Liquid Scintillation Vials, Caps Attached to Vials, Glass, Polyethylene Cone, 22-400, 20 mL**

NAME

Vial

TYPE

Wheaton

BRAND

DWK986546

SKU

<https://www.dwk.com/na/wheaton-liquid-scintillation-vials-caps-attached-to-vials-glass-polyethylene-cone-22-400-20-ml-986546>

LINK

20 mL Glass Vial with Polyethylene cone Caps

SPECIFICATIONS



## Equipment

**WHEATON® Shorty Vials clear with PTFE faced rubber lined cap**

NAME

vial

TYPE

WHEATON

BRAND

DWKW224607

SKU

<https://www.sigmaaldrich.com/US/en/product/aldrich/dwkw224607>

LINK



## Equipment

### Nutating Mixer

NAME

Mixer

TYPE

Fisherbrand

BRAND

88-861-043

SKU

<https://www.fishersci.com/shop/products/nutating-mixers-variable-speed/88861043><sup>LINK</sup>

16.3 × 11.5 × 10.7 in.(415 × 293 × 273 mm)

SPECIFICATIONS



## RECIPES

### PTxw: 0.05% Tween 20, 0.1% Triton X-100, 0.04% NaN<sub>3</sub> in PBS

Combine the following reagents and stir at Room temperature until fully dissolved, store at 4 °C .

	Reagent	Amount
	Milli-Q Water	up to 500 mL
	10X PBS	50 mL
	Triton X-100	500 uL
	Tween 20	250 uL
	5% Sodium azide	4 mL
	<b>Total</b>	<b>500 mL</b>



## Troubleshooting

### Before start

If performing the **Whole Mouse Brain Delipidation, Immunolabeling, and Expansion Microscopy** protocol, start with a brain that has been delipidated using the **Tetrahydrofuran and Dichloromethane Delipidation of a Whole Mouse Brain** and **Aqueous (SBiP) Delipidation of a Whole Mouse Brain** protocols.



## Immunolabeling of a Whole Mouse Brain

4w 3d

1 Start with a whole mouse brain perfused with 4% PFA, post-fixed, and stored in PBS in a 20 mL vial.

2 Using a 20 mL glass vial, wash the brain with PTxw, rotating at Room temperature Overnight

3 Transfer the brain to a small 4 mL vial containing diluted primary antibody in ~ 4 mL PTxw and incubate with rocking, a nutator or shaker, at Room temperature for 11 days.

- 10 µg/brain for primary antibody
- 20 µg/brain for primary conjugated antibody

1w 4d

4 Transfer the brain to a 20 mL glass vial and wash with PTxw, rocking on a nutator or shaker at Room temperature , and replacing solution for each step:

3d 11h

- PTxw for 01:00:00
- PTxw for 02:00:00
- PTxw for 04:00:00
- PTxw Overnight
- PTxw for 24:00:00
- PTxw for 24:00:00
- PTxw for 24:00:00 +

5 Transfer the brain to a small 4 mL vial containing diluted secondary antibody in ~ 4 mL PTxw and incubate with rocking at Room temperature for 11 days.

- 20 µg/brain for secondary antibody

6 Transfer brain to a 20 mL glass vial and wash with PTxw, rocking at Room temperature , and replacing solution for each step:

3d 11h

- PTxw for 01:00:00
- PTxw for 02:00:00
- PTxw for 04:00:00



- PTxw Overnight
- PTxw for 24:00:00
- PTxw for 24:00:00
- PTxw for 24:00:00 +

7 Wash the brain with 1X PBS, rocking at Room temperature Overnight .

16h