

May 24, 2023

Immunoprecipitation (IP)

DOI

dx.doi.org/10.17504/protocols.io.eq2ly79yelix9/v1

nguyen.tha¹

¹Laboratory of Michael Lazarou, Walter and Eliza Hall Institute of Medical Research, Parkville, Victoria, Australia

ASAP workspace



nguyen.tha

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2ly79yelix9/v1>

Protocol Citation: nguyen.tha 2023. Immunoprecipitation (IP). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.eq2ly79yelix9/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 07, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 76523

Keywords: immunoprecipitation, ASAPCRN, protocol details about immunoprecipitation, immunoprecipitation, protocol detail, protocol, ip

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000350

Abstract

This protocol details about immunoprecipitation using anti-HA magnetic beads.

Attachments








[635-1314.docx](#)

17KB





Materials

Buffers and reagents:

- **IP base buffer:** [M] 50 millimolar (mM) Tris-Cl ( 7.5 when cold), [M] 150 millimolar (mM) NaCl
- **Bead equilibration buffer:** IP base buffer supplemented with 0.1% Tween20
- **IP wash buffer:** IP base buffer supplemented with 0.1% TX-100 and 1x cOmplete, EDTA-free protease inhibitor cocktail
-  Pierce & Warriner Anti-HA Magnetic Beads **Thermo Fisher Catalog #88836**
-  Benzonase® Nuclease **Merck Millipore (EMD Millipore) Catalog #E1014-25KU**
-  cOmplete™ EDTA-free Protease Inhibitor Cocktail **Roche Catalog #4693132001**
- **Elution buffer:**  NuPAGE™ LDS Sample Buffer (4X) **Thermo Fisher Scientific Catalog #NP0007**

Note






















Elution buffer can be aliquoted and stored at  -20 °C or  -80 °C .

Troubleshooting



Procedures

3h 50m

- 1 Lyse cell pellets ( 5-7 mg) in  500 μ L IP lysis buffer containing IP base buffer supplemented with 1x cOmplete, EDTA-free protease inhibitor cocktail and  0.1 μ L of    benzonase and incubate samples  On ice for  00:30:00 . Mix the sample by inverting the eppies gently every 5 min. 30m
- 2 Wash anti-HA beads with  500 μ L of bead equilibration buffer. 
- 3 Repeat step 2 twice.
- 4 Centrifuge the cell lysates at max speed for  00:10:00 at  4 $^{\circ}$ C . 10m 
- 5 Carefully transfer cleared lysates into 2 ml eppies and take  50 μ L from each tube for "Input" samples.
- 6 Gently add  1000 μ L of IP base buffer containing 1x cOmplete, EDTA-free protease inhibitor cocktail to the rest of each sample to dilute out the detergent. 
- 7 Incubated the diluted cleared lysates with the anti-HA magenetic beads on a rotary mixer for  03:00:00 at  4 $^{\circ}$ C . 3h 
- 8 Collect beads on a magnetic rack and aspirate the unbounds.
- 9 Wash with  1 mL IP wash buffer. 
- 10 Repeat steps 7-8 another 4 times.

Note

For the last wash, make sure to remove all the liquid off the beads.



11

Elute with  25 μ L elution buffer by boiling at shaking at  99 °C for  00:10:00

10m