



Sep 23, 2023

Immunoprecipitation of NAP1-GFP

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl8pmndg2w/v1>

Elias Adriaenssens¹

¹Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8pmndg2w/v1>

Protocol Citation: Elias Adriaenssens 2023. Immunoprecipitation of NAP1-GFP. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.5jyl8pmndg2w/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 27, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 84128

Keywords: Immunoprecipitation, ASAPCRN, immunoprecipitation of nap1, gfp from hap1 cell, hap1 cell, immunoprecipitation, gfp this protocol, nap1, gfp, cell

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-000350

Marie Skłodowska-Curie MSCA Postdoctoral fellowship

Grant ID: 101062916

Abstract

This protocol describes immunoprecipitation of NAP1-GFP from HAP1 cells.

Attachments



[757-1925.pdf](#)

43KB



Materials

Materials

- NAP1-EGFP (Addgene xxx)
- Lipofectamine 2000 (Thermo Fisher)
- SDS-PAGE gels (NP0322BOX, Thermo Fisher)
- PageRuler Prestained protein marker (Thermo Fisher)

Lysis buffer

	A	B
	KCl	100 mM
	MgCl ₂	2.5 mM
	Tris-HCl pH 7.4	20 mM
	NP-40	0.50%

Fixation solution

	A
	40% ethanol
	10% acetic acid
	50% dH ₂ O

 Pierce & Warriner; Detergent Compatible Bradford Assay Kit **Thermo Fisher Catalog #23246**

 ChromoTek GFP-Trap® Agarose **Proteintech Catalog # gta-20**



Immunoprecipitation

1h

- 1 Transiently transfect FIP200 knockout HAP1 cells with pcDNA3.1 NAP1-EGFP (Addgene), or empty pcDNA3.1 vector as a negative control, using Lipofectamine 2000 (Thermo Fisher).

Note

This cell line was selected as FIP200 deletion results in TBK1 hyperactivation and thus increased NAP1 phosphorylation.

- 2 After 48 h, collect the cells by trypsinization and wash the cell pellet with PBS once. Then, lyse the cells in lysis buffer.

- 3 Lyse the samples for 00:20:00 On ice . Then, clear the cell lysates by centrifugation at 20000 x g, 4°C, 00:10:00 .

30m



- 4 Determine the protein concentrations of the cleared protein lysates with the Pierce Detergent Compatible Bradford Assay Kit (23246, Thermo Fisher).

- 5 For both samples, negative control and NAP1-EGFP lysates, incubate 12 mg of cell lysate Overnight with 20 µL or GFP-Trap agarose beads (GTA-20, Chromotek).

20m



- 6 In the morning, wash the samples three times in lysis buffer. Then, resuspend the beads in protein loading dye, supplemented with 100 mM DTT. Finally, boil the samples for 00:05:00 at 95 °C .

5m



- 7 Load the samples on 4-12% SDS-PAGE gels (NP0322BOX, Thermo Fisher) with PageRuler Prestained protein marker (Thermo Fisher).

- 8 After the run, stain the SDS-PAGE gel with Coomassie, fix for 00:10:00 in fixation solution, and destain it Overnight in dH2O.

15m





- 9 Cut the band corresponding to NAP1-EGFP from the gel with a fresh scalpel and submit for mass spectrometry analysis.