



Dec 01, 2025

Immunoprecipitation with GFP-Trap Beads

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg3k6b7v25/v1>

Devin M. Fuller^{1,2}, Thomas Melia^{1,2}

¹Department of Cell Biology, Yale University School of Medicine, New Haven, CT;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research N



Devin M. Fuller

Yale University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg3k6b7v25/v1>

Protocol Citation: Devin M. Fuller, Thomas Melia 2025. Immunoprecipitation with GFP-Trap Beads. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kqdg3k6b7v25/v1>

Manuscript citation:

Fuller Devin M, Wu Yumei, Schueder Florian, Rasool Burha, Nag Shanta, Korfhage Justin L, Garcia-Milian Rolando, Melnyk Katerina D, Bewersdorf Joerg, De Camilli Pietro, Melia Thomas J (2025) ATG2A engages RAB1A and ARFGAP1 positive membranes during autophagosome biogenesis eLife 14:RP107316

<https://doi.org/10.7554/eLife.107316.1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 09, 2025

Last Modified: December 01, 2025

Protocol Integer ID: 126468

Keywords: tagged protein, trap bead, protein, experiments with gfp, gfp

Funders Acknowledgements:

National Institutes of Health

Grant ID: R01 GM100930

National Institutes of Health

Grant ID: R35 GM153482

National Institutes of Health

Grant ID: DA018343

National Institutes of Health

Grant ID: F31 AG079606

National Institutes of Health

Grant ID: F31 DK136246

Aligning Science Across Parkinson's

Grant ID: ASAP-025173

Human Frontier Science Program

Grant ID: LT000056/2020-C

Abstract

A protocol to perform co-immunoprecipitation (CoIP) experiments with GFP tagged proteins.



Materials

Cell culture materials:

DMEM (Thermo Fisher Scientific, 11965-092)
FBS (Thermo Fisher Scientific, 16140-071)
Penicillin/Streptomycin (10,000 U/mL; Thermo Fisher Scientific, 15140122)
PBS (Thermo Fisher Scientific, 10010023)
Earle's Balanced Salt Solution (EBSS; Thermo Fisher Scientific, 24010043)

Transfection Reagents:

Opti-Mem (Thermo Fisher Scientific, 31985062)
Lipofectamine 3000 (Thermo Fisher Scientific, L3000008)

Lysis Buffer:

lysis buffer (10 mM Tris pH 7.5, 150 mM NaCl, 0.5 mM EDTA, 0.5% NP-40) with Protease Inhibitor Cocktail
Tris (American Bio, AB02000-05000)
NaCl (Sigma-Aldrich, S9888-25G)
EDTA (Sigma-Aldrich, 3690)
Nonidet P40 Substitute (Roche, 11754599001)
Complete mini EDTA Free (Roche, 11836170001)

Other reagents:

GFP-Trap beads (Proteintech, gtma-20)

Cell culture and treatments

10m

- 1 Culture the HEK293 cells at  37 °C in 5% CO₂ and DMEM containing 10% FBS,  1000 U/mL penicillin, a  1 mg/mL streptomycin.
- 2 For any given experiment, plate the cells on 15 cm plates at such density so as to be approximately 95% confluent at the time of lysis.
- 3 If transfecting the sample, do so when the cells are approximately 70-80% confluent. Prepare the following reagents for transfection with lipofectamine 3000:
 - 3.1 For transfected plasmids, each 15 cm plate requires  10 µg of each target DNA construct. Prepare a master mix consisting of  10 µg of DNA,  700 µL of opti-mem low serum medium, and  20 µL of P3000 reagent per well. Prepare an additional mix consisting of  700 µL of opti-mem and  20 µL of lipofectamine 3000 per well. Combine the DNA + P3000 and lipofectamine solutions together, and allow to incubate for  00:10:00 at  Room temperature .

10m

Cell lysis and sample preparation

- 4 Supplement lysis buffer (10 mM Tris pH 7.5, 150 mM NaCl, 0.5 mM EDTA, 0.5% NP-40) with Protease Inhibitor Cocktail (Roche) and chill  On ice .
- 5 Aspirate media from cells, add PBS, and scrape cells using a cell scraper (Corning)  On ice
- 6 Pipette the resuspended cells into Eppendorf tube.
- 7 Centrifuge at  500 x g, 4°C, 00:03:00 , aspirate the excess PBS.
- 8 Pipette 150 µL lysis buffer onto cell pellet, pipette mixture up and down 10 times with a P-200 pipette tip

**Note**

NOTE: Take care not to introduce bubbles.

9 Incubate Eppendorf tube On ice for 00:05:00 .

10 Determine protein concentration in sample using Bradford Assay.

Immunoprecipitation**2h 10m 5s**

11 Prepare GFP-Trap beads for the IP reaction:

11.1 Pipette 25 μL of the GFP-Trap bead slurry into eppendorf tubes for each reaction.

11.2 Pipette 500 μL of lysis buffer onto the bead slurry. Mix the tube by inverting five times.

11.3 Precipitate the beads using a magnetic rack. Aspirate the lysis buffer, including the buffer trapped at the top of the tube.

11.4 Pipette 500 μL of lysis buffer onto the bead slurry. Mix the tube by inverting five times.

11.5 Precipitate the beads using a magnetic rack. Aspirate the lysis buffer, including the buffer trapped at the top of the tube.

12 Pipette equal amounts of lysate (by protein level) onto the washed GFP-Trap beads. Adjust the volume up to 500 μL .

13 Incubate the samples at 4 $^{\circ}\text{C}$ for 02:00:00 , rotating.

2h



- 14 Perform the following washes at Room temperature :
- 14.1 Aspirate the flow through. Pipette 500 μ L of lysis buffer onto the bead slurry. Mix the tube by inverting five times. Precipitate the beads using a magnetic rack. Aspirate the lysis buffer, including the buffer trapped at the top of the tube.
- 14.2 Pipette 500 μ L of lysis buffer onto the bead slurry. Mix the tube by inverting five times. Precipitate the beads using a magnetic rack. Aspirate the lysis buffer, including the buffer trapped at the top of the tube.
- 14.3 Pipette 500 μ L of lysis buffer onto the bead slurry. Mix the tube by inverting five times. Precipitate the beads using a magnetic rack. Aspirate the lysis buffer, including the buffer trapped at the top of the tube.
- 14.4 Pipette 500 μ L of lysis buffer onto the bead slurry. Mix the tube by inverting five times. Precipitate the beads using a magnetic rack. Aspirate the lysis buffer, including the buffer trapped at the top of the tube.
- 14.5 Pipette 500 μ L of lysis buffer onto the bead slurry. Mix the tube by inverting five times. Precipitate the beads using a magnetic rack.
- 14.6 Briefly centrifuge the samples for 00:00:05 . Aspirate the buffer. 5s
- 15 Elute the protein by heating the samples at 95 °C in 2x LDS loading buffer for 00:10:00 . Pipette the eluted sample into a new eppendorf tube or directly in to the gel. 10m

Gel electrophoresis and immunoblotting

2h 10m 5s

- 16 Prepare samples at desired concentration, add 1 mM DTT and 1x LDS loading buffer.
- 17 Incubate samples at 95 °C for 00:05:00 .
- 18 During this incubation, prepare gel apparatus with NuPAGE™ Bis-Tris Mini Protein Gels, 4–12%, 1.0–1.5 mm (ThermoFisher) and 1x MOPS running buffer.



19 Load samples into gel and run until dye front reaches bottom (180 V, roughly 60 min).

20 Remove gel and set up transfer cassette with nitrocellulose membrane.

21 Transfer at 30-35 V for  01:30:00 in 1x transfer buffer.

22 Block membrane with 5% BSA in PBST for  01:00:00 at  22 °C .

23 Incubate membrane with primary antibodies in 5% BSA and 0.03% sodium azid in PBST  Overnight at  4 °C .

Note

NOTE: Optimal primary antibody incubation time and temperature can be determined empirically for a given primary antibody

24 Wash membrane with PBST for five minutes. Repeat a total of 3 times.

24.1 Wash membrane for  00:05:00 with PBST (1/3).

24.2 Wash membrane for  00:05:00 with PBST (2/3).

24.3 Wash membrane for  00:05:00 with PBST (3/3).

25 Incubate membrane with secondary antibodies (1:5,000) in 5% milk in PBST for  01:00:00 at  22 °C .

26 Wash membrane with PBST. Repeat a total of 3 times.



- 26.1 Wash membrane for  00:05:00 with PBST (1/3).
- 26.2 Wash membrane for  00:05:00 with PBST (2/3).
- 26.3 Wash membrane for  00:05:00 with PBST (3/3).
- 27 Image membranes using a BioRad VersaDoc imaging system.

Protocol references

https://www.ptglab.com/products/pictures/pdf/gtma_Manual_GFP-Trap_Magnetic_Agarose.pdf

Acknowledgements

This work was supported by grants from the National Institutes of Health (R01 GM100930 and R35 GM153482 to TJM; R01 GM151829 to JB; DA018343 to PDC), F31 AG079606 to DMF and F31 DK136246 to JLK. This research was also funded in part through Aligning Science Across Parkinson's (ASAP-025173 to TJM and PDC) through the Michael J. Fox Foundation for Parkinson's Research (MJFF) and the Howard Hughes Medical Institute (HHMI; PDC). FS acknowledges support from the Human Frontier Science Program (LT000056/2020-C). JB acknowledges support by the Wellcome Leap Foundation. Imaging was supported by the Yale Center for Cellular and Molecular Imaging (both the fluorescence and electron microscopy facilities). We also thank the MS & Proteomics Resource at Yale University for providing the necessary mass spectrometers and the accompany biotechnology tools funded in part by the Yale School of Medicine and by the Office of The Director, National Institutes of Health (S10OD02365101A1, S10OD019967, and S10OD018034). The funders had no role in study design, data collection and analysis, decision to publish, or preparation of the manuscript.