



Jan 22, 2024

Immunostaining

DOI

dx.doi.org/10.17504/protocols.io.e6nvwd7n7lmk/v1

Cole S Sitron¹, Victoria A Trinkaus¹, F Ulrich Hartl¹

¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry



Cole Sitron

Max Planck Institute of Biochemistry

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwd7n7lmk/v1>

Protocol Citation: Cole S Sitron, Victoria A Trinkaus, F Ulrich Hartl 2024. Immunostaining. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.e6nvwd7n7lmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 16, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 93738

Keywords: ASAPCRN, cellular epitopes by immunofluorescence microscopy, immunofluorescence microscopy, immunofluorescence, immunostaining, cellular epitope, antibody

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Abstract

This is a protocol that describes how to use antibody staining to detect cellular epitopes by immunofluorescence microscopy.

Attachments



[260-2491.docx](#)

16KB

Materials



Cell culture

25-50% confluent wells of cells growing on coverslips in 24-well plates.



Equipment

- Platform shaker
- Microcentrifuge

Buffers & Reagents

- 4% paraformaldehyde (diluted from  Pierce™ 16% Formaldehyde (w/v), Methanol-free **Thermo Fisher Catalog #28908**) in 1X PBS (diluted from  PBS (10X), pH 7.4 **Thermo Fisher Scientific Catalog #70011051**)
- 1X PBS
- 0.1 % Triton X-100 in PBS
- Blocking Buffer:









	A	B
	Milk powder (Sucofin)	5%
	Triton X-100	0.1%
	PBS	pH 7.4

- 2 drops/ml  NucBlue™ Fixed Cell ReadyProbes™ Reagent **Thermo Fisher Catalog #R37606**
-  Fluorescence Mounting Medium **Agilent Technologies Catalog #S302380-2**
- Appropriate primary and fluorophore-conjugated secondary antibodies
- Microscope slides

Troubleshooting



Fixation


















- 1 Aspirate media and gently wash each well with 1X PBS. 
- 2 Aspirate PBS and place  500 μL of 4% paraformaldehyde onto each well. 
- 3 Fix at  Room temperature for  00:10:00 . 
- 4 Remove 4% paraformaldehyde and add  1 mL 1X PBS. 

Note




At this point, the plate can be stored for several months at  4 °C .

Immunostaining

2h 40m

- 5 Permeabilize  1 mL 0.1% Triton X-100 and gently shake at  Room temperature for  00:05:00 . 
- 6 Block in  1 mL Blocking Buffer for  01:00:00 . 
- 7 Dilute primary antibodies into Blocking Buffer and add  300 μL of diluted antibodies to each well. 
- 8 Shake at  4 °C  Overnight . 
- 9 Wash 3x with  1 mL PBS, gently shaking at  Room temperature each time for  00:05:00 .  





10 During washes, centrifuge tubes of secondary antibodies at  16000 x g for  00:10:00 at  4 °C .

10m





11 Dilute secondary antibodies into Blocking Buffer 1:500.





12 Incubate coverslips in  300 µL diluted secondary antibody for 2-3 hr at  Room temperature with gentle shaking, protected from light.



13 Wash with  1 mL PBS for  00:05:00 with gentle shaking.

5m



14 Add  1 mL diluted NucBlue and incubate with gentle shaking for  00:05:00 , protected from light.


5m



15 Wash 2x with  1 mL PBS for  00:05:00 with gentle shaking.

5m





16 Allow fluorescence mounting medium to come to  Room temperature .

17 Add one drop of fluorescence mounting medium onto microscope slides for each coverslip.



18 Mount coverslips onto slides.

19 Let slides cure at  Room temperature  Overnight , protected from light.

5m

