

Jul 13, 2024

## Immunostaining- Fluorescent

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlk88y6l5r/v1>

daniel.dautan daniel<sup>1,2</sup>, Per Svenningsson<sup>1,2</sup>

<sup>1</sup>Department of Clinical Neuroscience, Karolinska Institutet, 171 76 Stockholm, Sweden;

<sup>2</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA



Jacquelyn Haytayan

Weill Cornell Medicine

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk88y6l5r/v1>

**Protocol Citation:** daniel.dautan daniel, Per Svenningsson 2024. Immunostaining- Fluorescent. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.j8nlk88y6l5r/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** June 05, 2024

**Last Modified:** September 24, 2024

**Protocol Integer ID:** 101197

**Keywords:** ASAPCRN, immunofluorescence, mouse, brain, immunofluorescent visualization, mouse brain section, brain,

**Funders Acknowledgements:**

**Aligning Science Across Parkinson's**

Grant ID: 020608

## Abstract

Staining of mouse brain sections with immunofluorescent visualization.



- 1 Wash freshly sectioned slices of tissue 2-3 times with 1X PBS to remove OCT.
- 2 Block with 5% Normal donkey/Goat serum in 0.3% TRITON X-100 in 1X PBS for  01:00:00 at  Room temperature . 1h
- 3 Wash sections 5 times in 1X PBS.
- 4 Following washes, transfer sections to a primary antibody solution containing a mix of primary antibody diluted in 0.3% TRITON X-100 in 1X PBS with 1% Normal Donkey/Goat serum  Overnight at  4 °C . 1h
- 5 The following day, wash sections 5 times in 1X PBS.
- 6 Transfer sections to secondary solution containing the adequate secondary antibody diluted in 0.3% TRITON X-100 in 1X PBS with 1% Normal Donkey/Goat serum for  04:00:00 at  Room temperature . 4h
- 7 Wash sections 3 times in 1X PBS.
- 8 Mount sections on microscope slides with medium (Vectashield with 4',6-diamidino-2-phenylindole (DAPI), H-1800-10).