

Jan 06, 2026

Immunostaining

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlk141wg5r/v1>

Roberta Marongiu¹, Sabina Marciano²

¹Weill Cornell Medical College; ²weill cornell medicine



Jacquelyn Haytayan

Weill Cornell Medicine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk141wg5r/v1>

Protocol Citation: Roberta Marongiu, Sabina Marciano 2026. Immunostaining . protocols.io

<https://dx.doi.org/10.17504/protocols.io.j8nlk141wg5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 06, 2026



Last Modified: January 07, 2026

Protocol Integer ID: 237149

Keywords: ASAPCRN, IHC, Immunostaining , Parkinson's, Gene Therapy , microglial marker iba1, multiplexed fluorescent detection of dopaminergic neuron, synuclein aggregate, microglia, study of neuroinflammatory, neurodegenerative features in parkinsonian model, neurodegenerative feature, neuroinflammatory, synuclein, pathological markers in aof, immunohistochemistry, pd tissue, multiplexed fluorescent detection, endogenous peroxidase, dopaminergic neuron, immunohistochemistry protocol, pathological marker, tyrosine hydroxylase

Funders Acknowledgements:

Aligning Science Across Parkinson's
Collaborative Research Network
Weill Cornell Medicine

Abstract

This immunohistochemistry protocol was developed to visualize pathological markers in AOF-PD tissue, focusing on tyrosine hydroxylase (TH), microglial marker Iba1, GFAP, and DAT with quenching of endogenous peroxidases to reduce background. The protocol enables multiplexed fluorescent detection of dopaminergic neurons, microglia, and pathological alpha-synuclein aggregates, facilitating the study of neuroinflammatory and neurodegenerative features in Parkinsonian models.

Materials

Solutions:

PBS-T Triton 0.2 % 150mL: 300uL Triton in 150mL PBS

Immunostaining: AOF-PD: TH, Iba1, GFAP, DAT

- 1 **1) Washes** Wash 3x in PBS-T, 10min each (0.2% Triton in PBS)
- 2 **Quench** Incubate in 3% H₂O₂ in PBST for 2hrs at RT in the DARK! (x3) (can be done the day before and sections can be stored in PBS ON at 4degrees)
- 3 **Washes** Wash 3x in PBS-T, 10min each (0.2% Triton in PBS)
- 3.1 **All washing steps are performed for 10min and on a shaker**
- 4 **Blocking** Incubate in 3% BSA/PBS-T for 1 hr at room temperature.
- 5 **5) Primary Antibody** Incubate in primary antibody (diluted in blocking solution) overnight.
- 5.1 **Goat anti-Iba1 Iba1, WAKO, cat. 01919741 1:1500**
- 5.2 **Rabbit anti-TH Abcam, cat. Ab113, 1:1000**
- 5.3 **Mouse anti-GFAP ThermoFisher cat. PA1-10019 1:1500**
- 5.4 **Rat anti-DAT Millipore Cat #MAB369 1:500**
- 6 **Washes** Wash 3x in PBS-T, 10min each
- 6.1 **All washing steps are performed for 10min and on a shaker**



7 **Secondary Antibody** Incubate in secondary antibody (diluted in blocking solution) 1hr at room temperature.

7.1 **Donkey anti-goat 594 1:1000**

7.2 **Donkey anti-rabbit 488 1:500**

7.3 **Donkey anti-mouse 405 1:1000**

7.4 **Donkey anti-rat 647 1:1000**

8 **Washes** Wash 3x in PBS-T, 10min each

8.1 **All washing steps are performed for 10min and on a shaker**