

Jan 07, 2026

Immunostaining

DOI

<https://dx.doi.org/10.17504/protocols.io.14egn1z16v5d/v1>

Roberta Marongiu¹, Sabina Marciano²

¹Weill Cornell Medical College; ²weill cornell medicine



Jacquelyn Haytayan

Weill Cornell Medicine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn1z16v5d/v1>

Protocol Citation: Roberta Marongiu, Sabina Marciano 2026. Immunostaining . **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.14egn1z16v5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 07, 2026



Last Modified: January 07, 2026

Protocol Integer ID: 237151

Keywords: ASAPCRN, immunofluorescence, endogenous peroxidase activity, primary antibodies against tyrosine hydroxylase, fluorescent secondary antibody, tyrosine hydroxylase, procedure with hydrogen peroxide quenching, staining procedure, detection of th, hydrogen peroxide quenching, primary antibody, reduced autofluorescence, secondary antibody, bsa before overnight incubation, fluorescent, antibody

Funders Acknowledgements:

Aligning Science Across Parkinson's
Collaborative Research Network
Weill Cornell Medicine

Abstract

This protocol describes an immunofluorescence staining procedure with hydrogen peroxide quenching to reduce endogenous peroxidase activity, optimized for detection of TH and Myc/Cherry signals. Tissue sections are permeabilized with PBS-T, quenched with 3% H₂O₂ in the dark, and blocked with BSA before overnight incubation with primary antibodies against tyrosine hydroxylase (TH) and Myc. Fluorescent secondary antibodies and NeuroTrace are then applied for visualization, with thorough PBS-T washes between steps to minimize background. The protocol supports robust, specific labeling with reduced autofluorescence and background.

Materials

Solutions:

PBS-T Triton 0.2 % 150mL: 300uL Triton in 150mL PBS

Before start

All washing steps are performed for 10min and on a shaker

Immunostaining: AOF-PD: myc/cherry, TH quench

- 1 **Washes** Wash 3x in PBS-T, 10min each (0.2% Triton in PBS)
- 2 **Quench** Incubate in 3% H₂O₂ in PBST for 2hrs at RT in the DARK! (x3) (can be done the day before and sections can be stored in PBS ON at 4degrees)
- 3 **Washes** Wash 3x in PBS-T, 10min each (0.2% Triton in PBS)
- 4 **Blocking** Incubate in 3% BSA/PBS-T for 1 hr at room temperature.
- 5 **Primary Antibody** Incubate in primary antibody (diluted in blocking solution) overnight.
 - 5.1 **Rabbit anti-TH Abcam, cat. Ab113, 1:1000**
 - 5.2 **Mouse anti-Myc Cell Signaling, cat. 2276 1:1000**
- 6 **Washes** Wash 3x in PBS-T, 10min each
- 7 **Secondary Antibody** Incubate in secondary antibody (diluted in blocking solution) 1hr at room temperature.
 - 7.1 **Donkey anti-rabbit 647 1:500**
 - 7.2 **Donkey anti-mouse 488 1:1000**
 - 7.3 **Neurotrace 455 (1:200)**
- 8 **Washes** Wash 3x in PBS-T, 10min each

