

Feb 16, 2025

Immunostaining Mouse Tissue

 Forked from [Immunostaining Mouse Brain Tissue or Neuronal Cultures](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg3q7ypv25/v1>

Chelsie Steele¹, Yujie Fan¹

¹Caltech



Chelsie Steele

Caltech

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg3q7ypv25/v1>

Protocol Citation: Chelsie Steele, Yujie Fan 2025. Immunostaining Mouse Tissue . **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.kqdg3q7ypv25/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 07, 2025

Last Modified: February 16, 2025

Protocol Integer ID: 117844

Keywords: ASAPCRN, immunostaining mouse tissue, protein expression in mouse brain, mouse tissue, immunostaining procedure, cultured neuron, fixation of tissue, mouse brain, immunocytochemical method, protein expression, mouse, tissue, brain, fixation

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-020495

Abstract

This protocol describes immunohistochemical and immunocytochemical methods to analyze protein expression in mouse brain. In the protocol, we describe the fixation of tissue or cultured neurons, the immunostaining procedure, and collecting images for analysis.

Materials

Buffers, Reagents, Materials

- 4% paraformaldehyde (diluted from  Pierce & Warriner; 16% Formaldehyde (w/v), Methanol-free Thermo Fisher Catalog #28908) in 1X DPBS
- 1X DPBS
- 0.1 % Triton X-100 in DPBS
- Blocking Buffer: Normal Donkey Serum -Jackson Immune Research Labs - Cat No NC9624464
- DAPI - ThermoFisher Cat No D1306
- Prolonged Diamond Antifade Mountant - ThermoFisher Cat No P36965
- Appropriate primary and fluorophore-conjugated secondary antibodies
- Microscope slides

Troubleshooting

Collection of fixed tissues slices from mice

1d 0h 10m

- 1 Anesthetize mouse, pin onto dissection tray and open chest cavity
- 1.1 Using a perfusion pump, pierce needle into left ventricle and sever the right atrium; immediately begin perfusing cold 1X PBS followed by 4% paraformaldehyde (PFA) in PBS
- 1.2 Dissect out brain or wanted tissue into 15 mL tube containing cold 4% PFA in PBS
- 1.3 Incubate  Overnight at  4 °C in 4% PFA in PBS 12h
- 1.4 Incubate  Overnight at  4 °C in 30% sucrose in PBS 12h
- 1.5 Tissue can be kept in long-term storage at  -80 °C ; when ready for IHC, slice 100 µm sections using a cryostat (Leica CM3050 S)

Equipment

Cryostat

NAME

Leica

BRAND

Cryostat

SKU

<https://www.leicabiosystems.com/histology-equipment/cryostats/>^{LINK}

- 1.6 Once slices are collected, IHC is performed on free-floating sections in cold 1X DPBS with 5% Azide

Tissue staining

18h

- 2 Wash samples three times for 10 minutes at  Room temperature in 1X DPBS + 0.1% Triton X-100 on a rocker 30m
- 3 Block and permeabilize for 2 hours at  Room temperature in 1X DPBS containing 4% normal goat serum (NGS) and 0.1% Triton X-100 on a rocker 2h
- 4 Prepare primary antibody solution by diluting antibody at desired concentration in 1X DPBS containing 2% NGS and 0.1% Triton X-100
- 5 Incubate in primary antibody overnight at  4 °C 12h
- 6 Wash samples three times for 10 minutes at  Room temperature in 1X DPBS+ 0.1% Triton X-100 on a rocker 30m
- 7 Prepare secondary antibody solution by diluting antibody at 1:1000 in 1X DPBS containing 2% NGS and 0.1% Triton X-100
- 8 Incubate in secondary antibody for 2 hours at  4 °C in the dark 2h
- 9 Wash samples three times for 10 minutes at  Room temperature in 1X PBS + 0.1% Triton X-100 on a rocker 30m
- 10 Add DAPI 1:1000 to free floating samples for 10 minutes 10m
- 11 Wash samples two times for 10 minutes at RT in 1X DPBS + 0.1% Triton X-100 on a rocker 20m
- 12 Mount free floating slices onto a slide.
 - For all sample types, mount using Prolonged Diamond Antifade Mountant
- 13 Store samples at  4 °C for short-term storage or  -20 °C for long-term storage



Imaging

- 14 Collect images of fixed tissue using a microscope:
 - For high resolution: use Leica DMI8 Inverted Microscope with Spinning Disk