



May 06, 2025

Immunostaining of Brain Sections (100µm and 40µm)

 [The Journal of Cell Biology](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.261ge876og47/v1>

Pia Rodriguez¹

¹Duke University



Pia Rodriguez

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ge876og47/v1>

External link: <https://doi.org/10.1083/jcb.202410130>

Protocol Citation: Pia Rodriguez 2025. Immunostaining of Brain Sections (100µm and 40µm). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.261ge876og47/v1>

Manuscript citation:

Salazar MPR, Kolanukuduru S, Ramirez V, Lyu B, Manigault G, Sejourne G, Sesaki H, Yu G, Eroglu C (2025) Mitochondrial fission controls astrocyte morphogenesis and organization in the cortex. The Journal of Cell Biology 224(10). doi: [10.1083/jcb.202410130](https://doi.org/10.1083/jcb.202410130)



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 06, 2025

Last Modified: May 06, 2025

Protocol Integer ID: 217754

Keywords: immunostaining, brain section

Funders Acknowledgements:

Howard Hughes Medical Institute

Aligning Science Across Parkinson's Initiative

Grant ID: ASAP-020607

BRAIN Initiative

Grant ID: U19-NS123719

Chan Zuckerberg Initiative, Neurodegeneration Challenge Network Collaborative

Grant ID: DAF2021-237435

Chan Zuckerberg Initiative, Neurodegeneration Challenge Network Collaborative

Grant ID: DAF2018-191999

HHMI Gilliam Fellowship

Grant ID: GT15076

Abstract

General Immunohistochemistry protocol



Prepare Solutions

- 1 Make TBST (0.2% Triton): 49ml 1x TBS + 1ml 10% TritonX (note: make solution fresh to begin your staining. Make sure that Triton bottle has not been open for more than 2-4 weeks)
- 2 Blocking and antibody solution: 10% goat serum in TBST
- 3 Note: All incubation and washing steps done with shaking.

Blocking and Primaries

- 4 Wash sections 3× 10 min in 1x TBS. 1ml of solution per well is sufficient for 4-6 sections per well in a 24 well plate.
- 5 Incubate in blocking solution (10% goat serum in TBST) for 1hr at room temperature.
- 6 Prepare primary antibody solution by diluting antibodies in 10% goat serum TBST solution. Vortex briefly to mix. Spin antibodies down at 5 mins, max speed. (This step is done to pellet any insoluble antibody aggregates that might cause background staining.)
- 7 Incubate in primary antibody solution at 4°C on shaker for at least 2 nights. Can incubate up to 3 nights.

Secondaries and Mounting

- 8 Incubate in secondary antibodies for 3hr at room temperature.
- 9 Wash in TBST 2 × 10 mins, 3rd rinse for 30 min.
- 10 Prepare secondary antibody solutions by diluting 1:200 in 10% goat serum TBST and spinning down as described with primary antibodies.
- 11 Wash in TBST 3× 10 mins.



- 12 One by one, transfer sections first to a container filled with 2/3 TBS, 1/3 ddH₂O. Mount the sections (3 per slide) with the help of a clean brush on glass slides. Allow sections to dry for approximately 5-10 mins.
- 13 Use 1 drop of VectaShield mounting media (with DAPI) per section. Coverslip and then seal with nail polish.
- 14 Important! Wait at least 2hrs before imaging. This time is required for the sections to settle and achieve accurate z-stack imaging.