



Aug 20, 2020

Improved protocol of ELISA for determining the serum concentration of IL-23 in humans.

DOI

dx.doi.org/10.17504/protocols.io.bjzakup2e

Angel A Justiz-Vaillant¹

¹University of the West Indies St. Augustine

University of the West In...

angel.vaillant@sta.uwi.e...



Angel A Justiz-Vaillant

University of the West Indies St. Augustine

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bjzakup2e

Protocol Citation: Angel A Justiz-Vaillant 2020. Improved protocol of ELISA for determining the serum concentration of IL-23 in humans.. [protocols.io https://dx.doi.org/10.17504/protocols.io.bjzakup2e](https://dx.doi.org/10.17504/protocols.io.bjzakup2e)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 20, 2020

Last Modified: August 20, 2020

Protocol Integer ID: 40706



Abstract

A previous protocol was contaminated with IL-17. This is an improved version of an ELISA for quantification of Interleukin 23. It activates the production of IL-17A, which is a pro-inflammatory cytokine, produced by a group of T helper cell known as T- helper17 lymphocytes in response to the said IL-23 [1,2].

Reference

1. Cho ML, Kang JW, Moon YM, Nam HJ, Jhun JY, Heo SB, Jin HT, Min SY, Ju JH, Park KS, Cho YG, Yoon CH, Park SH, Sung YC, Kim HY (May 2006). Journal of Immunology. 176 (9): 5652–61.
2. Nobee A, Justiz-Vaillant A, Akpaka PE, Poon-King P (2016) Levels of Interleukin 17 and 23 in Patients with Systemic Lupus Erythematosus (SLE) in Trinidad and Tobago. Immunochem Immunopathol 2: 115. doi: 10.4172/2469- 9756.1000115

Materials

MATERIALS

☒ Set of one 96-well filter plate with 2 plate sealers **Millipore Catalog #MX-PLATE**

☒ Interleukin 23 by IBL corporation Germany



- 1 Ninety-six well ELISA plates are coated with monoclonal anti-human antibodies to interleukin-23 (IL-23).
- 2 Patient serum samples are added to the plates.
- 3 The plate is incubate x 1.30 hour at RT.
- 4 The plate is washed 4 times with PBS-tween 20 buffer.
- 5 The wells are incubated with a biotin conjugated anti-human IL-23 for 1.30 hour at RT.
- 6 The plates are washed again as above.
- 7 To the plate a peroxidase-labeled streptavidin conjugate is added and incubated for 1 hour at RT.
- 8 After a further washing procedure a substrate solution reactive is added and allowed to produced a colored reaction in positive controls.
- 9 The level of IL-23 in the sample is proportional to the colored product developed.
- 10 The addition of 3M H₂SO₄ stops the reaction.
- 11 The absorbance is measured at 450 nm.
- 12 The IL-23 concentration can be calculated by generating an standard curve using lineal regression.