


Aug 27, 2025

In situ species-level sequencing of seawater and sediment

DOI

<https://dx.doi.org/10.17504/protocols.io.261ge5endg47/v1>



Laura C. Martinez¹, Nallely Magaña-Montiel¹, Ricardo Grande¹, Liliana Pardo-López¹

¹Instituto de Biotecnología, UNAM



Laura C. Martinez

Instituto de Biotecnología UNAM

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ge5endg47/v1>

Protocol Citation: Laura C. Martinez, Nallely Magaña-Montiel, Ricardo Grande, Liliana Pardo-López 2025. In situ species-level sequencing of seawater and sediment . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.261ge5endg47/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: June 18, 2024

Last Modified: August 27, 2025

Protocol Integer ID: 102054

Keywords: Marine water, Marine sediment, DNA extraction, 16S gene, Nanopore , dna from marine sediment, level sequencing of seawater, marine sediment, identification of the microorganism, sequencing data, 16s gene dna library, preparation of 16s gene dna library, microorganism, oceanographic research vessel, minion mk1c sequencer, oxford nanopore, water sample, nanopore, seawater, oxford nanopore technology, sediment, dna, subsequent sequencing, level sequencing, methodology for the extraction, tris buffer, situ species, extraction, sequence, gulf of mexico

Funders Acknowledgements:

CONAHCYT Master Scholarship

Grant ID: 1321139

LANCO-CONAHCyT

Grant ID: APOYOS LNC-2023-80

DGAPA-PAPIIT

Grant ID: IG200223

Abstract

This protocol outlines a methodology for the extraction of DNA from marine sediment and water samples collected in the Gulf of Mexico. It further details the preparation of 16S gene DNA libraries and their subsequent sequencing using the Oxford Nanopore Technologies (MinION Mk1C sequencer) aboard an oceanographic research vessel. It has been observed that the addition of water to the reactions may result in the fragmentation of the sequence into two parts; therefore, Tris buffer is used instead of water. The sequencing data can be used for the identification of the microorganisms present in the samples.

Image Attribution

Marine Biotechnology Laboratory - IBT-UNAM



Materials

PCR Barcoding Expansion 1–96 Oxford Nanopore Technologies Catalog #EXP-PBC096
LongAmp Taq 2X Master Mix - 100 rxns New England Biolabs Catalog #M0287S
ZymoBIOMICS Lysis Solution Zymo Research Catalog #D4300-1-150
Q5 High-Fidelity 2X Master Mix - 100 rxns New England Biolabs Catalog #M0492S
ZymoBIOMICS DNA Binding Buffer Zymo Research Catalog #D4300-2-100
ZymoBIOMICS HRC Prep Solution Zymo Research Catalog #D4300-7-30
ZymoBIOMICS DNA Wash Buffer 2 Zymo Research Catalog #D4300-4-60
UltraPure™ Low Melting Point Agarose Thermo Fisher Scientific Catalog #16520050
ZymoBIOMICS DNA Wash Buffer 1 Zymo Research Catalog #D4300-3-50
Zymo-Spin III-HRC Filters Zymo Research Catalog #C1058-50
Zymo-Spin III-F Filter Zymo Research Catalog #C1057-50
Zymo-Spin IICR Columns Zymo Research Catalog #C1078-50
ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm) Zymo Research

Protocol materials

- ✕ ZymoBIOMICS DNA Binding Buffer **Zymo Research Catalog #D4300-2-100**
- ✕ Zymo-Spin III-HRC Filters **Zymo Research Catalog #C1058-50**
- ✕ ZymoBIOMICS HRC Prep Solution **Zymo Research Catalog #D4300-7-30**
- ✕ Zymo-Spin III-F Filter **Zymo Research Catalog #C1057-50**
- ✕ Zymo-Spin IICR Columns **Zymo Research Catalog #C1078-50**
- ✕ ZymoBIOMICS DNA Wash Buffer 2 **Zymo Research Catalog #D4300-4-60**
- ✕ UltraPure™ Low Melting Point Agarose **Thermo Fisher Scientific Catalog #16520050**
- ✕ ZymoBIOMICS DNA Wash Buffer 1 **Zymo Research Catalog #D4300-3-50**
- ✕ PCR Barcoding Expansion 1–96 **Oxford Nanopore Technologies Catalog #EXP-PBC096**
- ✕ LongAmp Taq 2X Master Mix - 100 rxns **New England Biolabs Catalog #M0287S**
- ✕ ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm) **Zymo Research**
- ✕ Q5 High-Fidelity 2X Master Mix - 100 rxns **New England Biolabs Catalog #M0492S**
- ✕ ZymoBIOMICS Lysis Solution **Zymo Research Catalog #D4300-1-150**

Safety warnings

- ! The flow cell must be kept at a temperature range of 2-8 °C to ensure proper functionality.



Before start

Since this protocol is carried out on board a research vessel, it is necessary to ensure that all necessary materials and equipment are on board before setting sail.

It has been observed that the addition of water to the PCR reactions may result in the fragmentation of the sequence; therefore, in this protocol Tris buffer is used instead of water.



DNA extraction

29m

- 1 Add the sample to the lysis tube ZR BashingBead. Add  750 μ L of  ZymoBIOMICS Lysis Solution **Zymo Research Catalog #D4300-1-150** to the tube and cap.

	A	B
	Sample type	Maximun input
	Soil	250 mg
	Liquid samples	250 μ l

For seawater samples, first use a membrane filter with the desired pore size to collect enough cells for the DNA extraction.

Note

To ensure optimal homogenization, the membrane filters were initially shredded and then placed in Falcon tubes, prior to addition of lysis solution.
The sediment sample from 5 cm exhibited excessive viscosity, the addition of water before add the lysis buffer, was necessary.

The DNA extraction was made following the ZymoBIOMIC DNA Miniprep Kit protocol.

For more information on the filtration process, please refer to the protocol "Water and sediment sampling in the sea to discover marine bacteria"
(<https://www.protocols.io/edit/water-and-sediment-sampling-in-the-sea-to-discover-dx4w7qxe>)

- 1.1 Vortex for  00:15:00 .

15m

Transfer the lysis solution in the Falcon tube to the

 ZR BashingBeadTM Lysis Tubes (0.1 & 0.5 mm) **Zymo Research**



Note

The sediment samples do not require to be transfer to another tube.

- 2 Centrifuge the lysis tube at  13600 rpm, Room temperature, 00:01:00

1m

Note

13,600 rpm is the maximum speed of the microcentrifuge



Equipment

MiniSpin

NAME

Microcentrifuge

TYPE

Eppendorf

BRAND

5452000010

SKU

<https://www.eppendorf.com/es-es/Productos/Centrifugación/Microcentrifugas/MiniSpin-MiniSpinplus-p-PF-240750>

LINK



- 3 Transfer  400 μL supernatant to the  Zymo-Spin III-F Filter **Zymo Research Catalog #C1057-50** and placed it in a collection tube. Centrifuge at  10000 rpm, Room temperature, 00:01:00 . Discard the  Zymo-Spin III-F Filter **Zymo Research Catalog #C1057-50**

1m

Note

For the sediment samples 350 μL were transferred.

- 4 Add  1200 μL of  ZymoBIOMICS DNA Binding Buffer **Zymo Research Catalog #D4300-2-100** to the filtrate in the collection tube and mix well.
- 5 Transfer  800 μL of the mixture from step 4 to a  Zymo-Spin IICR Columns **Zymo Research Catalog #C1078-50** in a collection tube. Centrifuge at  13600 rpm, Room temperature, 00:01:00

1m



5.1 Discard the flow through from the collection tube and repeat step 5.

6 Add  400 μL   ZymoBIOMICS DNA Wash Buffer 1 **Zymo Research Catalog #D4300-3-50** to the  Zymo-Spin IICR Columns **Zymo Research Catalog #C1078-50** in a new collection tube . Centrifuge at  13600 rpm, Room temperature, 00:01:00 . Discard the flow through from the collection tube

7 Add  700 μL to the   Zymo-Spin IICR Columns **Zymo Research Catalog #C1078-50** in a collection tube. Centrifuge at  13600 rpm, Room temperature, 00:01:00 . Discard the flow through from the collection tube

8 Add  200 μL of   ZymoBIOMICS DNA Wash Buffer 2 **Zymo Research Catalog #D4300-4-60** to the  Zymo-Spin IICR Columns **Zymo Research Catalog #C1078-50** in a collection tube. Centrifuge at  13600 rpm, Room temperature, 00:01:00

9 Transfer the  Zymo-Spin IICR Columns **Zymo Research Catalog #C1078-50** to a clean 1.5 ml microcentrifuge tube.

9.1 Add  100 μL (50 μL minimum) of steril miliQ water directly to the column matrix and incubate for 1 minute

Note

The original protocol use ZymoBIOMICSTM DNase/RNase Free Water, but miliQ water works as well.

9.2 Centrifuge at  13600 rpm, Room temperature, 00:01:00 to elute the DNA. 



- 10 Place a  Zymo-Spin III-HRC Filters **Zymo Research Catalog #C1058-50** in a new collection tube.

Note

This step is to prepare the filter through which the DNA will pass.

- 10.1 Add  600 μ L of  ZymoBIOMICS HRC Prep Solution **Zymo Research Catalog #D4300-7-30** to the  Zymo-Spin III-HRC Filters **Zymo Research Catalog #C1058-50**

- 10.2 Centrifuge at  10000 rpm, Room temperature, 00:03:00

3m

- 11 Transfer the elute DNA (from step 9) to a prepared  Zymo-Spin III-HRC Filters **Zymo Research Catalog #C1058-50** (from step 10) in a clean 1.5 ml microcentrifuge tube.

- 11.1 Centrifuge at  13600 rpm, Room temperature, 00:03:00 . Do not discard the flow through because it is the DNA. This DNA is suitable for PCR.


3m



- 12 To verify DNA extraction run the product on an agarose gel using

 UltraPure™ Low Melting Point Agarose **Thermo Fisher Scientific Catalog #16520050**

[M] 0.5 % (v/v) at 70 V for 20 min.



Note

The electrophoresis equipment was placed on a special table to ensure that the movement of the vessel did not affect the gel.

12.1 Visualize the gel in a transilluminator.

PCR for 16S rRNA gene amplification

1h 22m

13 Mix the following reagents in a 0.2 ml PCR tube



It has been observed that the addition of water to the reaction mixture can result in the cleavage of DNA in two fragments. Therefore, it is recommended to add Tris buffer. In the following steps of this protocol, **Tris buffer will be used insted of water and will be specified in the reagents.**

A	B
Reagents	Volume
Master Mix	10 μ l
Primers	2 μ l



A	B
DNA	8 µl
Tris Buffer	Bring to 20 µl

The



Q5 High-Fidelity 2X Master Mix - 100 rxns **New England**
Biolabs Catalog #M0492S

was used for the amplification.

The primers used for the amplification were:

A	B
Primer	Sequence
27F	TTT CTG TTG GTG CTG ATA TTG CAG AGT TTG ATC MTG GCT CAG
1492R	ACT TGC CTG TCG CTC TAT CTT CCG GTT ACC TTG TTA CGA CTT

Note

- Between 10 and 20 ng of DNA are required for the reaction; however, the vessel lacks the necessary we do not have equipment for DNA quantification. For this reason, 8 µl of the extracted DNA were added.

- 13.1 Use the miniPCR equipment for the amplification.
Using the miniPCR app, set the following cycle conditions.

1h 22m

Software

miniPCR app for android

NAME

miniPCR

DEVELOPER

<https://www.minipcr.com/downloads/>

SOURCE LINK

Equipment

miniPCR® mini8 thermal cycler

NAME

Thermal cycler

TYPE

miniPCR®

BRAND

mini8

SKU

<https://www.minipcr.com/>

LINK




 01:22:00

	A	B	C
	Temperature	Time	No. of cycles
	98 °C	2 mins	1
	98 °C	30 sec	30
	55 °C	30 sec	
	72 °C	1.5 mins	
	72 °C	5 mins	1
	10 °C	Hold	See note below

Note

Due to limitations inherent to the miniPCR instrument to set at 10°C, it was necessary to place the PCR tubes were placed on ice at 10°C following the final cycle at 72°C.



DNA purification after 16S amplification

30m

14 Place the AMPure Beads (AXP) at Room temperature for 00:30:00 before use it

30m

14.1 Resuspend the AMPure XP Beads (AXP) by vortexing.

15 Transfer the DNA sample to a clean 1.5 ml Eppendorf DNA LoBind tube.

15.1 Use a 0.6X ratio of the resuspended the AMPure XP Beads to the PCR product and mix by flicking the tube.

In this case, 12 μL of the beads were added.

Note

Example: if the final volume is 20 μL (20×0.6), 12 μL of AMPure XP Beads will be add.

15.2 Mix by vortexing and incubate for 00:05:00 at room temperature.

5m

16 Prepare 500 μL of fresh 80% ethanol in nuclease-free water.

17 Spin down the sample and pellet on a magnet until supernatant is clear and colourless. Keep the tube on the magnet for 00:05:00 and pipette off the supernatant.


5m



Note

- Sometimes, it may be necessary to slightly rotate the tubes on the magnet to facilitate the faster production of pellets.
- To open the tubes, it is recommended that both hands be used, making sure the beads do not move (one hand holds the magnet with the tubes, while the other hand opens the tube).

18 Keep the tube on the magnet and **wash the beads with** 200 μL **of freshly prepared 80% ethanol** without disturbing the pellet. Remove the ethanol using a pipette and discard.



Note

Add the ethanol along the walls, not directly onto the beads. The beads should not dry in this step; if the washing takes too long, close the tubes.

18.1 Repeat this step (step 18).

19 Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for  00:02:00 .

2m

20 Remove the tube from the magnetic rack, resuspend the pellet in  15 μL of Tris buffer and mix.
Incubate for  00:02:00  Room temperature .

2m

21 Spin down and place the tube back on the magnet for  00:02:00 or until the eluate is clear and colourless.

2m

- 21.1 Keep the tube on the magnet, pipette off the supernatant and transfer to a clean 1.5 ml Eppendorf tube.
- 22 To verify the amplification of the 16S gene, run the PCR product on an agarose gel using

20m



 UltraPure™ Low Melting Point Agarose **Thermo Fisher Scientific Catalog #16520050**

[M] 0.5 % (v/v) at 70 V for  00:20:00 .

Library preparation: Barcoding PCR

- 23 For the library preparation the



 PCR Barcoding Expansion 1–96 **Oxford Nanopore Technologies Catalog #EXP-PBC096**

was used.

This kit contains the necessary barcodes to prepare the library.

Mix the following reagents in a 0.2 ml PCR tube

A	B
Reagent	Volume
PCR Barcode	1 µl
DNA	4 µl + 20 µl Tris Buffer
Master mix	25 µl
Final Volume	50 µl

For each sample add one PCR Barcode (We used the barcodes A12, B12, C12 and D12).

Mix by pipetting and spin down.

Note

As the PCR Barcoding Expansion kit recommend, the

 LongAmp Taq 2X Master Mix - 100 rxns **New England Biolabs Catalog #M0287S**

was used (Master Mix).



23.1 Amplify using the following cycle conditions

	A	B	C
	Temperature	Time	No. of cycles
	98 °C	2 mins	1
	98 °C	15 sec	13
	62 °C	15 sec	
	72 °C	1.5 mins	
	72 °C	5 mins	1
	10 °C	Hold	

24 Purify the barcoded DNA following the steps listed in the section "**DNA purification after 16S amplification**"

25 Prepare 1 µg of pooled barcoded libraries and adjust the volume to 50 µL with TRIS buffer.



Note

Due to the inability to quantify the DNA could not be quantified on the vessel, so **12 µL of each barcoded library (4 libraries) were added**, adjusting to a final volume of 50 µL with TRIS buffer (48 µL of DNA and 2 µL of TRIS buffer).

End-prep

30s

26 Thaw all reagents on ice.

26.1 Invert the reagent tubes to ensure they are well mixed but do not vortex the Ultra II End Prep Enzyme Mix.

30s

If the Ultra II End Prep Buffer have precipitate, allow the mixture to come to room temperature and pipette the buffer up and down until break up the precipitate, followed by vortexing the tube for  00:00:30 .

26.2 Spin down the tubes before opening

27 Mix the following reagents in a 0.2 ml PCR tube

A	B
Reagent	Volume
DNA	50 µL
Ultra II End-prep Reaction Buffer	7 µL
Ultra II End-prep Enzyme Mix	3 µL
Final Volume	60 µL

Mix by pipetting

27.1 Incubate at the following conditions:



A	B
Temperature	Time
20°C	5 mins
65°C	5 mins

The miniPCR cannot be set at 20°C so we use a water bath at the indicated temperature. The temperature was measured with a thermometer.



Using the "**Heat Block**" function of the miniPCR, the tubes were incubated at 65°C.

- 28 Purify the DNA following the same steps listed in the section "**DNA purification after 16S amplification**".

Note

The necessary volume of beads in this step is 36 μ l ($60 \times 0.6 = 36$).

- 28.1 After the last wash (step 19) resuspend the DNA in 62 μ l of TRIS buffer and follow the steps 20-22 (pipette off 60 μ l of the supernatant).

Adapter ligation and clean-up

20m 30s

- 29 Spin down the follow reagents: **Adapter Mix H (AMX H)** and **Quick T4 Ligase**. Place it on ice.
Thaw Ligation Buffer (LNB) at room temperature, spin down and mix by pipetting. Place on ice immediately after thawing and mixing.
Thaw the Elution Buffer (EB) at room temperature, mix by vortexing, spin down and place on ice.
- 30 Mix the following reagents in a clean 1.5 ml Eppendorf DNA LoBind tube. Add all the reagents in the order of the following table.



A	B
Reagent	Volume
DNA	60 µL
Ligation Buffer (LNB)	25 µL
NEBNext Quick T4 DNA Ligase	10 µl
Adapter Mix H (AMX H)	5 µl
Final volume	100 µl

After adding all the reagents, do not pipette to mix them, but gently tap with your finger.

31 Resuspend the AMPure XP beads by gentle agitation.

32 Transfer the DNA sample to a clean 1.5 ml Eppendorf tube.

32.1 Use a 0.4 ratio of the resuspended the AMPure XP Beads to the final preparation reaction and mix by gently moving the tube.

In this case the final volume is 100 µl, so 40 µl of beads were added.

32.2 Vortex and incubate for  00:05:00 at room temperature.

5m

33 Spin down to pellet the beads and place the tube on a magnet until the supernatant is clear and colorless. Incubate for  00:05:00 .

5m

33.1 Keep the tube on the magnet and remove the supernatant.

34 Keep the tube on the magnet and wash the beads with 200 µl of **Short Fragment Buffer (SFB)**. Mix gently and wait for  00:00:30 .

30s

Note

The Short Fragment Buffer is used to purified fragments of all sizes, for DNA fragments of 3 kb or longer, use the Long Fragment Buffer (LFB).



34.1 Pipette of the SFB and discard it.

34.2 Repeat this step (step 34).

35 Spin down and place the tube back on the magnet. Remove any remaining SFB. Let it dry for 2 minutes.

36 Remove the tube from the magnetic rack, resuspend the pellet in 15 µl of **Elution Buffer (EB)** and mix.

Incubate for  00:10:00 at  Room temperature .

10m

36.1 Spin down and place the tube back on the magnet for 2 minutes until the supernatant is clear and colourless.

37 Keep the tube on the magnet and transfer 13 µl of the supernatant to a clean 1.5 ml Eppendorf tube.

Priming and loading the SpotON flow cell

5m

38 Thaw the **Sequencing Buffer II (SBII)**, **Loading Beads II (LBII)**, **Flush Tether (FLT)** and one tube of **Flush Buffer (FB)** at  Room temperature before mixing the reagents by vortexing and spin down at  Room temperature .

39 To prepare the flow cell priming mix, add 30 µl of Flush Tether (FLT) directly to the tube of Flush Buffer (FB), and mix by vortexing at  Room temperature .

40 Open the MinION device lid andx slide the flow cell under the clip.



Equipment

MinION	NAME
Sequencer	TYPE
Oxford Nanopore Technologies	BRAND
MinION 1B / MinION 1C	SKU

Note

Press the flow cell firmly to ensure proper contact. **The yellow area of the cell should never be dry.**

- 41 Slide the flow cell priming port cover clockwise to open the priming port. After opening the priming port, check for a small air bubble under the cover.
- 41.1 Set a P1000 pipette to 200 µl and insert the tip into the priming port.
- 41.2 Turn the wheel until the dial shows 220-230 µl, to draw back 20-30 ul, or until you can see a small volume of buffer entering the pipette tip.

Note

Do not draw back more than 20-30 µl.
Visually check that there is continuous buffer from the priming port across the sensor array.

- 42 Slowly load 800 µl of the **priming mix** into the flow cell via the **Priming port**, avoiding the introduction of air bubbles. Wait for  00:05:00 . During this time, prepare the library for loading.
- 43 Mix the contents of the Loading Beads II (LBII) by pipetting.

5m

Note

The LBII settle very quickly, so mixed immediately before use.

44



In a new tube, prepare the library for loading as follows:

A	B
Reagent	Volume
Sequencing Buffer II (SBII)	37.5 µl
Loading Beads II (LBII)	25.5 µl
DNA library	See note below
Final volume	75 µl

Note

The recommended loading range is between 5-10 fmol adjusted to a final volume of 12 µl with the elution buffer. However, the DNA could not be quantified; therefore 12 µl of the DNA library were added undiluted.

It is important to load the library onto the flow cell immediately after adding the Sequencing Buffer II (SBII) and Loading Beads II (LBII).

45

Complete the flow cell priming by lifting the **SpotON sample port cover** to make the SpotON sample port accessible and **load 200 µl of the priming mix into the flow cell priming port** (not the SpotON sample port), avoiding the introduction of air bubbles.

**Note**

The SpotON sample port and the priming port are different. Be sure of add the priming mix into the priming port.

- 46 Add 75 µl of the prepared library to the flow cell via the **SpotON sample port in a dropwise fashion.**

Be sure each drop flows into the port before adding the next.

Note

Mix the prepared library by pipetting up and down just prior to loading

- 47 Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the priming port and replace the MinION device lid.

Start sequencing

- 48 To start the sequencing, the following kits were selected:
- Ligation Sequencing Kit (SQK-LSK110)
 - PCR Barcoding Expansion 1-96 (EXP-PBC096)

Flow Cell Cleaning

- 49 If the flow cell is going to be reuse, it is neccesary to whash it using the Flow Cell Wash Kit.

Mix the following reagents in a clean Eppendorf tube

	A	B
	Reagent	Volume
	Solution A	20 µl
	Solution B	380 µl
	Final volume	400 µl



- 50 The process is similar to the washing of the cell mentioned in the "**Priming and loading the SpotON flow cell**".
Open the port.
- 50.1 Set a P1000 pipette to 200 μ l, insert the tip into the priming port and remove between 15-20 μ l.
- 50.2 Add the 400 μ l of the mixed solution
- 50.3 Load 400 μ l of the **Storage buffer** after the washing if the flow cell is going to be storage. If another sample is to be run, only the washing buffer is added, and then the entire cell preparation protocol is repeated and the next sample is loaded.