

Aug 01, 2023

In-vitro GCase Activity Assay

DOI

dx.doi.org/10.17504/protocols.io.dm6gp3b7dvzp/v1

Tae-Un Han¹, sidranse¹

¹Medical Genetics Branch, National Human Genome Research Institute, National Institutes of Health, USA.



Tae-Un Han

NIH/NHGRI

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gp3b7dvzp/v1>

Protocol Citation: Tae-Un Han, sidranse 2023. In-vitro GCase Activity Assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.dm6gp3b7dvzp/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 01, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 85785

Keywords: ASAPCRN, GCase, GCase activity, gcase activity assay, gcase activity, containing protease inhibitor, protease inhibitor, assay condition, protein lysate, cell lysate, protein

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000458

Abstract

This is version 2 of in-vitro GCase activity assay.

We optimized the assay condition to enhance sensitivity and specificity.

Both cell lysates and protein lysates prepared in M-buffer containing protease inhibitor and 0.25% TritonX can use this protocol.

Attachments



GCase Activity Assay...

147KB

Materials

- Cell lysates and protein lysates prepared in M-buffer containing protease inhibitor and 0.25% TritonX
- 0.2M Na₂HPO₄
- 0.1M citrate
- protease inhibitor tablet (Roche cOmplet mini)
- Triton X-100 solution
- 384 well plate (flat bottom, black)
- 10 mM CBE (5 mg CBE in 3086 ul DMSO)
- DMSO
- aluminum foil
- sodium taurocholate powder
- 1 M 4-Methylumbelliferyl-B-D-glucoside (4-MU-G, 338 mg per 1ml DMF)
- 5 M NaOH
- glycine
- microplate reader

Troubleshooting

- 1 Mix 20 mL 0.2M Na₂HPO₄ (0.852 g in 30 mL water) with 14 mL 0.1M citrate (0.576 g in 30 mL water) to make M-buffer (5.4 , no additional pH adjustment required).
- 2 Dissolve one protease inhibitor tablet (Roche cOmplet mini) in 10 mL M-buffer , then add Triton X-100 solution to 0.25 % (v/v) (e.g. 25 μ L in 10 mL) and 0.2 Mass / % volume of sodium taurocholate to make active GCase buffer.
- 3 Set up the desired plate layout in a 384 well plate (flat bottom, black). There should be two sections, as each sample must be prepared and assayed both with and without CBE (CBE: GCase1 inhibitor).
- 4 Prepare 0.8 millimolar (mM) CBE by diluting 10 millimolar (mM) CBE in DMSO with the GCase buffer (GCase buffer:CBE = 92:8).
- 5 Prepare CBE-free carrier solution in the same volume of GCase buffer/DMSO (GCase buffer:DMSO = 92:8).
- 6 Pipette 10 μ L protein lysate diluted with GCase buffer into wells of a 384- well plate. Four replication sets should be run (sample concentration: 0.7 mg/mL ~ 1.2 mg/mL). Protein concentration should be adjusted to be similar between control and experiment groups using GCase buffer.
- 7 Add 5 μ L 0.8 mM CBE solution to the CBE-positive wells or the same volume of CBE-free carrier solution to the CBE-negative wells.
- 8 Cover the plate with aluminum foil and briefly centrifuge. Incubate shaking: 600 rpm, 37°C, 00:15:00 .
- 9 During incubation, prepare 4-Methylumbelliferyl-B-D-glucoside (4-MU) diluent by diluting 1 M 4-Mu (338 mg per 1ml DMF) with GCase buffer to a final concentration of 2.5 millimolar (mM) (1:400 dilution).
- 10 After the CBE incubation, spin down the plate and add 15 μ L assay buffer with 2.5 mM 4-MU to reach a total volume of 30 μ L in each well.
- 11 Cover the plate with aluminum and briefly centrifuge, and Incubate shaking: 450 rpm, 37°C, 01:00:00 .



- 12 Prepare stop solution by adding 4 mL 5M NaOH and 1.877 g glycine up to 25 mL in water. Final concentration of glycine is 1 Molarity (M) (pH 10.5)
- 13 After incubation, spin down the plate again and Add 30 µL stop solution to each well.
- 14 Read the 4-MU fluorescence with a microplate reader (Excitation: 365 nm; Emission: 449 nm; Cutoff: 435nm; 3 reads/well).
- 15 GCase activity in each protein lysate can be calculated as below.

$$\frac{\text{fluorescence of CBE} - \text{free sample} - \text{fluorescence of sample with CBE}}{\text{protein concentration (measured by BCA)}}$$