

Mar 30, 2023

In vitro GCase activity assay (total cell lysate)

 In 2 collections

DOI

dx.doi.org/10.17504/protocols.io.261ge3767l47/v1

Federico Bertoli¹, Michela Deleidi¹

¹German Center for Neurodegenerative Diseases (DZNE), Tübingen, 72076 Germany



Federico Bertoli

ASAP

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ge3767l47/v1>

Protocol Citation: Federico Bertoli, Michela Deleidi 2023. In vitro GCase activity assay (total cell lysate). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.261ge3767l47/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 28, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 79530

Keywords: In vitro GCase activity assay, total cell lysate, ASAPCRN, glucocerebrosidase, assay detects gba activity, lysosomal enzyme, gcase activity, glucosylceramide, gba1 inhibitor, total cell lysate, hydrolysis of glucosylceramide, cell lysate, reacts with cell lysate, glcCer, glucose, assay, lysate, sphingolipid, cell, hydrolysis

Abstract

Glucocerebrosidase is a lysosomal enzyme that catalyzes the hydrolysis of glucosylceramide (GlcCer), a membrane glyco-sphingolipid, to ceramide and glucose. This assay detects GBA activity by using a fluorogenic substrate that reacts with cell lysates previously treated with or without CBE (GBA1 inhibitor).

Attachments



[ggmvbqjbx.pdf](#)

615KB

Materials

Reagents

-  4-Methylumbelliferyl β -D-glucopyranoside **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M3633**
-  Condurotol-b-epoxide **Merck Millipore (EMD Millipore) Catalog #234599**
-  AMP-Deoxynojirimycin (CAS 216758-20-2) **Catalog #sc-223780**

▪ 1% Triton Base Buffer:

	A	B	C
	1% Triton Base Buffer	Final concentration	Amount
	Triton X-100	1%	0.5 mL
	5 M NaCl	150 mM	1.5 mL
	1 M HEPES pH 7.4	20 mM	1 mL
	0.5 M EDTA	1 mM	100 μ L
	1 M MgCl ₂	1.5 mM	75 μ L
	100% glycerol	10%	5 mL
	Milli-Q H ₂ O	n/a	41.825 mL

▪ 1% Triton extraction buffer:

	A	B	C
	1% Triton Extraction Buffer	Final concentration	Amount
	1% Triton Base Buffer	n/a	4.425 mL
	PIC	n/a	½ tablet



	A	B	C
	500 mM NaF	50 mM	500 μ L
	200 mM Na ₃ VO ₄	2 mM	50 μ L
	0.1 M PMSF	0.5 mM	25 μ L


■ **Mcllvaine Buffer:**

	A	B	C
	pH	0.2 M NaHPO₄ (mL)	0.1 M citric acid (mL)
	6.0	12.63	7.37

Troubleshooting



Sample Lysis

1 Suspend samples in  50 μL of 1% Triton extraction buffer.



2 Homogenize with a Dounce homogenizer for 25 strokes.

3 Rotate samples for  00:30:00 at  4 $^{\circ}\text{C}$.

30m




4 Centrifuge at  13500 x g ,  4 $^{\circ}\text{C}$ for  00:15:00 .

15m






5 Collect supernatants.

Substrate preparation

6 Add  20.30 mg 4-Methylumbelliferyl- β -D-glucopyranoside for  10 mL ddH₂O of substrate ( 6 millimolar (mM)).





7 Incubate at  55 $^{\circ}\text{C}$ and vortex every  00:05:00 until dissolved (approx.  00:30:00).

35m



8 Store at  4 $^{\circ}\text{C}$ until needed.

Sample preparation

9 Add the equivalent of  10 μg total protein in ddH₂O to reach a final  45 μL volume.




**Note**

For each sample





- 10 Add to each  25 μL McIlave Buffer  6 and mix it.

**Note**

For GBA2 inhibition,  5 nanomolar (nM) AMP-Deoxynojirimycin




- 11 Divide the overall  70 μL volume into two tubes ( 35 μL each).



- 11.1 Incubate one tube with  5 μL CBE  1 millimolar (mM) at  Room temperature for  00:30:00 .


30m



- 11.2 Incubate the other one with  5 μL ddH₂O at  Room temperature for  00:30:00 .

30m

**Enzymatic reaction**


- 12 Add  25 μL substrate to each reaction tube.



- 13 Incubate at  37 °C for  02:00:00 .



2h

**Measurement**

- 14 Take  10 μL of each reaction tube into a 96-well plate (in triplicate).





15 Add  90 μ L  0.2 Molarity (M) glycine  10.2 to each well to stop the reaction.



16 Measure fluorescence: Excitation 355nm, Emission 460nm.

Note

GBA1 activity is obtained by subtracting the background and GBA2 activity from the total GCase activity.