

Aug 26, 2023

In vitro GTPase activity

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g74qwqgwz/v1>

Xinbo Wang^{1,2}, Pietro De Camilli^{1,2}

¹1. Departments of Neuroscience and of Cell Biology, Howard Hughes Medical Institute, Program in Cellular Neuroscience, Neurodegeneration and Repair, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

²2. Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815



xinbo.wang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g74qwqgwz/v1>

Protocol Citation: Xinbo Wang, Pietro De Camilli 2023. In vitro GTPase activity. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.q26g74qwqgwz/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 19, 2022

Last Modified: May 31, 2024

Protocol Integer ID: 68879

Keywords: GTPase, LRRK2, In vitro, ASAPCRN, gtpase activity testing of purified Irrk2, gtpase activity, gtpase activity testing, purified Irrk2, gtpase activity this protocol details method, Irrk2

Abstract

This protocol details methods for the *in vitro* GTPase activity testing of purified LRRK2.

Attachments



[iut9bvk9p.docx](#)

16KB

Materials

-  EnzChek® Phosphate Assay Kit **Thermo Fisher Catalog #E6646**
-  96-well Clear Flat Bottom Polystyrene TC-treated Microplate Individually Wrapped with Low Evaporat **Corning Catalog #3595**
- Microplate reader (Synergy H1; BioTek)



In vitro GTPase activity

1h 16m

- 1 Set up the reaction mixtures in a 96-well plate with  5 μL 20 $\times$ reaction buffer ( 1 Mass Percent Tris-HCl,  20 millimolar (mM) MgCl_2 ,  7.5 and  2 millimolar (mM) sodium azide),  200 micromolar (μM) 2-amino-6-mercapto-7-methylpurine riboside, 0.1 U purine nucleoside phosphorylase, and  9 micromolar (μM) LRRK2 protein or  0.8 micromolar (μM) Dynamin1 or buffers for the control in separate wells. 

Note

Note: For best measurement results, we usually use a total volume of  80 μL -  100 μL .

- 2 Preincubate samples in the microplate reader for  00:30:00 at  37 $^{\circ}\text{C}$. 

30m

- 3 Add  0.5 millimolar (mM) GTP (final concentration) to initiate the reaction. 

- 4 Immediately begin reading absorbance at  360 nm , every  00:01:00 over  00:45:00 at  37 $^{\circ}\text{C}$. 

46m

- 5 For data analysis, subtract the last values determined before GTP was added from the corresponding values for the experimental reaction.