



Dec 06, 2024

In vitro human tyrosinase inhibitory assay (human melanoma cell lysate based)

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygxw97dv8j/v1>

Gengxuan Shi¹

¹Institute of Biomedicine and Glycomics, Griffith University



Gengxuan Shi

Institute of Biomedicine and Glycomics, Griffith University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygxw97dv8j/v1>

Protocol Citation: Gengxuan Shi 2024. In vitro human tyrosinase inhibitory assay (human melanoma cell lysate based). [protocols.io https://dx.doi.org/10.17504/protocols.io.kxygxw97dv8j/v1](https://dx.doi.org/10.17504/protocols.io.kxygxw97dv8j/v1)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: December 03, 2024

Last Modified: December 06, 2024

Protocol Integer ID: 113708

Keywords: human tyrosinase inhibitory assay, human melanoma cell lysate, simple hstyr assay, based tyr assay, inhibitory assay, pigmented human melanoma cell line, known hstyr inhibitor, human tyrosinase, human melanoma cell line, hstyr inhibitor, tyr inhibitor, tyr assay, hstyr inhibition, skin hyperpigmentation, mushroom tyr, kojic acid, assay, tyrosinase, inhibitor, hstyr, using cell lysate, cell lysate, human tyr, melanoma

Abstract

Tyrosinase (TYR) inhibitors are required to treat skin hyperpigmentation. Currently live cell-based TYR assays and mushroom TYR *in vitro* assays are the common methods used to screen for TYR inhibitors. However, these methods are either time consuming and expensive or are not human TYR (*hsTYR*) specific. Here, we describe a simple *hsTYR* assay using cell lysate prepared from pigmented human melanoma cell lines that takes less than 3 hours to complete after collecting cell pellets. We confirmed the assay is species specific by using a known *hsTYR* inhibitor, kojic acid, as a positive control, while arbutin, which inhibits mushroom TYR, but not *hsTYR*, was not effective. This assay is a simple method to confirm *hsTYR* inhibition before conducting follow-up studies in live biological models.

Materials

Cell lines and materials

MM418C1 is a clone of the MM418 melanoma cell line that originated from the primary lesion of the cutaneous surface (Clark *et al.*, 1994 & Fechner *et al.*, 1994).

Reagents

L-Glutamine, penicillin/streptomycin, trypsin, RPMI 1640 medium (powder) were purchased from Life Technologies (CA, USA). Bovine serum albumin (BSA), 3,4-Dihydroxy-L-phenylalanine (L-Dopa), sodium phosphate, kojic acid, Triton X-100 and T75 flask were purchased from Sigma Chemicals (NSW, Australia). Fetal bovine serum (FBS) was purchased from Bovogen Biologicals (VIC, Australia). Trypan blue 0.4% solution was purchased from MP biomedical (VIC, Australia). The DC protein assay kit was purchased from Bio-Rad (NSW, Australia)

Safety warnings

- ⚠ Proper PPE must be worn all the time.
Risk assessment for handling chemicals and equipment use must be undertaken.



In vitro human tyrosinase inhibitory assay (human melanoma cell lysate based)

2h 25m

- 1 Seed human melanoma MM418C1 cells (or equivalent melanoma cell line) in a T75 cell culture flask and grow with RPMI 1640 medium, which contains 10% (v/v) fetal bovine serum (FBS), 200 mM L-glutamine, 100 U/mL penicillin, and 100 g/mL streptomycin (Sigma chemicals, Australia) until 80% confluent.  
- 2 Mix  10 μL of cell suspension with  10 μL of 0.4% (w/v) trypan blue solution, then apply  10 μL of the mixture onto a clean hemocytometer and count the cells. 10m 
- 3 Collect 3×10^6 cells (can be stored at  -80°C for up to at least  672:00:00 without loss of tyrosinase activity) by trypsinization and centrifugation at  300 x g, 25°C , 00:05:00 . 10m 
- 4 Lyse the cell pellet using  600 μL Triton X-100 lysis buffer (1% (v/v) Triton X-100, 50 mM sodium phosphate buffer,  6.8 , Sigma chemicals, Australia) and sonicate for  00:00:30  On ice using Epishear Sonicator (CA, USA). 5m 
- 5 Collect the supernatants by centrifugation at  7500 x g, 4°C , 00:30:00 . 30m 
- 6 Determine the protein concentration of the clarified cell lysate with the DC protein assay kit (Bio-Rad, Australia) 15m 
- 7 For each sample,  20 μL of the potential tyrosinase inhibitors (1 mg/mL) were added to a 96-well plate in triplicate, followed by  160 μL of 15 mM L-3,4-dihydroxyphenylalanine(L-Dopa) solution in 0.1 M sodium phosphate buffer (0.049 M $\text{Na}_2\text{HPO}_4 \cdot 7\text{H}_2\text{O}$, 0.051 M $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$). Finally,  20 μL of the clarified cell lysate was added to each well. 10m  
- 7.1  20 μL of 1 mg/mL kojic acid +  160 μL L-DOPA +  20 μL cell lysate is used as the positive control, 



20 µL of 1% (v/v) Triton X-100 lysis buffer + 160 µL L-DOPA + 20 µL cell

lysate is used as the blank.

- 8 Measure the change in absorbance at $\lambda 475$ nm over 01:00:00 using a CLARIOstar multi-mode plate reader (BMG Labtech, Germany) and normalize the data by using the protein concentration obtained in step 6.

1h



- 9 Tyrosinase activity (% remaining) can be calculated using the following equation:

5m



$$\text{Tyrosinase activity (\%)} = [(A1-A2) / (A3-A4)] * 100\%$$

L-DOPA solution and sodium phosphate buffer were used as a blank (A1)

L-DOPA and cell lysate without test sample were used as a negative control (A2)

L-DOPA and cell lysate with the test sample were the test reactions (A3)

L-DOPA and the test sample without cell lysate were used as a negative control (A4)

Protocol references

Clark, J., Grabs, A. J., Parsons, P. G., Smithers, B. M., Addison, R. S., & Roberts, M. S. (1994). Melphalan uptake, hyperthermic synergism and drug resistance in a human cell culture model for the isolated limb perfusion of melanoma. *Melanoma Research*, 4(6), 365-370.

Fechner, G. A., Michel, J., Sturm, R. A., Jacobs, J. J., & Parsons, P. G. (1994). Reduction of DNA synthesis, pigment synthesis, pigmentation gene mRNA and resistance to UVB in human melanoma cells treated with analogues of a histamine (H2) agonist. *Biochemical pharmacology*, 48(1), 121-130.

Acknowledgements

Griffith University Postgraduate Research Scholarship (GUPRS)

Griffith University International Postgraduate Research Scholarship (GUIPRS)

Professor Peter Parsons from QIMR Berghofer Research Institute (Brisbane, Australia) for providing the human melanoma cell line