



Aug 25, 2023

In vitro LRRK2 autophosphorylation

DOI

<https://dx.doi.org/10.17504/protocols.io.81wgb6m91lpk/v1>

Xinbo Wang^{1,2}, Pietro De Camilli^{1,2}

¹1. Departments of Neuroscience and of Cell Biology, Howard Hughes Medical Institute, Program in Cellular Neuroscience, Neurodegeneration and Repair, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

²2. Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815



xinbo.wang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgb6m91lpk/v1>

Protocol Citation: Xinbo Wang, Pietro De Camilli 2023. In vitro LRRK2 autophosphorylation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.81wgb6m91lpk/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: August 18, 2022

Last Modified: May 31, 2024

Protocol Integer ID: 68881

Keywords: LRRK2, Autophosphorylation, In vitro, ASAPCRN, Irrk2 autophosphorylation assay, Irrk2 autophosphorylation this protocol details method, Irrk2 autophosphorylation, Irrk2, assay

Abstract

This protocol details methods for the *in vitro* LRRK2 autophosphorylation assay.

Attachments



[iuubbv9p.docx](#)

18KB

Materials

✂ Prescission Protease **Genscript Catalog #Z02799**

Glutathione beads (GE Healthcare, 17075601),

✂ Amicon Ultra-15 Centrifugal Filter Unit **Merck MilliporeSigma (Sigma-Aldrich) Catalog #UFC901024 & UFC903024**

✂ Slide-A-Lyzer[®]; MINI Dialysis Device, 10K MWCO, 0.1 mL **Thermo Fisher Catalog #69572** ,

✂ Anti-LRRK2 (phospho T1357) antibody **Abcam Catalog #ab270606**

Solutions to prepare:

10x Kinase buffer:

	A	B
	Tris-HCl (pH7.5)	200 mM
	MgCl ₂	75 mM
	EGTA	1 mM

Dialysis buffer

	A	B
	HEPES (7.4)	20 mM
	NaCl	150 mM
	MgCl ₂	2.5 mM
	Glycerol	5%
	DTT	2 mM
	GDP	20 μM



In vitro LRRK2 autophosphorylation

- 1 Set up the reaction mixture in a 1.7 mL Eppendorf tube with  1.4 mL purified LRRK2 protein, 1x kinase buffer with  1 millimolar (mM) ATP and  0.01 U/ μ L GST-Prescission Protease (to remove the Flag tag).

Note

Note: the LRRK2 protein used in this experiment was obtained by elution from the anti-FLAG M2 resin as described in the LRRK2 purification protocol.

- 2 Incubate samples  Overnight at  4 °C .

- 3 Add Glutathione beads to remove GST-Prescission Protease.

- 4 Concentrate samples by centrifugal filters and dialyze  Overnight at  4 °C against dialysis buffer.

- 5 Check autophosphorylation by Western blotting using a LRRK2 phospho-specific (pT1357) antibody.

- 6 Determine protein concentration by SDS-PAGE using Bovine Serum Albumin (BSA) as standard and used without freezing in the liposome tubulation experiments.