

Sep 23, 2023

In vitro phosphatase assay

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov1qyepgr2/v1>

Elias Adriaenssens¹

¹Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna (AT)



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1qyepgr2/v1>

Protocol Citation: Elias Adriaenssens 2023. In vitro phosphatase assay. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.ewov1qyepgr2/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 27, 2023



Last Modified: May 31, 2024

Protocol Integer ID: 84093

Keywords: In vitro phosphatase assay, ASAPCRN, phosphatase assay, phosphatase, assay this protocol, assay

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-000350

Marie Skłodowska-Curie MSCA Postdoctoral fellowship

Grant ID: 101062916

Abstract

This protocol describes *in vitro* phosphatase assay.

Attachments



[756-1924.pdf](#)

49KB

Materials

Materials

- 6 well plates
- MnCl_2

 Pierce & Warriner; Detergent Compatible Bradford Assay Kit **Thermo Fisher Catalog #23246**

 Lambda Protein Phosphatase - 100,000 units **New England Biolabs Catalog #P0753L**

Lysis buffer

	A	B
	HEPES pH 7.4	50 mM
	NaCl	150 mM
	MgCl_2	2.5 mM
	DTT	2 mM
	NP-40	0.50%
	protease inhibitor cocktail	



In vitro phosphatase assay

5m

- 1 Seed HAP1 wild-type or FIP200 knockout cells in 6 well plates and grow until confluency.
- 2 Collect cells by trypsinization and pellet by centrifugation at  300 x g, 4°C, 00:05:00  .

- 3 After a PBS wash to remove the remaining cell medium, resuspend the cell pellets in lysis buffer. 
- 4 Lyse the samples for  00:20:00  On ice . Clear the cell lysates by centrifugation at  20000 x g, 4°C, 00:10:00 . 
 
- 5 Determine the protein concentrations of the cleared protein lysates with the Pierce Detergent Compatible Bradford Assay Kit (23246, Thermo Fisher).
- 6 For both samples, wild-type and FIP200 knockout lysates, incubate  100 µg of cell lysate with  5 µL of 10x NEBuffer for Protein MetalloPhosphatases (P0753, New England Biolabs) and  5 µL of  10 millimolar (mM) of MnCl₂ to make a total reaction volume of  50 µL . 
- 7 Add  1 µL of Lambda Protein Phosphatase (NEB) to the reaction and incubate the samples at  30 °C for the indicated time.  
- 8 Terminate the phosphatase reactions by the addition of 6x Protein Loading dye and heat inactivation at  95 °C for  00:05:00 . 

- 9 Analyze the samples by western blot analysis as described above. 

Protocol references

<https://international.neb.com/protocols/2017/12/28/typical-reaction-protocol-for-lambda-protein-phosphatase>