



Mar 21, 2025

In Vitro Phosphatase Assay

 [Nature Cell Biology](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvm9b95l3p/v1>

Elias Adriaenssens¹

¹Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvm9b95l3p/v1>

External link: <https://doi.org/10.1038/s41556-025-01712-y>

Protocol Citation: Elias Adriaenssens 2025. In Vitro Phosphatase Assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.yxmvm9b95l3p/v1>

Manuscript citation:

Adriaenssens, E., Schaar, S., Cook, A.S.I. *et al.* Reconstitution of BNIP3/NIX-mitophagy initiation reveals hierarchical flexibility of the autophagy machinery. *Nat Cell Biol* **27**, 1272–1287 (2025). <https://doi.org/10.1038/s41556-025-01712-y>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 10, 2025

Last Modified: March 21, 2025

Protocol Integer ID: 120012

Keywords: ASAPCRN, vitro phosphatase assay this protocol detail, vitro phosphatase assay, phosphatase assay, assay

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-000350

Marie Skłodowska-Curie MSCA Postdoctoral fellowship

Grant ID: 101062916

Abstract

This protocol details the in-vitro phosphatase assay.

Materials

Lysis Buffer:

	A	B
	Tris- HCl, pH 7.4	20 mM
	NaCl	150 mM
	EDTA	1 mM
	EGTA	1 mM
	Triton X-100	1%
	Protease-inhibitor cocktail, but no phosphatase inhibitors	

1x phosphatase buffer:

	A	B
	Tris- HCl, pH 7.4	50 mM
	NaCl	100 mM
	DTT	2 mM
	MnCl ₂	1 mM

-  Pierce & Warriner; Detergent Compatible Bradford Assay Kit **Thermo Fisher Catalog #23246**



Procedure

1h 5m

- 1 To verify the activity of our recombinantly purified λ phosphatase, collect HeLa cells by trypsinisation and wash the cell pellet with PBS once before cells were lysed in lysis buffer.



Lysis buffer:

	A	B
	Tris- HCl, pH 7.4	20 mM
	NaCl	150 mM
	EDTA	1 mM
	EGTA	1 mM
	Triton X-100	1%
	Protease-inhibitor cocktail, but no phosphatase inhibitors	

- 2 Lyse the cells for 00:20:00 On ice before cell lysates were cleared by centrifugation at 20000 x g, 4°C, 00:10:00 .

30m



- 3 Determine the protein concentrations of the cleared protein lysates with the Pierce Detergent Compatible Bradford Assay Kit (23246, Thermo Fisher) and incubate the equal amounts with or without home-made λ phosphatase in 1x phosphatase buffer diluted from a 10x stock.



1x phosphatase buffer:

	A	B
	Tris- HCl, pH 7.4	50 mM
	NaCl	100 mM
	DTT	2 mM
	MnCl ₂	1 mM



- 4 For this, protein samples were, where necessary, supplement with lysis buffer to a total volume of  45 μL and supplement with  5 μL of 10x phosphatase buffer.
- 5 The total amount of protein lysate depends on the abundance of the target protein you want to verify with western blotting. However, for us, this typically ranges between  50 μg -  500 μg in total.
- 6 Incubate the samples for  00:30:00 at  30 $^{\circ}\text{C}$ before reactions were terminated by the addition of protein loading dye and heat inactivation for  00:05:00 at  95 $^{\circ}\text{C}$.
- 7 Analyse the samples by SDS-PAGE and western blotting.

35m



8

Note

Note, some detergents can inhibit the phosphatase reaction, in that case opening the cells with a 26G needle can serve as a substitute in following buffer (50 mM HEPES pH 7.4, 100 mM NaCl, 2 mM DTT).