

Aug 03, 2025

In-vitro Thermodynamic Solubility

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygxqdmzv8j/v1>

Nick Lynch^{1,2}

¹Curlew Research; ²ASAP Discovery Consortium

ASAP Discovery



Nick Lynch

Curlew Research, ASAP Discovery Consortium

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygxqdmzv8j/v1>

Protocol Citation: Nick Lynch 2025. In-vitro Thermodynamic Solubility. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kxygxqdmzv8j/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 25, 2025

Last Modified: August 03, 2025

Protocol Integer ID: 218908

Keywords: ADME, DMPK, drug discovery, solubility, physical chemical, thermodynamic, thermodynamic solubility assay, thermodynamic solubility, saturation solubility, true solubility, dissolved compound, given solvent, formulation, solvent under specific condition, crucial assay, assay, vitro

Disclaimer

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgement, advice, diagnosis, or treatment. Any action you take or refrain from taking, using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held

Abstract

A thermodynamic solubility assay measures the maximum amount of a compound that can dissolve in a given solvent under specific conditions at equilibrium. It's a crucial assay for understanding a compound's true solubility, particularly in formulations and early drug development. This method helps determine the saturation solubility, where the dissolved compound is in equilibrium with the undissolved solid.

Safety warnings

- ! Always wear appropriate PPE for this protocol
Refer to Material Safety Data Sheets for additional safety and handling information.

Summary

- 1 Thermodynamic solubility: Dry solid compound are used, where particle size and crystal morphology influence the amount solubilized.

Results provide important information for PK formulations, pH adjustment and cosolvent being two of the most used methods to solubilize poorly soluble drug compounds during later stages of development.

Sample Preparation

- 2 The thermodynamic solubility assay is performed using 1 mg of compound, accurately weighed into a 1.5 mL glass vial
- 3 To each sample, 1 mL of 0.1 M phosphate buffer (pH 7.4) is added
- 4 The vials are then incubated in a Thermomixer at 21°C and 700 rpm for 24 hours to allow the compound to reach equilibrium solubility

Assay

- 5 Following incubation, the supernatant is separated from undissolved solids via a two-step centrifugation process
- 6 After the second centrifugation, 5 µL of the clear supernatant is diluted with 495 µL of 30% methanol/70% water and mixed thoroughly to produce "Solution X."
- 7 Both the remaining supernatant and Solution X are used for analytical preparation.
- 7.1 For LC-MS/MS analysis, samples are prepared directly into a 96-well plate.
- 7.2 Two analytical solutions are made for each sample. Solution A consists of 5 µL of Solution X, 395 µL of 30% methanol/70% water, and 100 µL of internal standard (IS).



- 8 Solution B is prepared by mixing 25 μL of the original supernatant with 375 μL of 30% methanol/70% water and 100 μL IS.
- 9 A standard curve is also generated using four concentrations:
Standard 1 (1.25 nM),
Standard 2 (12.5 nM),
Standard 3 (125 nM), and
Standard 4 (1250 nM).
- 10 All samples and standards are injected onto the LC-MS/MS system, and the compound concentrations in the samples are quantified by comparing their responses to those of the known standards.

Data Analysis

- 11 For each test compound injection:

Response ratio = Test peak area / Internal standard peak area

Test compound response ratio is converted to concentration (μM) from the standard curve.

Protocol references

https://cdn.prod.website-files.com/659db4dde7f13f395f49157c/65f05d2d5345728df57c4cbf_Concept%20Life%20Sciences_LE_SOLUBILITY%20_2_.pdf

Acknowledgements

Grateful to Concept Life Sciences for supplying the original protocol summary