



Aug 05, 2024

In vivo microdialysis for striatal DA release

 In 2 collections

DOI

dx.doi.org/10.17504/protocols.io.n92ld8q1xv5b/v1

Ariadna Laguna¹, Miquel Vila¹

¹VHIR-CIBERNED-ASAP

Vilalab Public

Nuria



Miquel Vila

VHIR-CIBERNED-ASAP

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ld8q1xv5b/v1>

Protocol Citation: Ariadna Laguna, Miquel Vila 2024. In vivo microdialysis for striatal DA release. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n92ld8q1xv5b/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 05, 2024

Last Modified: August 05, 2024

Protocol Integer ID: 104683

Keywords: vivo microdialysis for striatal da release, vivo microdialysi, local effects of nomifensine, nomifensine, net inhibitor

Disclaimer

The **protocols.io** team notes that research involving animals and humans must be conducted according to internationally-accepted standards and should always have prior approval from an Institutional Ethics Committee or Board.

Abstract

To assess local effects of nomifensine (DAT and NET inhibitor) on striatal DA release in microdialysis experiments

Troubleshooting

Preparation of reagents

- 1 Nomifensine (DAT and NET inhibitor) is dissolved in artificial cerebrospinal fluid (aCSF in mM: NaCl, 125; KCl, 2.5; CaCl₂, 1.26 and MgCl₂ 1.18)
- 2 Concentrated solutions (1 mM; pH adjusted to 6.5–7 with NaHCO₃ when necessary) are stored at -80°C and working solutions are prepared daily by dilution in aCSF.

Mouse surgery

- 3 One concentric dialysis probe (Cuprophane membrane; 6000 Da molecular weight cut-off; 1.5 mm-long) is implanted in the striatum (AP, 0.5; ML, -1.7; DV, -4.5 in mm) of isoflurane-anesthetized mice
- 4 Microdialysis experiments are performed in freely-moving mice 24h after surgery.

Microdialysis

- 5 Probes are perfused with aCSF at 1.5 µL/min
- 6 Following an initial 100-min stabilization period, 5 or 7 baseline samples are collected (20 min each) before local drug application
- 7 Nomifensine is administered by reverse dialysis at 10 and 50µM (uncorrected for membrane recovery) and then successive dialysate samples are collected

Dopamine levels determination

- 8 The concentration of DA in dialysate samples is determined by HPLC coupled to electrochemical detection (+0.7 V, Waters 2465), with 3-fmol detection limit.
- 9 The mobile phase containing 0.15 M NaH₂PO₄·H₂O, 0.9 mM PICB8, 0.5 mM EDTA (pH 2.8 adjusted with orthophosphoric acid), and 10 % methanol is pumped at 1 ml/min (Waters 515 HPLC pump)
- 10 DA is separated on a 2.6 mm particle size C18 column (7.5 × 0.46 cm, Kinetex, Phenomenex) at 28°C.



Data representation

- 11 Microdialysis data is expressed as femtomoles per fraction (uncorrected for recovery) and are shown as percentages of basal values (individual means of 5-7 pre-drug fractions).