

Jun 06, 2023

iNDI Papain Dissociation protocol for iNeurons Version 1

DOI

dx.doi.org/10.17504/protocols.io.ewov1oddylr2/v1

Erika Lara Flores¹, Joel Reyes², Mark Cookson¹, Michael Ward¹

¹National Institutes of Health, NIH Center for Alzheimer's and Related Dementias (CARD);

²National Institutes of Health, National Institute of Neurological Disorders and Stroke (NINDS)

Neurodegeneration Met...

iNDI Protocol Developm...



Erika Lara Flores

NIH/NIA

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1oddylr2/v1>

Protocol Citation: Erika Lara Flores, Joel Reyes, Mark Cookson, Michael Ward 2023. iNDI Papain Dissociation protocol for iNeurons Version 1. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.ewov1oddylr2/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: December 13, 2022

Last Modified: June 06, 2023

Protocol Integer ID: 73921

Keywords: papain dissociation protocol for neurons version, papain dissociation protocol, proteolytic enzyme, cell dissociation, neurons for scRNAseq, papain, neurons version, other protease, effective than other protease, derived neuron, neuronal culture

Abstract

Proteolytic enzymes are widely used in cell dissociation and papain has proved less damaging with some tissues and neuronal cultures and more effective than other proteases. Here we describe a papain dissociation protocol that has been used on our iPSC-derived neurons for scRNAseq and FACS assays readouts.

Materials

⊗ Papain Dissociation System **Worthington Biochemical Corporation Catalog #LK003153**

⊗ TrypLE™ Express Enzyme (1X) no phenol red **Thermo Fisher Scientific Catalog #12604013**

PBS

PBS-EDTA 0.05mM

PBS-0.04% BSA (For scRNAseq cell preparation, if applicable)

Troubleshooting



Preparation of Reagents

- 1 **Equilibrate EBSS/phenol red** with 95%O₂:5%CO₂ in incubator for several hours of overnight.
- 2 **Reconstitute the Ovomucoid mixture (Papain inhibitor-PI)** by adding  32 mL of EBSS and allow contents to dissolve. This will yield solution at an effective concentration of 10mg of ovomucoid inhibitor and 10mg of albumin per mL. Reconstitute for the first use, then store  1 mL at  4 °C .
- 3 **Papain solution:** Add  10 mL of TrypLE™ Express Enzyme (1X), no phenol red to the vial and place it in  37 °C bath for 10 minutes, this vial has 10 units/mL.

Note

When working with i3Neurons (WTC11 line), dissolve the papain with 20 mL of TrypLE™ Express Enzyme (1X).

- 4 **DNase I vial D2:** Add 500ul of EBSS media and mix gently (2000 units/mL).
- 5 **Trituration Medium:**
 30 mL of Neuronal Maturation medium (the same medium used for your differentiation)
 30 µL of Rock inhibitor (1x)
1 vial of DNase I
- 6 **Papain Inhibitor media:**  9 mL of Trituration medium +  1 mL of papain inhibitor.

Dissociation protocol for iNeurons at day 28 of differentiation

5m

- 7 Aspirate culture medium and wash gently with PBS (calcium/magnesium free).



- 8 Aspirate PBS and wash 2 times with PBS-0.5mM EDTA
- 9 Aspirate PBS-EDTA and add to one well of 6 well plate  1 mL of Papain solution (10 units/ml) to cells, and incubate at  37 °C for  00:05:00 .

5m

Note

Incubation time can vary and go up to 30 minutes and it depends on cell line, confluency and maturation day.

Incubation times for i3Neurons (WTC11):

Day 7: 3 minutes

Day 10: 5 minutes

Day 14: 10 minutes

- 10 Aspirate papain solution when the cell bodies have an halo morphology like for and EDTA split. If the incubation is longer (up to 30 minutes) the cells will detach into papain solution, then continue with the next step.
- 11 Add  1 mL of trituration medium per well of a 6-well plate to dissociate cells and pipette medium over the plate up and down to detach cells and dissociate them as single-cell suspension.
- 12 Collect cells into a sterile conical 15 mL tube and adjust volume to 5 mL with PBS or trituration medium.
- 13 Centrifuge cells  00:05:00 at 200 - 300 x g at  Room temperature .
- 14 Remove supernatant and resuspend the cell pellet gently with  1 mL of Papain inhibitor medium to wash off the papain.
- 15 Centrifuge cells  00:05:00 at 200 - 300 x g at  Room temperature
- 16 Aspirate supernatant and resuspend cell pellet in  1 mL of Neuronal Maturation Medium.

5m

5m



- 17 After resuspending the cell pellet, count cells and use as necessary.

Note

For scRNAseq 10X genomics assay resuspend in cold PBS-0.04%BSA.
For FACS assay resuspend in cold PBS-0.5mM EDTA.