

Aug 02, 2023

Indirect Proximity Ligation Assay (PLA) - Brightfield

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l6x94klqe/v1>

Ayse Ulusoy¹, Rita Pinto-Costa¹, Angela Rollar¹, Donato Di Monte¹

¹DZNE



Ayse Ulusoy

German Center for Neurodegenerative Diseases

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l6x94klqe/v1>

Protocol Citation: Ayse Ulusoy, Rita Pinto-Costa, Angela Rollar, Donato Di Monte 2023. Indirect Proximity Ligation Assay (PLA) - Brightfield. protocols.io <https://dx.doi.org/10.17504/protocols.io.bp2l6x94klqe/v1>

Manuscript citation:

doi: 10.1126/sciadv.abn0356

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: August 02, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 85848

Keywords: post-translational modification, alpha-synuclein, oxidative stress, ASAPCRN, brightfield indirect proximity ligation assay, indirect proximity ligation assay, nitration of mitochondrial enzyme, nitration of protein, mitochondrial enzyme, synuclein, protein modification, dependent ligation, protein, protein interaction, subunit ndufb8, protein complex formations within cell, nitrated alpha, powerful molecular technique, specific antibody, pla protocol

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000420

Abstract

Indirect Proximity Ligation Assay (PLA) is a powerful molecular technique used to detect and visualize protein-protein interactions, protein modifications, and protein complex formations within cells or tissues. This method is based on the principles of proximity-dependent ligation and utilizes specific antibodies to detect the nitration of proteins on free-floating brain sections. Here we describe the PLA protocol that we routinely use in our laboratory to detect nitrated alpha-synuclein and nitration of mitochondrial enzymes such as SOD2 and the mitochondrial complex 1 subunit NDUF8.



Materials

Duolink InSitu PLA probe anti-rabbit PLUS kit: DUO92002-100RXN

Duolink InSitu PLA probe anti-mouse MINUS kit: DUO92004-100RXN

Duolink InSitu Detection Reagents Brightfield kit: DUO92012-100RXN

Duolink wash buffer A: DUO82047-20L

Histomount mounting medium: Life Technologies 008030

Antibodies:

mouse anti-3-NT: 1:250; ab61392, Abcam

rabbit anti-human alpha-synuclein (clone MJFR1): 1:4000, ab138501, Abcam

rabbit anti-SOD1: 1:1000; ADI-SOD-110, Enzo Life Sciences

rabbit anti-NDUFB8: 1:300; 14794, Proteintech



Day 1

1

Pick 35um cut brain sections and transfer them to  1.5 mL Eppendorf tubes

Note: all incubation and wash steps are performed by shaking Eppendorf tubes at 250rpm (e.g., thermomixer)

2

Wash 2x  00:05:00 with Tris-HCl

5m

3

Peroxidase quenching with BloxAll solution

use  300 µL and incubate for 30min

4

Wash 3x  00:05:00 with wash buffer A (see materials)

5m

5

Antigen-retrieval with citrate buffer (+tween, pH = 6) at 95°C for  00:04:00

4m

Note: Heat the solution to 95°C before adding on the samples

6

Wash 3x  00:05:00 with wash buffer A

5m

7

Blocking: Incubate in Duolink Blocking Solution ( 300 µL) at 37°C for  01:00:00

1h

8

Primary Antibody Incubation: Dilute antibodies in Duolink Antibody Diluent ( 200 µL

µL/tube) and incubate  Overnight room temperature.

for nitrated alpha-synuclein: mouse anti-3-NT and rabbit anti-h-alpha-synuclein

for nitrated SOD2: mouse anti-3-NT and rabbit anti-SOD2

for nitrated NDUFB8: mouse anti-3-NT and rabbit anti-NDUFB8

see materials for dilutions and catalog numbers

5m



Day 2

10m

- 9 Wash 3x  00:10:00 with wash buffer A 10m
- 10 **PLA probe incubation:** 1h 30m
Prepare anti-mouse MINUS and anti-rabbit PLUS probes following manufacturer's instructions, add solution on sections and incubate at 37°C for  01:30:00
- 11 Wash 3x  00:10:00 with wash buffer A 10m
- 12 **Ligation:** 1h 15m
Dilute the Duolink Ligation Buffer 1:5 in high-purity water
Add the ligase (diluted 1:40) just before incubation (keep cold on freezer block!)
add the ligation solution on sections and incubate at 37°C for  01:15:00
- 13 Wash 3x  00:05:00 with wash buffer A 5m
- 14 **Amplification:** 2h 15m
Dilute the Duolink Amplification Buffer 1:5 in high-purity water
Add the polymerase (diluted 1:80) just before incubation (keep cold on freezer block!)
Add the amplification solution to the sections and incubate at 37°C for  02:15:00
- 15 **Detection:** 1h 10m
Dilute the Duolink Detection Brightfield Stock 1:5 in high-purity water
add solution to the sections and incubate at RT for  01:10:00
- 16 Wash 2x  00:05:00 with wash buffer A 10m



Wash 1x  00:05:00 with Tris-buffered saline (no detergent)

17 **Developing:**

Dilute the Duolink Brightfield Substrate solutions in high-purity water:

§ Substrate A (1:70)

§ Substrate B (1:100)

§ Substrate C (1:100)

§ Substrate D (1:50)

Incubate sections at RT for 1-5min (depending on desired intensity outcome)

18 Wash 3x  00:05:00 with wash Tris-buffered saline

5m

19 Mount the samples on slides and let them air-dry

20 **Nuclear staining** (optional): Incubate for  00:01:00

11m

Wash the staining off with running tap water for  00:10:00

21 Dehydrate of samples in 1min each: 70% ethanol, 95% ethanol, 2× 100% ethanol 2x xylene

22 Coverslip samples with Histomount mounting medium

Protocol references

doi: 10.1126/sciadv.abn0356