



Aug 02, 2023

## Indirect Proximity Ligation Assay (PLA) - Fluorescence

DOI

[dx.doi.org/10.17504/protocols.io.261ged36dv47/v1](https://dx.doi.org/10.17504/protocols.io.261ged36dv47/v1)

Ayse Ulusoy<sup>1</sup>, Rita Pinto-Costa<sup>1</sup>, Angela Rollar<sup>1</sup>, Donato Di Monte<sup>1</sup>

<sup>1</sup>DZNE



Ayse Ulusoy

German Center for Neurodegenerative Diseases

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.261ged36dv47/v1>

**Protocol Citation:** Ayse Ulusoy, Rita Pinto-Costa, Angela Rollar, Donato Di Monte 2023. Indirect Proximity Ligation Assay (PLA) - Fluorescence. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.261ged36dv47/v1>

**Manuscript citation:**

doi: 10.1126/sciadv.abn0356

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 02, 2023

**Last Modified:** May 31, 2024

**Protocol Integer ID:** 85854

**Keywords:** post-translational modification, alpha-synuclein, oxidative stress, ASAPCRN, fluorescence indirect proximity ligation assay, indirect proximity ligation assay, nitration of mitochondrial enzyme, nitration of protein, synuclein, mitochondrial enzyme, protein, subunit ndufb8, protein complex formations within cell, protein modification, protein interaction, dependent ligation, nitrated alpha, powerful molecular technique, specific antibody

**Funders Acknowledgements:**

**ASAP**

Grant ID: ASAP-000420

## Abstract

Indirect Proximity Ligation Assay (PLA) is a powerful molecular technique used to detect and visualize protein-protein interactions, protein modifications, and protein complex formations within cells or tissues. This method is based on the principles of proximity-dependent ligation and utilizes specific antibodies to detect the nitration of proteins on free-floating brain sections. Here we describe the PLA protocol that we routinely use in our laboratory to detect nitrated alpha-synuclein and nitration of mitochondrial enzymes such as SOD2 and the mitochondrial complex 1 subunit NDUF8.



## Materials

Blocking solution: Vector Lab SP-6000-100

Duolink In Situ Wash Buffers, Fluorescence: DUO82049-20L

Duolink InSitu PLA probe anti-rabbit PLUS kit: DUO92002-100RXN

Duolink InSitu PLA probe anti-mouse MINUS kit: DUO92004-100RXN

Duolink In Situ Detection Reagents Red: DUO92008-100RXN

Duolink *In Situ* Mounting Medium with DAPI: DUO82040-5ML

### **Antibodies:**

mouse anti-3-NT: 1:250; ab61392, Abcam

rabbit anti-human alpha-synuclein (clone MJFR1): 1:4000, ab138501, Abcam



rabbit anti-SOD1: 1:1000; ADI-SOD-110, Enzo Life Sciences

rabbit anti-NDUFB8: 1:300; 14794, Proteintech

## Troubleshooting



## Day 1

- 1 Pick 35um cut brain sections and transfer them to  1.5 mL Eppendorf tubes:  
Note: all incubation and wash steps are performed by shaking Eppendorf tubes at 250rpm (e.g., thermomixer)
- 2 Wash 2x  00:05:00 with Tris-HCl 5m
- 3 Wash 3x  00:05:00 with wash buffer A (see materials) 5m
- 4 **Antigen-retrieval** with citrate buffer (+tween, pH = 6) at 95°C for  00:04:00 4m  
Note: Heat the solution to 95°C before adding on the samples
- 5 Wash 3x  00:05:00 with wash buffer A 5m
- 6 **Blocking:** Incubate in Duolink Blocking Solution (  300 µL ) at 37°C for  01:00:00 1h
- 7 **Primary Antibody Incubation:** Dilute antibodies in Duolink Antibody Diluent (  200 µL ) and incubate  Overnight room temperature. 5m  
for nitrated alpha-synuclein: mouse anti-3-NT and rabbit anti-h-alpha-synuclein  
for nitrated SOD2: mouse anti-3-NT and rabbit anti-SOD2  
for nitrated NDUFB8: mouse anti-3-NT and rabbit anti-NDUFB8  
**see materials for dilutions and catalog numbers**

## Day 2

- 8 Wash 3x  00:10:00 with wash buffer A 10m



- 9 **PLA probe incubation:** 1h 30m

Prepare anti-mouse MINUS and anti-rabbit PLUS probes following manufacturer's instructions, add solution on sections and incubate at 37°C for 01:30:00
- 10 Wash 3x 00:10:00 with wash buffer A 10m
- 11 **Ligation:** 1h 15m

Dilute the Duolink Ligation Buffer 1:5 in high-purity water  
Add the ligase (diluted 1:40) just before incubation (keep cold on freezer block!)  
add the ligation solution on sections and incubate at 37°C for 01:15:00
- 12 Wash 3x 00:05:00 with wash buffer A 5m
- 13 **Amplification:** This step is light-sensitive! incubation is performed using aluminum foil to protect samples from light. 2h 30m

Dilute the Duolink Amplification Buffer (red) 1:5 in high-purity water  
Add the polymerase (diluted 1:80) just before incubation (keep cold on freezer block!)  
Add the amplification solution to the sections and incubate at 37°C for 02:30:00
- 14 Wash 2x 00:05:00 with wash buffer B 6m

Wash 1x 00:01:00 with wash 0.1x buffer B
- 15 **Additional immunohistochemical staining (optional):** if immunohistochemical counter-staining is going to be performed it is important to avoid the use of detergents in wash buffers and perform long incubations in cold.
- 16 Mount the samples on slides and let them air-dry for 5 mins
- 17 Coverslip samples with Duolink *In Situ* Mounting Medium with DAPI, store samples in the dark at 4°C

## Protocol references

doi: 10.1126/sciadv.abn0356