

Mar 19, 2025

Induced Pluripotent Stem Cells (iPSCs) Differentiation into iNeurons

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2lyxz3wgx9/v1>

Patricia Yuste-Checa¹, F Ulrich Hartl¹

¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry, Martinsried, Germany



Patricia Yuste-Checa

MPI Biochemistry

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2lyxz3wgx9/v1>

Protocol Citation: Patricia Yuste-Checa, F Ulrich Hartl 2025. Induced Pluripotent Stem Cells (iPSCs) Differentiation into iNeurons. [protocols.io https://dx.doi.org/10.17504/protocols.io.eq2lyxz3wgx9/v1](https://dx.doi.org/10.17504/protocols.io.eq2lyxz3wgx9/v1)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 11, 2025

Last Modified: March 19, 2025

Protocol Integer ID: 124471

Keywords: ASAPCRN, stemdiff forebrain neuron differentiation kit, stemdiff forebrain neuron maturation kit, induced pluripotent stem cell, stemdiff smadi neural induction kit, pluripotent stem cell, differentiation into ineuron, stem cell technology, type neuron, neuron, ineuron, ipsc

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Abstract

This protocol details how to efficiently differentiate induced pluripotent stem cells (iPSCs) into iPSC-derived forebrain-type neurons using the commercial kits, STEMdiff SMADi Neural Induction Kit (08581, Stem Cell Technologies, monolayer protocol), STEMdiff Forebrain Neuron Differentiation Kit (08600, Stem Cell Technologies) and STEMdiff Forebrain Neuron Maturation Kit (08605, Stem Cell Technologies).

Materials

- Matrigel (Corning) or Geltrex (Thermo Fisher Scientific)
- ROCK inhibitor (Y27632)
- iNeurons in Neurobasal Plus medium (Thermo Fisher Scientific) supplemented with:
 1. B-27 Plus Supplement 1x (Thermo Fisher Scientific)
 2. [M] 0.5 millimolar (mM) GlutaMAX (Thermo Fisher Scientific)
 3. 100 U/ml penicillin,  100 μ L streptomycin sulfate (Thermo Fisher Scientific).
- Anti-MAP2 antibody (AB554, Merck)
-  beta-3 Tubulin Monoclonal Antibody (TU-20) Thermo Fisher Scientific Catalog #MA1-19187
-  STEMdiff™ SMADi Neural Induction Kit 1 Kit STEMCELL Technologies Inc. Catalog #8581
-  STEMdiff™ Forebrain Neuron Differentiation Kit STEMCELL Technologies Inc. Catalog #08600
-  STEMdiff™ Forebrain Neuron Maturation Kit STEMCELL Technologies Inc. Catalog #08605



iPSC Culture

- 1 Coat 6-well plates with Matrigel or Geltrex by dilution 1/100 in F-12 DMEM medium (for example by diluting  100 μ L in 10 mL F-12 DMEM). Add  1 mL /well (6-well plate) and incubate at least 1 hour in the CO₂ incubator ( 37 °C).



Note

After the incubation time, you can keep the plates at  4 °C for several weeks. Make sure the liquid does not evaporate by adding more F-12 DMEM after coating was completed.

- 2 Thaw quickly a cryotube containing iPSCs in a water bath at  37 °C . Add the cells to 5 mL F-12 DMEM and centrifuge at  300 x g, 00:04:00 .
- 3 Remove the supernatant, resuspend the pellet in mTeSR or mTeSR plus medium supplemented with  10 micromolar (μ M) ROCK inhibitor with a Pasteur plastic pipette to avoid formation of single cells, place cell suspension in the Matrigel-coated well, move the plate softly in several back-and-forth and side-to-side motions to distribute the cells evenly and incubate in the CO₂ incubator ( 37 °C , 5% CO₂).

4m



Note

- iPSCs grow in colonies and are routinely passaged as small cell aggregates or clumps. Generation of single cells is stressful for the cells and can lead to genetic aberrations. Therefore, a Pasteur plastic pipette or a 1000 μ L pipette (wide openings, if using 1000 μ L pipette, do not pipette the cells more than 4-5 times) should be used when routinely passaging iPSCs to prevent production of single cells.
- ROCK inhibitor (Y27632) is an apoptosis inhibitor that should be added to the cells when the cells are likely to be stressed (e.g. when thawing or splitting with Accutase to generate single cells, see below).

- 4 Do full-medium change according to medium manufacturer's instructions (mTeSR every day, mTeSR plus every other day).



- 5 Once confluent, split the cells using ReLeSR (Stem Cell Technologies) or similar (PBS-EDTA) to obtain cellular aggregates and avoid making single cells for maintenance.

Generation of Neural Progenitor Cells (NPCs, Monolayer Protocol)

- 6 Remove the medium from the wells.
- 7 Add  1 mL Accutase to make single cells and incubate for  00:05:00 in the CO₂ incubator ( 37 °C , 5% CO₂).    5m
- 8 Add  1 mL F-12 DMEM supplemented with  10 micromolar (μM) Y-27632 (ROCK inhibitor) (no need of ROCK inhibitor if you are quick) and collect the cells with the 1000 μL pipette. 
- 9 Count viable cells using Trypan blue staining.
- 10 Centrifuge the cells at  300 x g, 00:04:00 .  4m
- 11 Aspirate the supernatant and resuspend the pellet with STEMdiff Neural Induction medium supplemented with SMADi and  10 micromolar (μM) Y-27632 (ROCK inhibitor) to a final concentration of 1×10⁶ cells/mL. 
- 12 Add  2 mL of cell suspension per Matrigel-coated well of a 6-well plate (2×10⁶ cells/well), move the plate softly in several back-and-forth and side-to-side motions to distribute the cells evenly and incubate in the CO₂ incubator ( 37 °C , 5% CO₂).   
- 13 Do a full-medium change with STEMdiff Neural Induction Medium supplemented with SMADi every day for the next 6 days.
- 14 On day 7 (Day in vitro (DIV) 7), remove the medium from the well, add  1 mL Accutase and incubate for  00:05:00 in the CO₂ incubator ( 37 °C , 5% CO₂).    5m



- 15 Add  1 mL F-12 DMEM supplemented with  10 micromolar (μM) Y-27632 (ROCK inhibitor) (no need of ROCK inhibitor if you are quick) and collect the cells with the 1000 μL pipette. 
- 16 Centrifuge the cells at  300 x g, 00:04:00 . 

- 17 Aspirate the supernatant and resuspend the pellet with  2 mL STEMdiff Neural Induction Medium supplemented with SMADi and  10 micromolar (μM) Y-27632 (ROCK inhibitor). 
- 18 Count viable cells using Trypan blue staining.
- 19 Dilute cells to a final concentration of 1×10^6 cells/mL with STEMdiff Neural Induction Medium supplemented with SMADi and 10 μM Y-27632 (ROCK inhibitor).
- 20 Add  2 mL of cell suspension per Matrigel-coated well of a 6-well plate (2×10^6 cells/well), move the plate softly in several back-and-forth and side-to-side motions to distribute the cells evenly and incubate in the CO₂ incubator ( 37 °C , 5% CO₂).   
- 21 Do a full-medium change with STEMdiff Neural Induction Medium supplemented with SMADi every day for the next 6 days.
- 22 Repeat splitting on DIV 13 (steps 14-20).
- 23 On DIV 20, split cells as previously described (steps 14-18).
- 24 Centrifuge cells and freeze neural progenitor cells (NPCs) using commercial freezing media or STEMdiff Neural Induction Medium supplemented with SMADi and 10% DMSO ($2\text{--}4 \times 10^6$ cells/cryotube). 

**Note**

On DIV 20 split, you can also plate the cells directly to continue the differentiation as described in the next section.

Generation of iPSC Derived-Forebrain-Type Neurons from NPCs

- 25 Coat Poly-L-Ornithine/Laminin (PO/L) wells. Coat wells of a 6-well plate with  500 μ L poly-Ornithine solution (0.01% in PBS) during  04:00:00 or more hours at  Room temperature . Aspirate solution, then add  500 μ L Laminin solution (10-20 μ g/mL Laminin in PBS) per well of a 6-well plate during 4 hours to  Overnight (better) at  Room temperature .

4h

**Note**

1. To coat a 12-well plate use  300 μ L -  350 μ L PO and L and for a 24-well plate use  200 μ L -  250 μ L PO and L.
2. You can keep the coated plates at  4 $^{\circ}$ C for several weeks. Make sure the liquid does not evaporate by adding more PBS after coating was completed.

- 26 Thaw NPCs quickly in a water bath, add cells to 5 mL of F-12 DMEM media and centrifuge at  300 x g, 00:04:00 .



- 27 Discard supernatant, resuspend the pellet with STEMdiff Neural Induction Medium supplemented with SMADi to a final concentration of 1×10^6 cells/mL.

- 28 Add  2 mL of cell suspension per PO/L-coated well of a 6-well plate (2×10^6 cells/well), move the plate softly in several back-and-forth and side-to-side motions to distribute the cells evenly and incubate in the CO₂ incubator ( 37 $^{\circ}$ C , 5% CO₂).



- 29 Do full-medium change with STEMdiff Forebrain Neuron Differentiation Medium for the next 5 days.



30 On day 6 perform the final plating of the cells. Remove medium from the well, add  1 mL Accutase and incubate for  00:05:00 in the incubator ( 37 °C , 5% CO₂).

5m



31 Add  1 mL F-12 DMEM and collect the cells with the 1000 µL pipette.



32 Centrifuge cell at  300 x g, 00:04:00 .

4m



33 Aspirate the supernatant, resuspend the pellet with 4 mL STEMdiff Forebrain Neuron Maturation Medium and pass the cells through a strainer to make single cells.

34 Count viable cells using Trypan blue staining.

35 Dilute cells to the desired final concentration with STEMdiff Forebrain Neuron Maturation Medium, dispense them in PO/L-coated wells, move the plate softly in several back-and-forth and side-to-side motions to distribute the cells evenly and incubate in the CO₂ incubator ( 37 °C , 5% CO₂).



Note

Dispense 5×10^5 cells per well of a 6-well plate, 2.5×10^5 cells per well of a 12-well plate and 1×10^5 cells per well of a 24-well plate.

36 Do half-medium change every other day with STEMdiff Forebrain Neuron Maturation Medium supplemented with  5 micromolar (µM) 5-Fluorouracil/Uridine (Merck) to stop growth of non-differentiated cells.



37 After 8 days in STEMdiff Forebrain Neuron Maturation Medium, maintain iNeurons in Neurobasal Plus medium supplemented with B-27 Plus Supplement 1x, 0.5 mM GlutaMAX and 100 U/ml penicillin,  100 µL streptomycin sulfate. Perform half-medium change every 2-3 days.



Note

Neuronal differentiation may be assessed by immunofluorescence using the anti-MAP2 antibody, 1/500 dilution and the anti- β -3-tubulin, 1/100 dilution or similar.