



Jan 22, 2024

Induction of aggregation in alpha-synuclein-expressing cells by treatment with preformed fibrils (PFFs)

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2lyjbbplx9/v1>

Cole S Sitron¹, Victoria A Trinkaus¹, F Ulrich Hartl¹

¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry



Cole Sitron

Max Planck Institute of Biochemistry

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2lyjbbplx9/v1>

Protocol Citation: Cole S Sitron, Victoria A Trinkaus, F Ulrich Hartl 2024. Induction of aggregation in alpha-synuclein-expressing cells by treatment with preformed fibrils (PFFs). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.eq2lyjbbplx9/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 16, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 93747

Keywords: ASAPCRN, synuclein aggregates in cell, synuclein aggregate, synuclein through treatment, synuclein, fibrils of alpha, preformed fibril, misfolding of intracellular alpha, fibril, cells by treatment, intracellular alpha

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Abstract

This protocol details how to efficiently produce alpha-synuclein aggregates in cells by templating the misfolding of intracellular alpha-synuclein through treatment with preformed fibrils of alpha-synuclein (PFFs).

Attachments



[957-2488.docx](#)

16KB

Materials

Cell culture

Cells expressing alpha-synuclein

Reagents

- Appropriate cell culture medium
-  Lipofectamine™ 3000 Transfection Reagent **Thermo Fisher Scientific Catalog #L3000008**
-  Opti-MEM™ Reduced Serum Medium **Thermo Fisher Scientific Catalog #31985062**
-  PBS, pH 7.2 **Thermo Fisher Scientific Catalog #20012068**
- [M] 5 mg/mL alpha-synuclein preformed fibrils (PFFs) (see [dx.doi.org/10.17504/protocols.io.btynnpve](https://doi.org/10.17504/protocols.io.btynnpve) for purification protocol)

Equipment

- BioRuptor Plus sonicator (Diagenode cat. no. B01020001) (or equivalent)

Equipment	
Bioruptor® Plus sonication device	NAME
Sonication device	TYPE
Bioruptor®	BRAND
B01020001	SKU
https://www.diagenode.com/en/p/bioruptor-plus-sonication-device	LINK



Induction of alpha synuclein aggregation

15m 10s

- 1 Seed cells into plates such that they are ~25% confluent the following day.

Note

For our HEK cells, the correct seeding density is ~250 cells/mm², which is equivalent to 100,000 cells in a 12-well dish.

- 2 On the following day, warm PBS and OptiMEM to  37 °C , thaw an aliquot of PFFs, and (if using) allow lipofectamine to reach  Room temperature .

Note

Lipofectamine greatly increases the ability of PFFs to nucleate intracellular alpha synuclein aggregation, perhaps by altering the endocytic route that PFFs use to enter the cell. It is therefore not recommended to use lipofectamine if the route of PFF entry is a concern in your experiment. Lipofectamine is additionally excessively toxic and should therefore be avoided in some cell types.

- 3 Dilute the thawed PFFs 1:20 into an eppendorf containing PBS. 

Note

The final volume must be between  100 µL and  300 µL for proper sonication in the Bioruptor Plus.

- 4 Sonicate the PFFs in the Bioruptor Plus on high for 25 cycles of  00:00:05 on and  00:00:05 off at  4 °C .

10s

Note

Sonication breaks up the PFFs into smaller, more nucleation-competent fibrils.



5 Make a master mix of PFF and PBS solutions.

5.1  1 μL of sonicated & diluted PFFs should be added per 10 mm^2 of plate area (5 ug alpha-synuclein/10 mm^2).



Note

For a 12-well plate, this would be  40 μL of sonicated & diluted PFFs.

5.2 Master mixes should contain a ratio of  50 μL OptiMEM :  20 μL sonicated & diluted PFFs (or PBS as a control) :  3 μL of lipofectamine 3000.



Note

For a 12-well plate, 100 OptiMEM :  40 μL sonicated & diluted PFFs/PBS :  6 μL lipofectamine 300 should be used.

5.3 If using lipofectamine, first incubate the lipofectamine 3000 in OptiMEM for  00:05:00 before adding PFFs or PBS. Gently vortex upon addition of lipofectamine 3000.

5m



6 After addition of PFFs/PBS, gently vortex master mix solutions and incubate for  00:10:00 at  Room temperature .

10m



7 Add master mix solutions dropwise to cells, gently vortexing before adding to each well.





Protocol references

1. Luk KC, Song C, O'Brien P, Stieber A, Branch JR, Brunden KR, Trojanowski JQ, Lee VM. Exogenous alpha-synuclein fibrils seed the formation of Lewy body-like intracellular inclusions in cultured cells. *Proc Natl Acad Sci U S A*. 2009 Nov 24;106(47):20051-6. doi: 10.1073/pnas.0908005106. Epub 2009 Nov 5. PMID: 19892735; PMCID: PMC2785290.
2. Volpicelli-Daley LA, Luk KC, Lee VM. Addition of exogenous α -synuclein preformed fibrils to primary neuronal cultures to seed recruitment of endogenous α -synuclein to Lewy body and Lewy neurite-like aggregates. *Nat Protoc*. 2014 Sep;9(9):2135-46. doi: 10.1038/nprot.2014.143. Epub 2014 Aug 14. PMID: 25122523; PMCID: PMC4372899.