




Mar 12, 2025

iNeuron RNA-Extraction and RT-qPCR

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov1drryvr2/v1>



Benjamin O'Callaghan^{1,2}, Helene Plun-Favreau^{1,2}

¹Department of Neurodegenerative Diseases, UCL Queen Square Institute of Neurology, London, UK;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, USA, 20815



Benjamin O'Callaghan

UCL Queen Square Institute of Neurology

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1drryvr2/v1>

Protocol Citation: Benjamin O'Callaghan, Helene Plun-Favreau 2025. iNeuron RNA-Extraction and RT-qPCR. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.ewov1drryvr2/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: March 04, 2025

Last Modified: March 12, 2025

Protocol Integer ID: 123793

Keywords: ASAPCRN, RT-qPCR, PINK1, iNeuron, qpcr measurements of pink1 mrna, rna from ineuron, ineuron rna, pink1 mrna, ineuron, qpcr protocol for the extraction, rna, qpcr measurement, qpcr protocol

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP 000478

Abstract

Protocol for the extraction of RNA from iNeurons and RT-qPCR measurements of *PINK1* mRNA.

Materials

- Monarch Total RNA Miniprep Kit (New England Biolabs, T2010S)
- UltraPure DNase/RNase-Free Distilled Water (Invitrogen, 10977035)
- SuperScript IV Reverse Transcriptase (Invitrogen, 18090050)
- dNTP Mix (10 mM each) (Invitrogen, R0192)
- Random Hexamers (50µM) (Invitrogen, N8080127)
- RNaseOUT Recombinant Ribonuclease Inhibitor (Invitrogen, 10777019)
- Power SYBR Green PCR Master Mix (Applied Biosystems, 4367659)
- Gene Specific Primers (Sigma)
- NanoDrop One Spectrophotometer (Thermo Fisher Scientific)
- QuantStudio 7 Flex Real-Time PCR System (Applied Biosystems, RRID:SCR_020245)

RNA Extraction

- 1 iNeuron culture and initial collection for RNA extraction described in "iNeuron Differentiation and Culture in N2B27 vs BrainPhys for Immunofluorescence and Biochemistry Assessments of Mitophagy" (<https://www.protocols.io/private/E5851436F8F11EF86990A58A9FEAC02>). RNA is extracted using the Monarch Total RNA Miniprep Kit (New England Bioscience) as per the manufacturers instructions and with inclusion of gDNA binding column and the optional on-column DNase treatment. RNA is eluted from the column in  30 μL of Nuclease free H₂O
- 1.1 RNA is then quantified by 260nm absorbance using a NanoDrop One Spectrophotometer.

Reverse Transcription

- 2 100ng of RNA is reverse transcribed in a  10 μL reaction with  2.5 U/ μL SuperScript IV reverse transcriptase,  2.5 micromolar (μM) random hexamers,  0.5 millimolar (mM) dNTPs,  5 millimolar (mM) DTT and  2 U/ μL RNaseOUT.
A 100ng reverse transcription negative control reaction consisting of the same reaction mixture above but without reverse transcriptase is also conducted for each sample. The reverse transcription and negative control reactions are thermocycled at:

	Temperature	Time	Cycles
	25°C	10min	x1
	50°C	40min	
	80°C	10min	

Reverse Transcription Thermocycle

- 2.1 The cDNA product and reverse transcription negative control reactions are then diluted such that  100 ng of reverse transcribed RNA would be in a  120 μL final volume (i.e. add  10 μL reaction to  110 μL Nuclease Free H₂O =  0.83 ng/ μL).



RT-qPCR

- 3 3.33ng of cDNA ( 4 μ L of the [M] 0.83 ng/ μ l dilution) is used in each well of a  10 μ L RT-qPCR reaction using Power Sybr Green PCR master mix, [M] 500 nanomolar (nM) gene-specific primer pairs and the QuantStudio 7 Flex Real-Time PCR System.
- Each sample and gene target reaction is conducted with at least 2x technical replicates, with the same volume of reverse transcription negative control conducted alongside in a separate well for each gene target sample pair.

RT-qPCR Primer Pairs

PINK1_F: 5'- GTGGAACATCTCGGCAGGTT-3'

PINK1_R: 5'- CCTCTCTTGGATTTTCTGTAAGTGAC-3'

UBC_F: 5'-CTGGAAGATGGTCGTACCCTG-3'

UBC_R: 5'- GGTCTTGCCAGTGAGTGTCT-3'.

- 3.1 Ct values are determined automatically using the Quantstudio Real-time PCR Software and relative *PINK1* mRNA expression levels are calculated using the $2^{-\Delta\Delta C_t}$ method with *UBC* as the house-keeping gene.

Protocol references

Soutar MPM, Melandri D, O'Callaghan B, Annuario E, Monaghan AE, Welsh NJ, D'Sa K, Guelfi S, Zhang D, Pittman A, Tratzuni D, Verboven AHA, Pan KS, Kia DA, Bictash M, Gandhi S, Houlden H, Cookson MR, Kasri NN, Wood NW, Singleton AB, Hardy J, Whiting PJ, Blauwendraat C, Whitworth AJ, Manzonì C, Ryten M, Lewis PA, Plun-Favreau H. Regulation of mitophagy by the NSL complex underlies genetic risk for Parkinson's disease at 16q11.2 and MAPT H1 loci. Brain. 2022 Dec 19;145(12):4349-4367. doi: [10.1093/brain/awac325](https://doi.org/10.1093/brain/awac325). PMID: 36074904; PMCID: PMC9762952.