


Aug 24, 2022

Infrared PAR-CLIP

DOI

<https://dx.doi.org/10.17504/protocols.io.261ge4yxyv47/v1>



svetlana.lebedeva¹

¹BIMSB (Berlin Institute for Medical Systems Biology, Berlin)



svetlana lebedeva

bimsb

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ge4yxyv47/v1>

Protocol Citation: svetlana.lebedeva 2022. Infrared PAR-CLIP. [protocols.io](https://dx.doi.org/10.17504/protocols.io.261ge4yxyv47/v1)
<https://dx.doi.org/10.17504/protocols.io.261ge4yxyv47/v1>

Manuscript citation:

Control of immediate early gene expression by CPEB4 repressor complex-mediated mRNA degradation. Poetz F., Lebedeva S., Schott J., Lindner D., Ohler U., Stoecklin G. Genome Biology, 2022.



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 22, 2021

Last Modified: August 24, 2022

Protocol Integer ID: 54430

Keywords: CLIP, RNA-binding proteins, Immunoprecipitation, UV crosslinking, 4-thiouridine, mediated mrna degradation, mrna degradation, sites of rna, binding protein, rna, wide identification of rna, mrna, interactions at nucleotide resolution, control of immediate early gene expression, cpeb4 repressor complex, immediate early gene expression, gene expression, protein interaction, transcriptome, clip crosslinking, infrared par, genome biology, nucleotide resolution, microRNA, gene

Funders Acknowledgements:

DFG

Grant ID: SPP1935

Disclaimer

This protocol is based on multiple contributions from PAR-CLIP, iCLIP and easyCLIP. I acknowledge Harm Wessels, Antje Hirsekorn and Nico Kastelic for ideas, contributions and discussions on this protocol.

Abstract

Crosslinking and immunoprecipitation (CLIP) protocol based on 4-thiouridine incorporation and sequencing to determine binding sites of RNA-binding proteins (RBPs).

This protocol has been used for the manuscript "Control of immediate early gene expression by CPEB4 repressor complex-mediated mRNA degradation" by Poetz, Lebedeva et al. , Genome Biology 2022.

Based on the following protocols.

- easy CLIP

Porter DF, Miao W, Yang X, Goda GA, Ji AL, Donohue LKH, Aleman MM, Dominguez D, Khavari PA. easyCLIP analysis of RNA-protein interactions incorporating absolute quantification. Nat Commun. 2021 Mar 10;12(1):1569. doi: 10.1038/s41467-021-21623-4. PMID: 33692367; PMCID: PMC7946914.

- PAR-CLIP

Hafner M, Landthaler M, Burger L, Khorshid M, Hausser J, Berninger P, Rothballer A, Ascano M Jr, Jungkamp AC, Munschauer M, Ulrich A, Wardle GS, Dewell S, Zavolan M, Tuschl T. Transcriptome-wide identification of RNA-binding protein and microRNA target sites by PAR-CLIP. Cell. 2010 Apr 2;141(1):129-41. doi: 10.1016/j.cell.2010.03.009. PMID: 20371350; PMCID: PMC2861495.

- iCLIP

Huppertz I, Attig J, D'Ambrogio A, Easton LE, Sibley CR, Sugimoto Y, Tajnik M, König J, Ule J. iCLIP: protein-RNA interactions at nucleotide resolution. Methods. 2014 Feb;65(3):274-87. doi: 10.1016/j.ymeth.2013.10.011. Epub 2013 Oct 25. PMID: 24184352; PMCID: PMC3988997.

How to analyze the data:

- For the general pipeline to map CLIP-seq reads, see https://github.com/ohlerlab/clip_pipeline or https://github.com/slebedeva/CLIP_mapping. For specific pipeline to map the data from this protocol, look into - https://github.com/slebedeva/CPEB4_public.
- For the Ohler lab tool to call CLIP peaks, see <https://github.com/ohlerlab/omniCLIP> or <https://github.com/slebedeva/omniCLIP>.
- For the Ohler lab deep learning tool to predict RBP binding sites based on CLIP peaks see <https://github.com/ohlerlab/DeepRiPe>.

Materials

Primers

	A	B
	5' adapter	/5AzideN/GTTCAGAGTTCTACAGTCCGACGATC[CTGATC]rNrNrNrNrNrN
	3' adapter	/5rApp/NN NNT GGA ATT CTC GGG TGC CAA GGA AAA AAA AAA AA/iAzideN/AAA AAA AAA AAA /3Bio/
	small RNA RT Primer	GCCTTGGCACCCGAGAATTCCA
	small RNA PCR F Primer	AATGATACGGCGACCACCGAGATCTACACGTTTCAGAGTTCTACAGTCCGA
	small RNA PCR R Primer I1	CAAGCAGAAGACGGCATACGAGATCGTGATGTGACTGGAGTTCCTTGGCACCCGA GAATTCCA
	small RNA PCR R Primer I2	CAAGCAGAAGACGGCATACGAGATACATCGGTGACTGGAGTTCCTTGGCACCCGA GAATTCCA
	small RNA PCR R Primer I3	CAAGCAGAAGACGGCATACGAGATGCCTAAGTGACTGGAGTTCCTTGGCACCCGA GAATTCCA
	small RNA PCR R Primer I4	CAAGCAGAAGACGGCATACGAGATTGGTCAGTGACTGGAGTTCCTTGGCACCCGA GAATTCCA
	small RNA PCR R Primer I5	CAAGCAGAAGACGGCATACGAGATCACTGTGTGACTGGAGTTCCTTGGCACCCGA GAATTCCA
	small RNA PCR R Primer I6	CAAGCAGAAGACGGCATACGAGATATTGGCGTGACTGGAGTTCCTTGGCACCCGA GAATTCCA
	small RNA PCR R Primer I7	CAAGCAGAAGACGGCATACGAGATGATCTGGTGACTGGAGTTCCTTGGCACCCGA GAATTCCA
	small RNA PCR R Primer I8	CAAGCAGAAGACGGCATACGAGATTCAAGTGACTGGAGTTCCTTGGCACCCGA GAATTCCA

Optional: primers on your genomic DNA that make 140bp and 160bp product.

RNase free buffers:

Buy or autoclave with 0.1% DEPC

Optional:  5g DEPC (Diethylpyrocarbonate) **G-Biosciences Catalog #RC-041**

- RNase free water
- Tris-Cl pH 7.4 or 7.5
- Tris-Cl pH 6.5



- NaCl 5M
- EDTA 0.5M
- MgCl₂ 1M

Chemicals:

- IRDye 680RD DBCO Infrared Dye, 0.5mg 929-50005 LI-COR
- IRDye 800CW DBCO Infrared Dye, 0.5mg 929-50000 LI-COR
- RNA Clean and Concentration Kit Zymo
- Gel extraction kit Zymo
- Ethanol pure
- Culturing medium for respective cells
- PBS for cell culture
- Antibody against RBP of interest or tag
- Illumina small RNA barcoded PCR primers5x
- 5x TBE
- 6x DNA Gel loading dye (any)
- midori or any other agarose gel dye

Optional (for RNase concentration adjustment): Urea, Decade Markers, ³²P-γATP

☒ Acrylamide/bis-acrylamide, 40% solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A7802 Sigma**

☒ APS **Fisher Scientific Catalog #BP179-25** or any other

☒ TEMED **Bio-Rad Laboratories Catalog #1610801**

☒ 4-thiouridine (4sU) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T4509**

☒ IGEPAL-CA630 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #I3021 SIGMA-ALDRICH**

☒ Sodium deoxycholate **P212121**

☒ SDS, 10% Solution **Life Technologies Catalog #AM9822 (we use 20%)**

☒ Tween 20 **Bio-Rad Laboratories Catalog #170-6606-MSDS**

☒ DTT **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D0632**

☒ Protein G Magnetic Beads - 1 ml **New England Biolabs Catalog #S1430S**

☒ RNase I - 5,000 units **New England Biolabs Catalog #M0243S**

☒ RNasin(R) RNase Inhibitor, 2,500u **Promega Catalog #N2111**

☒ Suprase-In RNase Inhibitor **Thermofisher Catalog #AM2694**

☒ cOmplete™, Mini, EDTA-free (Protease Inhibitor) **Roche Catalog ##11836170001**

☒ T4 RNA Ligase 1 (ssRNA Ligase) - 1,000 units **New England Biolabs Catalog #M0204S**

☒ T4 RNA Ligase 2 truncated K227Q **New England Biolabs Catalog #M0351S**



✂ T4 Polynucleotide Kinase - 500 units **New England Biolabs Catalog #M0201S**

✂ NUPAGE LDS sample buffer (4x) **Thermo Fisher Scientific Catalog #NP0007**

✂ NuPAGE 4-12% Bis-Tris gel 1.0 mm 10 well **Thermo Fisher Scientific Catalog #NP0321BOX**

✂ Blotted nitrocellulose membrane **LI-COR Catalog #926-31090/926-31092** or any other 0.45 µm nitrocellulose membrane (use 0.22 µm for small proteins)

✂ PageRuler®; Prestained NIR Protein Ladder **Thermo Fisher Catalog #26635**

✂ 1x NUPAGE MOPS SDS running buffer (20x) **Thermo Fisher Scientific Catalog #NP0001**

✂ NuPage transfer buffer **Life Technologies Catalog #NP00061**

✂ Proteinase K **Thermo Fisher Scientific Catalog #EO0491**

✂ Oligo dT25 Magnetic Beads - 5 ml **New England Biolabs Catalog #S1419S**

✂ Superscript II Reverse Transcriptase **Thermo Fisher Scientific Catalog #18064071** (or III)

✂ SYBR Gold Nucleic Acid Gel Stain **Catalog # S-11494**

✂ 100bp DNA Ladder, 250ul (50 lanes) **Promega Catalog #G2101** or any other 100bp ladder

✂ High Resolution Agarose (For Nucleotides < 1kb) **Gold Biotechnology Catalog #A-202**

✂ NEBNext High-Fidelity 2X PCR Master Mix **NEB Catalog #M0541S**

Equipment and consumables:

- Pipettes
- Filter tips RNase free
- Eppendorf tubes RNase free
- Falcon tubes
- Cooled centrifuge for Falcon tubes
- Centrifuge for Eppendorf tubes
- Magnet
- Water bath
- UV crosslinker with 365nm lamp
- Ice and ice tray that fits into the crosslinker
- Cell culture hood and incubator
- Cell culture dishes, flasks, pipettes, cell scrapers
- NuPAGE gel running chamber
- NuPAGE cassettes and combs
- Filter paper (thick or thin depending on your preference)
- Transfer device (preferably for wet transfer)
- Infrared gel/membrane scanner (such as GE Typhoon)
- Printer

- scalpel
- optional: agarose gel cutting helpers
- Thermomixer
- Overhead tube rotator
- cold room
- PCR machine
- Qubit for nucleic acid quantification (DNA high sensitivity)
- Bioanalyzer or TapeStation with respective kits (DNA 1000)
- Agarose gel chamber, tray and comb
- power supply

Equipment

new equipment

NAME

Bioanalyzer 2100 instrument

BRAND

G2939BA

SKU

with RNA 6,000 Pico LabChip kit

SPECIFICATIONS

Equipment

Qubit

NAME

Flurometer

TYPE

Invitrogen

BRAND

Q33228

SKU

<https://www.thermofisher.com/order/catalog/product/Q33228>^{LINK}



Equipment

4200 TapeStation System

NAME

Electrophoresis tool for DNA and RNA sample quality control.

TYPE

TapeStation Instruments

BRAND

G2991AA

SKU

<https://www.agilent.com/en/product/automated-electrophoresis/tapestation-systems/tapestation-instruments/4200-tapestation-system-228263>

LIN
K

Buffer recipes:

Safety information

Wear mask if using powdered detergents!

Important notes:

1. Do not trust DTT in NEB buffers! Always add fresh DTT from 1M stock.
2. RNase inhibitors: in the lysis buffer, use an inhibitor that does not inhibit RNaseI (rRNAsin) and after IP, use inhibitors that inhibit all RNases (Suprase).

- iCLIP lysis buffer

	final concentration	for 500 ml
	50 mM Tris-HCl, pH 7.4	25 ml of 1M
	100 mM NaCl	10 ml of 5M
	1% Igepal CA-630	5 ml of 100%



	final concentration	for 500 ml
	0.1% SDS	2.5 ml of 20%
	0.5% sodium deoxycholate	25 ml of 10%

Use 500ml filter bottle to sterile filter. Store at 4°C. On the day of usage add:

- 1/100 volume of Protease Inhibitor Cocktail Set III (Roche EDTA-free) or 1 tablet of EDTA-free Pierce Inhibitor per 25ml
- RNase inhibitor: for tissues 1/1000 volume of ANTI-RNase, for cell lines 1/1000 rRNasin

- iCLIP high salt wash buffer

	final concentration	for 250 ml
	50 mM Tris-HCl, pH 7.4	12.5 ml of 1M
	1 M NaCl	50 ml of 5M
	1 mM EDTA	0.5 ml of 0.5M
	1% Igepal CA-630	2.5 ml of 100%
	0.1% SDS	1.25 ml of 20%
	0.5% sodium deoxycholate	12.5 ml of 10%

Use 250ml filter bottle to sterile filter. (We will need less of wash buffer than of the lysis buffer). Final NaCl concentration is 1M for anti-tag antibodies, it can be changed to 0.66M for antibodies against endogenous protein.

- iCLIP PNK-Tween buffer

	final concentration	for 10 ml	for 500 ml
	20 mM Tris-HCl, pH 7.4	0.2 ml of 1M	10 ml
	10 mM MgCl ₂	0.1 ml of 1M	5 ml
	0.2% Tween-20	20 µl of 100%	1 ml

Sterile filter large volume, use RNase free water and reagents if preparing small volume.

Note: I am lazy and heat up Tween a bit so I can pipet 100%. Feel free to prepare less concentrated solutions :-)

- iCLIP de-phosphorylation PNK pH 6.5 buffer



	final concentration	for 100 µl of 5x buffer
	350 mM Tris-HCl, pH 6.5	35ul of 1M
	50 mM MgCl ₂	5ul of 1M
	5 mM dithiothreitol	always add fresh 1 µl of 1M per 200ul 1x reaction

DTT has to be added fresh.

- easyCLIP Biotin IP buffer

	final concentration	for 10 ml
	0.1M Tris-HCl, pH 7.5	1 ml of 1M
	1M NaCl	2 ml of 5M
	0.1% Tween-20	10 µl of 100%
	1mM EDTA	20 µl of 0.5M

- easyCLIP NT2 buffer

	final concentration	for 10 ml
	50 mM Tris-HCl, pH 7.5	500 µl of 1M
	150 mM NaCl	300 µl of 5M
	1 mM MgCl ₂	10 µl of 1M
	0.5% NP-40	5 µl of 10%

- ProtK-SDS buffer

	final concentration	for 10 ml
	100 mM TrisCl pH 7.5	1 ml of 1M
	50 mM NaCl	100 µl of 5M



	final concentration	for 10 ml
	1 mM EDTA	20 μ l of 0.5M
	0.2% SDS	100 μ l of 20%

Important notes:

- 1. Do not trust DTT in NEB buffers! Always add fresh DTT from 1M stock.**
- 2. RNase inhibitors: in the lysis buffer, use an inhibitor that does not inhibit RNaseI (rRNAsin) and after IP, use inhibitors that inhibit all RNases (Suprase).**

Safety warnings

 Protect your eyes and open skin when working with UV-transilluminator.

Use mask when preparing detergent solutions from powder.

Infrared dye labeling of the adapter oligos

5h

- 1 Order the following oligos from IDT:

3' adapter:

/5rApp/NN NNT GGA ATT CTC GGG TGC CAA GGA AAA AAA AAA AA/iAzideN/ AAA
AAA AAA AAA /3Bio/

5' adapter:

/5AzideN/GTTCAGAGTTCTACAGTCCGACGATC[**CTGATC**]rNrNrNrNrNrNrN

Choose the smallest possible scale, and add RNase free HPLC purification. The amount of oligo and dye will be enough for several hundred CLIP reactions.

The bold sequence in the 5' adapter is an optional barcode.

There are 2 possibilities of multiplexing:

1. Pool several IPs before gel purification, and proceed with the library prep. In this case, order as many 5' adapters with different barcodes as the number of samples you plan to pool. There is no need to label them all with the dye, it is enough to have one adapter labeled and the rest unlabeled.

2. Standard Illumina multiplexing at the level of library PCR. In this case, order the 5' adapter **without the barcode**. If you have only one adapter with the barcode, make sure to **run at least 20% PhiX** on the flowcells with the CLIP sequencing, if you do not have other samples on the same run. Otherwise the clustering will fail due to the same sequence in the beginning of the read!

Resuspend the oligos in $X \mu\text{l}$ PBS (to final concentration of $100 \mu\text{M}$).

- 2 Resuspend the dye in the dark tube in $42.9 \mu\text{l}$ PBS.



Be careful!!!! The lyophilized pellet electrostatically sticks to the plastic tip and flies away super easily! - I did it under the cell culture hood and transferred $\sim 2 \times 20 \mu\text{l}$ to dark tubes.

- 3 Perform the labeling reaction:

2h



A	B
Dye in PBS	20 µl
100µM oligo	86 µl
PBS	to 200 µl (94 µl)

Incubate 37 °C 02:00:00

Either freeze the other half of the oligo and dye at -20°C and use within one week or use all the reagents on the single day to label all available oligos.

I label 3'adapter with 800nm dye and 5'adapter with 680nm dye. But you can also reverse it :-)

4 Purify the product:

Follow **Zymo RNA Clean and Concentration Kit protocol**. Use 2 high capacity columns per 100nmoles oligo. Wash with 85% EtOH.

I got around 100 µl ~ 50 µM oligo.

Aliquot and store working aliquot at -20 °C . I store main stocks at -80 °C

Growing, labeling and crosslinking cells

- 5 Starting material varies depending on the expression of your RBP of interest, the quality of your antibody and labeling time and efficiency. For the case of our test RBP, ELAVL1 (HuR) with the 3A2 antibody (sc-5261) in HeLa or HEK293 cells, 5× 15cm plates (~ 100M cells) should be enough to see a good signal after 2h labeling with 200 µM 4-thiouridine (4SU).

Important: longer labeling with higher 4SU concentrations could lead to cell toxicity.

- 5.1 Add 4SU diluted in DMSO to 1M or 0.5M concentration directly to the cell growth medium of the ~80% confluent cells.

Example:

1 plate = 20 ml medium = 8 µl of 0.5M 4SU



For multiple plates, I collect some medium from all plates, dilute the 4SU and distribute the medium again. We use conditioned medium.

5.2 During the labeling time, prepare crosslinking bench. You will need:

- the 365nm UV crosslinker
- an ice tray that can be put inside the crosslinker
- larger ice trays to put the plates on
- an ice box with a large bottle of cooled PBS (does not have to be sterile) and pre-labeled, cooled 50ml falcon tubes
- cell scrapers
- liquid waste bottle
- a centrifuge for Falcon tubes cooled down to 4°C
- a bucket of liquid nitrogen

Note: it is also possible to proceed with the lysis directly, if you only have a small sample and a long day. Usually I freeze the cell pellets to proceed later because crosslinking can take several hours depending on the number of plates.

5.3 After labeling, take the cells out of the incubator. Pour off the medium, replace with ice cold PBS and stack the plates onto the ice tray. Proceed in pairs (crosslinker fits 2 plates):

- take the lid off, pour off the PBS
- place onto the ice tray inside the crosslinker without lid
- crosslink (I use 150 mJ/cm²) (takes a minute or two)
- pour 5-10 ml of PBS onto each plate and scrape the cells off quickly into the falcons
- keep falcons on ice until all plates are processed

5.4  4000 rpm, 4°C, 00:10:00 Centrifuge the cells at 4000 rpm (top speed) for 10 min. Pour off the PBS and flash freeze the pellets in the liquid nitrogen. Store pellets at -80°C until use.

10m

Immunoprecipitation

30m

6 Couple 100 µl ProteinG beads with 10 µg of your antibody per large scale IP reaction. Wash the beads once with PBST (PBS-Tween).

Rotate the beads with antibody in 2x volume PBS-Tween for  00:30:00 at

 Room temperature

Cool down the Falcon centrifuge to  4 °C

Preheat water bath at  37 °C

30m



7 Meanwhile, lyse the cells (it will take time).

Prepare lysis buffer on ice: count with ~ 10-15 volumes per cell pellet. (I usually need ~50 ml buffer for 4 samples.)

Add protease and RNase inhibitor (see buffer recipes).

7.1 Get the pellets from the -80°C onto ice. First, pour ~10V of lysis buffer on top of all pellets. Then start pipeting gently one by one until the pellet has been dissolved. Avoid making too many bubbles. Then let them sit on ice for another 00:10:00

40m

If there are many pellets, I sometimes let them rotate overhead in the cold room to dissolve. It can take up to  00:30:00 .

Additionally, you can pass the lysate through a syringe several times to make sure that everything is lysed.

7.2 4000 rpm, 4°C, 00:10:00

10m

7.3 Transfer the cleared lysate into a new Falcon tube.

3m

If the volume is more than 10ml, I would split it into several 14ml tubes to be able to quickly heat up the whole volume.

Add RNase I (1:1000) on ice. Ideally, you would need to adjust the RNase concentration and digestion time for your BP by using radioactively labeled RNA and running it on a 15% TBE-Urea gel to estimate the size. Alternatively, you could either run the RNA gel after 3'adapter ligation and subtract the adapter length, or visualize the RNA size on a Tapestation or Bioanalyzer. However, I have not tried it and cannot tell how sensitive those methods are for the very small amounts of the IP RNA.

Bring the ice bucket and timer to the water bath. Put all the tubes at once into the water and digest exactly  00:03:00 . I quickly turn the tubes over about every 45 sec.

Immediately place the tubes back on ice (even better: ice water) and leave to cool for at least 3 min.

If desired, collect some of the input sample for protein extraction to control IP efficiency.

A note on the RNA input sample: I have found that RNA extraction from the digested lysate does not result in good libraries. I would rather recommend a separate, full RNA-seq library preparation from the same cells.



- 7.4 Wash once the beads coupled with the antibody and add the beads to the lysates. Rotate overhead at  4 °C for  02:00:00 .

2h

Washing and adapter ligation

20m

- 8 Collect the beads onto the magnet.
If desired, save some of the supernatant for protein extraction to control the IP efficiency.
- 8.1 Wash the beads 3 times with the high salt wash buffer, followed by 2 times of PNK-Tween washes.
- 8.2 Dephosphorylate the RNA.

20m

This step is from the iCLIP protocol and is needed to open circular 2',3' cyclic phosphate bond which is left behind by RNaseI to make the 3'-OH end available for ligation.

Component	amount, µl for 20 µl reaction
Water, RNase-free	15
5x pH6.5 buffer	4
SupersIn	0.5
T4 PNK	0.5
fresh DTT	to 5 mM

Incubate in Thermomixer :  37 °C  00:20:00  1100 rpm

- 8.3 Wash the beads 2x with PNK-Tween buffer.
- 8.4 **Ligate the 3' adapter.**

We connect the 3'OH end of the IP-RNA to the 5'-App end of the adapter with the help of truncated, mutated RNA ligase 2.



Component	amount, μl for 20 μl reaction
Rnase-free water	11
3' adapter, labeled, from step 4 (diluted to 1 pmol/ μl)	1
PEG8000 50%	4
10x T4 RNA ligase buffer	2
Rnl2TrK227Q	1
SuperaIn	1
fresh DTT	to 5 mM

16 °C

Overnight

1100 rpm 15 sec on / 15 sec off

8.5 Wash the beads 2x with PNK-Tween buffer.

8.6 5'end phosphorylation

20m

Converts 5'-OH end left by RNaseI to 5'-phosphate.

Component	amount, μl for a 20 μl reaction
T4 PNK	2
SuperaIn	1
Rnase-free water	13
10x PNK buffer	2
10mM ATP	2
fresh DTT	to 5 mM

37 °C

00:20:00

1100 rpm

8.7 Wash the beads 2x with PNK-Tween buffer.



8.8 Ligate the 5' adapter

2h

Component	amount, μ l for 20 μ l reaction
T4 Rnl1 buffer 10x	2
PEG8000, 50%	10
10mM ATP	2
5' adapter from step 4 (diluted to 1 pmol/ μ l)	1
Supersasin	0.5
T4 RNA ligase 1 (normal or high concentration)	1
Rnase-free water	4
fresh DTT	to 5 mM

25 °C

02:00:00

1100 rpm 15 sec on / 15 sec off

8.9 Wash the beads 2x with PNK-Tween buffer.

8.10 Resuspend the beads in 40 μ l 1x LDS NuPAGE sample buffer.

5m

70 °C

00:05:00

1100 rpm

Running the gel, transfer and imaging

9 Run the gel

1h 30m

Prepare NuPAGE chamber with 1x MOPS running buffer and 4-12% Bis-Tris gel.

Using 1mm 10well combs, I run each sample on 2 adjacent lanes alongside with 5 μ l / lane of the NIR (near infrared) prestained protein ladder.

If you have time, reduce the voltage to avoid distortion of bands (I use 120V for

01:30:00).

While the gel is running, prepare transfer.

10 Transfer

Notes:

I prefer wet overnight transfer, because it is more complete. Proteins crosslinked to RNA transfer poorly with fast transfer such as iBlot. However, if your RBP is very abundant, like ELAVL1, also incomplete transfer can give enough material to prepare a good library.

For the same reason I use transfer buffer without methanol (to prevent gel shrinkage, which leads to even poorer recovery of crosslinked complexes). I add small amount of SDS to facilitate transfer.

Always use nitrocellulose membrane.

Wet transfer:

- Prepare ~ 2L of 1x NuPAGE transfer buffer with 0% methanol and 0.0375% SDS.
- Cut 2 thick or 6 thin filter papers per each gel to be transferred.
- Cut a membrane of the size of the gel.
- Soak all papers, transfer pads and membrane in the transfer buffer.
- Take out the gel, open up the cassette and remove only one side of the cassette (I use this trick to stabilize the gel and avoid tearing it).
- Soak the gel lying on the half-cassette in the transfer buffer.
- Fill up the chamber with transfer buffer.
- Assemble the sandwich starting from the bottom: mnemonic rule is **"gel to black"**
 1. **Red** half-sandwich
 2. Transfer pad (sponge)
 3. Thick filter paper (or 3 thin papers)
 4. Membrane
 5. Take the gel still lying on plastic half-cassette, turn it over, I usually try to have the half with the slit, then I can use the cassette opening tool to push the gel through the slit off the cassette, while holding it over the sandwich.
 6. You can use a roll like the one that comes with iBlot to remove air bubbles between gel and membrane
 7. Thick filter paper (or 3 thin papers)
 8. Transfer pad (sponge)
 9. **Black** half-sandwich

Place sandwich(es) into the chamber minding again the black-to-black rule.

Transfer at constant voltage of 30V in the cold room  4 °C  Overnight

11 Imaging

Next morning disassemble the cassette. I keep the membrane on the filter papers to avoid drying and image it wet. If the membrane is partially dry, the scan will have bright dry-wet boundaries and difficult to assess.

Scan both IR channels (short, 680nm and long, 800nm). NIR ladder will be short.

Use Fiji software to merge the channels with false colors. (Image → Color → Merge channels...)

Two colors have to give a distinct color when combined. Locate the merged band (for green and red, it will be yellow) which will indicate the complex ligated to both adapters.

Software

FIJI (Image J)

NAME

NIH

DEVELOPER

<https://fiji.sc/>


SOURCE LINK

As an example, I used blue color for short and red for long channel, so we need to look for a purple band.

We look around the size of our protein plus the size of the ≥ 20 nt RNA (~6kDa) plus the 2 adapters (~12kDa+~15kDa).


Example: 100kDa RBP $\Rightarrow$ look for ~135-140 kDa band.

I imaged the gel before transfer, where one can see different by-products: unligated adapters, adapter dimers, RBP-RNA complexes ligated to only one of the two adapters, and sometimes higher molecular weight complexes. Also NuPAGE loading dye lights up in the short channel.



gel before transfer

After transfer the adapters pass through the membrane, so we mostly see single and double ligated complexes with sometimes higher molecular weight ones:



membrane after transfer

Print out the channel where the double ligated band is seen best in original size and high contrast black and white.



RNA isolation and library preparation

2h 13m

- 12 Use printed out black and white papers to cut out the double-ligated band. I use a white light table (or a tablet under a big square plate lid), put the printout on top, then a transparent film treated with RNase away, and the membrane.

Cut the band out with the scalpel and place it into an RNase-free Eppendorf tube.

12.1 Proteinase K digestion

45m

This step will remove the RBP crosslinked to the RNA, leaving behind only one or few amino acids crosslinked to a 4SU nucleotide.

- Add 375 µl ProtK-SDS buffer to the membrane in the tube.
- Add 25 µl Proteinase K

55 °C

00:45:00

1100 rpm

12.2 Purification of the RNA on oligo-dT beads

20m

Follows easyCLIP protocol.

- Use 20 µl oligo-dT beads per sample in a DNA LoBind tube
- Wash the beads with BIB buffer and leave in 600 µl BIB
- Transfer the Proteinase K extract to the beads
- Overhead rotation in the cold room 4 °C 00:20:00
- Wash 1x BIB, 2x NT2 buffer and 3x with ice cold PBS with added RNase inhibitor

12.3 Reverse transcription

1h 8m

Use standard Illumina small RNA RT primer: GCCTTGGCACCCGAGAATTCCA

I use a magnet for PCR tubes for this step, but one could also use normal 1.5ml tubes and heat in the thermocycler.

- Prepare double number of PCR tubes as your number of samples
- Heat up the PCR block to 95°C
- Prepare the Superscript kit on ice



- After the last wash, resuspend oligo-dT beads in a mix of 12.4 µl RNase-free water + 3µl RT primer
- Heat 🔥 95 °C ⌚ 00:03:00 , immediately transfer to the magnet and quickly collect the supernatant into a fresh PCR tube
- Place the tubes on ice for ⌚ 00:05:00
- Add 6.7 µl of the master mix:

component	amount, µl per 1 reaction
5x Superscript buffer	4
10 mM dNTPs	1
0.1M DTT	1
Superscript III	0.5
Superscript III	0.5
Rnase-free water	1

🔥 50 °C ⌚ 01:00:00

12.4 Library PCR

Use Illumina small RNA barcoded primers (1 sample = 1 barcode).

You can also order them yourself (see materials).

Use the same forward primer and different barcode index reverse primer for each sample to be multiplexed.

Use NEBNext High Fidelity 2x PCR Master Mix.

12.5 Pilot PCR

Perform small scale PCR to determine the number of cycles.

Depending on the scarcity of the sample, I perform 3-5 pilot reactions.

The mix for 5 reactions:



	Component	amount for 5 reactions, μ l
	cDNA	5
	small RNA forward primer	1
	small RNA reverse barcode (index) primer	1
	NEBNext PCR MM	25
	water	18

Pipet the mix onto 5 PCR tubes 10 μ l each.

12.6 Perform the following PCR program:

5m 5s

(place samples in cycler when T > 80°C)

 00:02:00  95 °C

N cycles of:

 00:00:20  95 °C

 00:00:30  60 °C

 00:00:15  72 °C

Final extension:

 00:02:00  72 °C

for different amount of cycles (depending on your starting material and RBP expect different results).

For example, I use 17, 20, 23, 26, 29 if I expect less material. Or use 14, 16, 18, 20, 22 if expecting a lot.

Take out the tubes one by one after respective number of cycles and keep on ice.



12.7 Add 2 μ l 6x DNA loading dye to each tube.

3m

Run 6% TBE PAGE gel with a 100bp or another low molecular weight DNA ladder.

Stain the gel for  00:03:00 in 1x SYBR Gold (stock is 10,000x) in 1x TBE.

Determine the optimal number of cycles.

The desired productive library size is between 157 bp (expecting 20nt long RNA insert) and 187bp (50 nt insert). Because we rely on RT readthrough, unlike iCLIP, we try to sequence ideally reads of length around 30nt, which are long enough for mapping but short enough to give us single nucleotide resolution.

The adapter dimer is 137bp.

The full PCR product sequence (sense) will be:


```
AATGATACGGCGACCACCGAGATCTACAC GTTCAGAGTTCTACAGTCCGACGATC CTGATC
NNNNNNNN [CLIP_RNA 20-50nt] NNNNTGGAATTCTCGGGTGCCAAGGAACTCCAGTCAC
[XXXXXX] ATCTCGTATGCCGTCTTCTGCTTG
```

N's denote unique molecular identifiers (UMIs): random bases which purpose is:

1. Reduce RNA ligase biases
2. Identify unique ligation events and separate them from PCR duplicates.

X's denote Illumina index used for multiplexing samples on the same sequencing lane.

Be aware that depending on the RNaseI digestion efficiency, you may not get a very sharp band.



For this library, I would use 21 or 22 cycles.

Usually it is not recommended to go higher than 22-23 cycles, because low amount of starting material and high duplication rates will waste a lot of sequencing depth, and it will be difficult to analyze the scarce data. If your library comes up very late, try to increase the amount of starting material, overexpress the RBP, if possible, or use a more efficient antibody (in case your RBP is not tagged, make a FLAG-tagged cell line).

12.8 Assemble final library amplification:

component	amount, μ l
rest cDNA	15-17



component	amount, μ l
NEBNext 2x PCR MM	25
small RNA forward primer	1
small RNA reverse barcode (index) primer	1
water	to 50

Run the same PCR program for the chosen number of cycles.

12.9 Gel clean up.

We use 2.5% high resolution agarose/TBE gel for final library cleanup.

It can be that your library will be almost invisible on this gel - do not panic and proceed, you do not need very much for a good sequencing run!


To guide purification, you can design genomic DNA primers which give a PCR product of ~160bp and ~140bp. Run them alongside your samples and cut out around 160bp, avoiding the 140bp.

- pour thick-welled 2.5% agarose gel
- leave enough wells for several loads of marker and ideally empty lanes in between samples
- add 6 μ l of 6x DNA loading dye to PCR
- run on the gel alongside the size markers (check the separation of 140 and 160bp band to decide when to stop)
- excise the 160bp band (also when invisible) trying to avoid the 140bp band
- extract library using Zymo gel extraction kit (elute in 11 μ l)

12.10 Library validation and sequencing

- Determine library concentration with Qubit (it will be low, 1-10ng/ μ l)
- Run the samples on Bioanalyzer or Tapestation to confirm that you cut out correct size.

Here is an example trace of a successful library run on a D1000 Tapestation:
There is a peak of 160bp or more, and no adapter dimer or 140bp.



- Sequence the samples on Illumina NEXTSeq500 or comparable machine. Include PhiX if running CLIP samples alone.
- See the description on suggested pipelines to analyze the data.

Troubleshooting

13

Problem	Solution
No double ligated band is seen on the membrane after transfer	Troubleshoot possible steps at which material loss could have occurred (see below)
1. Not enough material	Try to increase the amount of input lysate at least 10-fold. Try to overexpress the RBP.
2. Immunoprecipitation is not working	Collect input, supernatant in IP samples in protein buffer for Western blot assessment of immunoprecipitation efficiency. Note that crosslinking rate of RBPs on average is only ~1%, so do not expect to see a Western blot signal where infrared ligated complexes are. If IP is not efficient, try a different antibody or tag the RBP of interest endogenously. FLAG tag usually works well.
2a. RBP is not in the soluble fraction	If you cannot detect your RBP in the input, try checking by Western blot the insoluble pellet after lysate centrifugation at step 7.2. Possibly you may need to change lysis buffer or add sonication if your RBP is very chromatin-bound.
3. Transfer is not working	Image the gel before and after transfer to the membrane. If the double-ligated band

Problem	Solution
	is still seen in the gel after transfer, try to increase transfer efficiency by reducing methanol, increasing SDS in transfer buffer and transferring for longer time. Try to increase RNase concentration to get shorter fragments. Alternatively, you can try the original PAR-CLIP protocol by cutting the band directly out of the gel and electroeluting the crosslinked complexes directly out of the gel piece.
4. Only adapter dimer and 3' end ligation product is seen on the gel, but not 5' or double ligation product	DTT in ligation buffer deteriorated. This can happen even you aliquot fresh NEB buffer and use aliquots only once. Always add fresh DTT diluted from 1M stock to 5 mM.
No library amplification in pilot PCR	
1. Not enough material	Extend pilot PCR up to 35-40 cycles. If the band of correct size is detected, try to increase the amount of input lysate at least 10-fold.
2. Not efficient Proteinase K digestion	Do use SDS-containing Proteinase K buffer. I had problems with digestion step when using urea-containing buffer.
3. Oligo-dT purification is not working	Check supernatants from the oligo-dT bead purification on an RNA urea-PAGE gel.
4. Reverse transcription is not working	Check your RT enzyme by performing an independent RT reaction with a positive control.
5. PCR is not working	Include a positive control to check the PCR master mix. Use the same primers for pilot PCR as for the main PCR. Check that your primers are from the correct kit (Illumina small RNA kit). (No joke: in case of NEXTflex kits, both small RNA and mRNA forward primers are called simply "Universal forward primer", and if it once gets into the wrong box, you will get no amplification).

To troubleshoot phosphorylation, ligation, RT and PCR steps, you can use a positive control from the small RNA kit or order yourself a small RNA oligo of length ~20nt. PCR product from these reactions can subsequently be used as a size marker to purify the real library.

Also, Ule lab put together a very [**detailed forum**](#) about troubleshooting all possible steps of CLIP. They are very friendly and responsive, do not hesitate to ask!

