

Nov 03, 2025

Version 1

Integrated Diagnostic Protocol for Syphilis Using Immunochromatographic HIV/Syphilis Combo, RPR, and TPHA Tests Detection of Antibodies anti *Treponema pallidum* V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.261gek1wjg47/v1>

Faramalala Ramamonjiarisoa, MD^{1,2,3}, HDR RAZAFIARIMANGA Zara Nomentsoa^{4,5}, Dr. Ramaroson Roseline PhD⁶

¹Doctoral Candidate, Doctoral School of Life and Environmental Sciences (SVE);

²Research Team: EAD 7 – Immunology, Immunopathology, Immunodiagnostics;

³Medical Biologist, Research and Training Laboratory in Medical Biology (LBM), Faculty of Medicine, University of Antananarivo, Madagascar;

⁴Associate Professor – Director of Documentation, International Publications, and Editing, University of Antananarivo;

⁵President of the Faculty College of Teachers, Faculty of Sciences, University of Antananarivo;

⁶Specialization in Fundamental and Applied Biochemistry, Faculty of Sciences, University of Antananarivo, Madagascar



framamonjiarisoa

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261gek1wjg47/v1>

Protocol Citation: Faramalala Ramamonjiarisoa, MD, HDR RAZAFIARIMANGA Zara Nomentsoa, Dr. Ramaroson Roseline PhD 2025. Integrated Diagnostic Protocol for Syphilis Using Immunochromatographic HIV/Syphilis Combo, RPR, and TPHA Tests Detection of Antibodies anti Treponema pallidum. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.261gek1wjg47/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 31, 2025

Last Modified: November 03, 2025

Protocol Integer ID: 231230

Keywords: Treponema pallidum, CTK Biotech, Syphilis, HIV, RPR, TPHA, Rapid labs, ISO 15189, CLSI, WHO, Serological Testing, Diagnostic Algorithm, integrated diagnostic protocol for syphili, syphilis combo test for initial screening, syphilis combo test, tpha tests detection of antibody, standardized diagnostic workflow for syphili, syphilis combo, antibodies for both hiv, accurate syphili, treponema pallidum hemagglutination assay, syphili, confirming treponemal infection, using immunochromatographic hiv, hiv diagnosis, simultaneous hiv screening, tpha tests detection, immunochromatographic hiv, treponemal infection, specific antibody, syphilitic individual, treponema pallidum this protocol, antibody, hiv, antigen, serum bind to these antigen, treponema pallidum, specific antibodies in the plasma, rpr detects reagin, integrated diagnostic protocol, cholesterol particles in the test, diagnostic accuracy, overall diagnostic process

Abstract

This protocol presents a standardized diagnostic workflow for syphilis using three complementary laboratory methods: (1) Immunochromatographic HIV/Syphilis Combo test for initial screening, (2) Rapid Plasma Reagin (RPR) test for detecting non-treponemal antibodies, and (3) Treponema pallidum Hemagglutination Assay (TPHA) for confirming treponemal infection. The combined approach improves diagnostic accuracy and facilitates simultaneous HIV screening. Immunochromatographic HIV/Syphilis Combo test is designed to simultaneously detect antibodies for both HIV and syphilis, providing a rapid and efficient initial screening method. RPR detects reagin, an antibody-like substance found in syphilitic individuals and sometimes in those with other conditions. The reagin binds to the cardiolipin-coated cholesterol particles in the test, causing flocculation. TPHA uses avian erythrocytes coated with T.pallidum antigens. Specific antibodies in the plasma or serum bind to these antigens, causing erythrocytes to agglutinate, forming a characteristic pattern in the test well. Non-specific reactions are prevented by absorbents. The integration of these methods enhances the overall diagnostic process, ensuring both accurate syphilis and HIV diagnosis.

Guidelines

This protocol complies with IFU, ISO 15189 (Medical Laboratories – Requirements for Quality and Competence), CLSI EP05/EP09/EP12 guidelines, and WHO recommendations for sexually transmitted infection testing and quality assurance in medical laboratories.

Materials

- OnSite® HIV/Syphilis Ab Combo Rapid Test (CTK Biotech)
- BD Macro-Vue™ RPR Card Test kit
- Rapid Labs TPHA kit
- Micropipettes and sterile tips
- Clean glass slides or test cards
- Timer, centrifuge, and mixing rotator
- Personal protective equipment (gloves, lab coat, goggles)
- Sample types: serum, plasma, or whole blood

Safety warnings

- ! Perform procedures under Biosafety Level 2 (BSL-2) conditions. Wear gloves, laboratory coat, and protective eyewear. Dispose of waste according to local biosafety and environmental regulations. Treat all biological specimens as potentially infectious. Dispose of used materials in accordance with institutional biosafety and national biomedical waste guidelines. Sample Handling: Acceptable specimens: Serum or plasma (per kit IFU). Store samples at 2–8 °C up to 7 days or ≤ -20 °C for long-term storage. Avoid repeated freeze-thaw cycles.

Ethics statement

The study protocol was reviewed and approved by the Ethics Biomedical Committee of the University of Antananarivo. All procedures were conducted in accordance with national research ethics guidelines. Written informed consent was obtained from all participants. Calibrate pipettes and microplate readers. Store reagents at 2–8 °C and allow to reach room temperature before use. Verify internal quality controls (positive/negative sera) before performing

Before start

The study protocol was reviewed and approved by the Ethics Biomedical Committee of the University of Antananarivo. All procedures were conducted in accordance with national research ethics guidelines. Written informed consent was obtained from all participants. Calibrate pipettes and microplate readers. Store reagents at 2–8 °C and allow to reach room temperature before use. Verify internal quality controls (positive/negative sera) before performing



Before Start

- 1 Calibrate pipettes and microplate readers. Store reagents at 2–8 °C and allow to reach room temperature before use. Verify internal quality controls (positive/negative sera) before performing.

Safety Warnings

- 2 Perform procedures under Biosafety Level 2 (BSL-2) conditions. Wear gloves, laboratory coat, and protective eyewear. Dispose of waste according to local biosafety and environmental regulations.
Treat all biological specimens as potentially infectious. Dispose of used materials in accordance with institutional biosafety and national biomedical waste guidelines. Sample Handling: Acceptable specimens: Serum or plasma (per kit IFU). Store samples at 2–8 °C up to 7 days or ≤ -20 °C for long-term storage. Avoid repeated freeze-thaw cycles.

Laboratory procedure

- 3 Allow reagents and samples to reach room temperature (20–25°C).
- 4 Add 50 μ L of specimen into the test well.
- 5 Add two drops (≈ 80 μ L) of sample diluent.
- 6 Read results at 15 minutes.
 - 6.1 - One control line: Negative.
 - 6.2 - Control + test line: Positive.
 - 6.3 - No control line: Invalid.
- 7 Equilibrate reagents and samples to room temperature.

- 8 Place 50 μ L of serum on a test circle.
- 9 Add 18 μ L of antigen suspension.
- 10 Rotate the card for 8 minutes at 100 ± 2 rpm.
- 11 Examine for macroscopic flocculation.
- 12 For quantitative determination, perform serial twofold dilutions.
- 13 Dilute serum 1:20 with sample diluent.
- 14 Add 25 μ L of diluted sample to each well.
- 15 Add 75 μ L of Test Cells to test wells and Control Cells to control wells.
- 16 Mix gently and incubate for 45–60 minutes at 15–30°C.
- 17 Interpret results:
 - 17.1 - Reactive: diffuse mat formation.
 - 17.2 - Non-reactive: compact button formation.

Interpretation of Results



- 18 Step 1: Perform HIV/Syphilis Combo rapid test.
- If reactive positive or no → confirm with RPR.
Step 2: If RPR positive or no → perform TPHA.
Final interpretation:
- RPR (+) / TPHA (+): Active syphilis.
- RPR (–) / TPHA (+): Past/treated infection.
- RPR (+) / TPHA (–): Biological false positive.

Quality Control

- 19 Include positive and negative controls in each run. Do not use expired reagents. Repeat the test if control results are invalid. Participate in external quality assurance programs periodically.

Data Storage and Confidentiality

- 20 LIMS Database: All raw data and test results are stored securely in the Laboratory Information Management System (LIMS), ensuring traceability and compliance with ISO 15189 standards. Confidentiality: Patient data are handled confidentially, and all personal information is stored securely in compliance with ethical guidelines and privacy laws.

Limitations

- 21 Early primary infection may yield false negatives.
Biological false positives may occur in pregnancy, autoimmune diseases, or viral infections.
Cross-reactivity with yaws and other treponematoses possible.
Qualitative tests require clinical correlation.

Workflow Algorithm

- 22 Sample Reception → Verification of identification → Serum/Plasma preparation → Initial Rapid Screening (CTK Biotech Syphilis HIV Ab PLUS Combo)
RPR BD Macro-Vue® Screening → If Reactive or no
TPHA Rapid Labs Confirmation
Consolidate all results → Quality Control verification → Biological validation
Final report issue



Protocol references

1. CTK Biotech. OnSite® HIV/Syphilis Ab Combo Rapid Test IFU
2. BD Diagnostics. Macro-Vue™ RPR Card Test IFU
3. Rapid Labs. TPHA 100–500T Rev 9, CE Single, 2023
4. WHO Laboratory Manual for STI Diagnosis, 2019
5. Larsen SA et al., J Clin Microbiol, 1981
6. Manufacturer Instructions for Use (IFU) for each specific kit
7. CLSI documents (e.g., EP05, EP12) for precision and qualitative test evaluation.
8. CDC. Sexually Transmitted Infections and Viral Hepatitis testing recommendations.