



May 04, 2024

# Intracellular and extracellular MHC-II immunostaining and microscopy on peritoneal macrophages

DOI

[dx.doi.org/10.17504/protocols.io.bp2l62191gqe/v1](https://dx.doi.org/10.17504/protocols.io.bp2l62191gqe/v1)

Rebecca Wallings<sup>1</sup>

<sup>1</sup>University of Florida



Rebecca Wallings

University of Florida

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l62191gqe/v1>

**Protocol Citation:** Rebecca Wallings 2024. Intracellular and extracellular MHC-II immunostaining and microscopy on peritoneal macrophages. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bp2l62191gqe/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working



**Created:** May 03, 2024

**Last Modified:** May 04, 2024

**Protocol Integer ID:** 99196

**Keywords:** extracellular mhci in peritoneal macrophage, microscopy on peritoneal macrophage, extracellular mhc, extracellular mhci, peritoneal macrophage, ii immunostaining, using evos fluorescent microscope, evos fluorescent microscope, microscopy

## Abstract

Staining of intracellular and extracellular MHCII in peritoneal macrophages in 96 well plate using EVOS fluorescent microscope and analysis on CellProfiler

## Troubleshooting

- 1 After harvesting and plating pMacs from mouse, plate in 96-well plate at 50,000 cells/well in 200uL as described in pMac harvesting protocol. Once cells have adhered, aspirate media and replace with fresh growth media containing vehicle or 100U of IFN $\gamma$  for 18
- 2 Retrieve cells from 37 degrees. Aspirate media and wash 3x DPBS+/+
- 3 Incubate cells with 25 ug/mL of APC-MHC-II (Biolegend, RRID:AB\_313329) in DPBS+/+ containing 1:100 FcR blocking reagent, final volume 50uL, for 30 minutes at room temperature, protected from light
- 4 Wash cells 3 x DPBS+/+
- 5 fix cells by incubation in 4% PFA for 10 minutes at room temperature, final volume of PFA 50uL, and then washed 3x with DPBS+/+
- 6 Permeabilize cells with 100uL permeabilization buffer (eBiosciences, #88-8824-00) on ice for 15 minutes
- 7 Spike in 25ug/mL of PE-610-MHC-II (Biolegend, RRID:AB\_2574618) into permeabilization buffer and incubate cells for 30 minutes at room temperature protected from light.
- 8 Wash cells 3 x DPBS+/+
- 9 Incubate cells in 1  $\mu$ g/ml DAPI (Invitrogen, RRID:AB\_2629482) for 10 minutes at room temperature in DPBS+/+, final volume 50uL.
- 10 Image cells on an EVOS™ M7000 (Invitrogen) at 20 x magnification.
- 11 Perform image analysis using Cellprofiler 4.2.5 (RRID:SCR\_007358). Using the 'IdentifyPrimaryObject' module in Cellprofiler, identify icMHC and exMHC in their respective channels and quantify MFI and calculate Ex:IcMHCII from these quantified MFI values.