

Oct 26, 2023

Version 2

## 🌐 Introducing GolgiTAG to Cells and Immunoprecipitation of Golgi V.2

📖 [Proceedings of the National Academy of Sciences of the United States of America](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvrjdjrogmk/v2>

Rotimi Y. Fasimoye<sup>1</sup>, Dario R. Alessi<sup>1</sup>


<sup>1</sup>Medical Research Council Protein Phosphorylation and Ubiquitylation Unit, School of Life Sciences, University of Dundee, Dow Street, Dundee DD1 5EH, UK

Team Alessi



**Rotimi Fasimoye**

University of Dundee



### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.6qpvrjdjrogmk/v2>

**External link:** <https://doi.org/10.1073/pnas.2219953120>



**Protocol Citation:** Rotimi Y. Fasimoye, Dario R. Alessi 2023. Introducing GolgiTAG to Cells and Immunoprecipitation of Golgi. protocols.io <https://dx.doi.org/10.17504/protocols.io.6qpvrjdjrogmk/v2> Version created by **Rotimi Fasimoye**

**Manuscript citation:**

Fasimoye R, Dong W, Nirujogi RS, Rawat ES, Iguchi M, Nyame K, Phung TK, Bagnoli E, Prescott AR, Alessi DR, Abu-Remaileh M, Golgi-IP, a tool for multimodal analysis of Golgi molecular content. Proceedings of the National Academy of Sciences of the United States of America 120(20). doi: [10.1073/pnas.2219953120](https://doi.org/10.1073/pnas.2219953120)

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 25, 2023

**Last Modified:** October 26, 2023

**Protocol Integer ID:** 89870

**Keywords:** golgi, golgi apparatus, anti-HA magnetic beads, immunoprecipitation, IP, TMEM115, ASAPCRN, purification of golgi, methods for golgi purification, golgi purification, golgi from the mammalian hek293 cell line, integral membrane protein of the golgi complex, immunoprecipitation of golgi golgi apparatus, golgitag to cell, purifying golgi, golgi golgi apparatus, endoplasmic reticulum, essential to the secretory pathway, cell for secretion, golgi, introducing golgitag, golgi complex, secretory pathway, tagged integral membrane protein, other organelle, vesicle, other organelles within the cytoplasm, organelle, purification, mammalian hek293 cell line, main transportation hub of the cell, various mammalian cell line, protein, cell, various cell line

## Abstract

Golgi apparatus is essential to the secretory pathway of the cell and by extension helps to maintain cellular homeostasis. It is the main transportation hub of the cell, where molecules destined for other organelles within the cytoplasm or outside the cell for secretion are packaged into vesicles. The purification of Golgi using available techniques is time-consuming and laborious, as it requires various centrifugation steps, sometimes in gradient buffers, which reduce yield, compromise its integrity, and increase the chance of contamination, especially from endoplasmic reticulum. Furthermore, depending on the cell type or tissue, the methods for Golgi purification differs. Here, we present a protocol for rapidly purifying Golgi from the mammalian HEK293 cell line, using high affinity anti-HA magnetic beads. This method relies on the immunoprecipitation, in phosphate buffer, of HA-tagged integral membrane protein of the Golgi complex, TMEM115. The HA-tagged TMEM115 expressing vector is packaged into a lentivirus, therefore various mammalian cell lines can be transduced, giving stable expression levels. Our protocol is fast, approximately 10 min, and can be used on various cell lines and tissues without any modification. The Golgi purified using this method are highly enriched, intact, contaminant-free and, depending on solubilisation buffer, could be used for various downstream applications, such as proteomics and metabolomics.

In this new version, the **last step** contains a supplemental video with extra context and tips, as part of the ASAP Protocol Particulars, featuring conversations with protocol authors

## Attachments




Golgi stack showing ...

1.1MB

## Guidelines

Fig. 1



### Figure 1: GolgiTag immunoprecipitation.

(a) The sketch diagram of a Golgi cisternae showing the integral protein TMEM115 (green) with 3xHA tag (red) at the C-terminus. Both the N- and the C-terminal are cytosolic. (Note: Diagram is not to scale).

(b) The flowchart of the Golgi Immunoprecipitation. Cells are lysed in KPBS using Dounce homogeniser, before Golgi is captured using magnetic anti-HA beads. (Note: Diagrams are not to scale)

# Materials

## Materials

### 1. Cell lines

1. HEK293FT for virus packaging and propagation (Invitrogen™. Catalog# R70007)
2. HEK293 (ATCC Catalog# CRL-1573)

### 2. Plasmids

1. pLJC5-KOZAK-TMEM115-3HA (DU68534 available at MRCPPU depository at MRCPPUreagents@dundee.ac.uk). This is the GolgiTag expression construct
2. pLJC5-KOZAK-3HA-Empty (DU70022 available at MRCPPU depository at MRCPPUreagents@dundee.ac.uk). This is the ControlTag expression construct
3. pVSVG. Lentivirus envelope plasmid. Lenti-X HTX Packaging system (Clontech. Catalog# 631247).
4. pGag/Pol. Lentivirus Gag/Pol plasmid. Lenti-X HTX Packaging system (Clontech. Catalog# 631247).

### 3. Media and Reagents

1. Growth Media: Dulbecco's Modified Eagle's Medium (DMEM) (GIBCO 11960-085); 10% Foetal Bovine Serum (FBS) (Sigma F7524 Batch# BCBW6817); 1% L-Glutamine (GIBCO 25030024); 1% PenicillinStreptomycin (GIBCO 15140122).
2. Selection Media: Growth Media with 2ug/ml Puromycin (Sigma P9620).
3. Transfection media: OptiMem (GIBCO 31985-062)
4. Dulbecco's phosphate-buffered saline (PBS) (GIBCO 14190169)
5. KPBS Buffer: 136mM KCL, 10 mM KH<sub>2</sub>PO<sub>4</sub>. Adjust to pH 7.25 with KOH. (Note On the day of use, add Roche cOmplete protease inhibitor cocktail tablet (REF# 11873580001) and Roche PhosSTOP tablet (REF# 04906837001))
6. Lysis Buffer: 50mM HEPES-KOH pH 7.4; 40mM NaCl; 2mM EDTA; 1.5mM NaVO<sub>4</sub>; 50mM NaF; 10mM NaPyrophosphate; 10mM NaBetaGlycerophosphate; 1% TritonX100. (Note On the day of use, add Roche cOmplete protease inhibitor cocktail tablet (REF# 11873580001))
7. Pierce™ BCA Protein Assay Kit (Ref# 23227. Lot# VA294738)
8. Linear polyethylenimine (PEI Max 40K. Polyscience #24765)
9. Polybrene Infection/Transfection reagent (Millipore TR1003G)
10. NuPAGE 4xLDS sample buffer (Invitrogen. Lot #1941674)
11. 2-Mercaptoethanol (Sigma-Aldrich M6250)

### 4. Equipment

1. Homogenizer, 2-ml vessel (VWR, Catlog# 89026-386) and plain plunger (VWR, Catlog#. 89026-398)
2. Orbiter (The Belly Dancer Shaker; IBI Scientific, model # BDRAA115S)
3. DynaMag™ Spin Magnet. (Thermofisher scientific. Catalog #12320D)
4. Incubator with FPI-sensor system and display controller MB1 (BINDER GmbH. Model: CB150. Power Output: 1.40kW, 230V, 6.1 Amp). This incubator has CO<sub>2</sub> and O<sub>2</sub> control.
5. Microcentrifuge with thermostat (VWR Micro Star 17R. S/N 42209232. REF# 521-1647)

### 5. Consumables

1. Pierce™ Anti-HA magnetic beads (Thermofisher. Catalog # 88837)
2. 10cm and 15cm tissue culture Petri Dishes(ThermoFisher. Catalog# 172931 and 168381 respectively).
3. 1.5ml low binding Eppendorf tubes (Sarstedt. REF # 72.706.600).
4. 15ml CELLSTAR® tubes (Greiner bio-one. Catalog# 188271).



5. 50ml CELLSTAR® tubes (Greiner bio-one. Catalog# 227261).
6. Standard 1ml and 200µl Pipette tips (Greiner bio-one. Catalog# 686271 and 685261 respectively).
7. Syringe filter (0.45µm. Sartorius, Item # ST16537-Q)
8. Syringes (10ml) (Medicina. REF# IVS10. LOT# 19111004).
9. Disposable Cell Lifter (FisherBrand. Catalog # 08-100-240)

## Safety warnings

⚠ Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.

### Safety information

Handle lentivirus under sterile condition in a CATEGORY 2 biological safety cabinet.



## Packaging GolgiTag and ControlTag construct into lentivirus system

2d 0h 35m

1

### Safety information

This is done under sterile condition in a CATEGORY 2 biological safety cabinet

2 Grow 2 dishes of HEK293FT cells to 50-60% confluency in Growth media in 10cm Petri Dish.

### Note

One dish is needed to generate the GolgiTag lentivirus and the other is needed for the ControlTag lentivirus.

3 Prepare a transfection mix to generate a lentivirus expressing the **GolgiTag** in 1.5ml Eppendorf tube containing:

- 3.8 µg pGag/Pol plasmid
- 2.2 µg pVSVG plasmid
- 6 µg pLJC5-KOZAK-TMEM115-3HA plasmid
- 300 µL OptiMem

### Note

We purify plasmids using a QIAGEN HiSpeed® Plasmid Maxi kit [Lot# 166034460] following manufactures protocols and ensure sterile reagents are used and mixtures prepared in tissue culture hood to avoid contamination.

4 Prepare a transfection mix to generate a lentivirus vector expressing **ControlTag** in 1.5ml Eppendorf tube containing:

- 3.8 µg pGag/Pol plasmid
- 2.2 µg pVSVG plasmid
- 6 µg pLJC5-KOZAK-3HA-Empty plasmid
- 300 µL OptiMem



## Note

We purify plasmids using a QIAGEN HiSpeed® Plasmid Maxi kit [Lot# 166034460] following manufactures protocols and ensure sterile reagents are used and mixtures prepared in tissue culture hood to avoid contamination.

- 5 Prepare PEI mixture in 1.5ml Eppendorf tube.
  - 40  $\mu$ L 1mg/ml PEI Max 40K dissolved in distilled water
  - 600  $\mu$ L OptiMem
- 6 Incubate each mixture (Step 3, 4, 5) separately for ~ 00:05:00 at Room temperature .
- 7 Add 310  $\mu$ L PEI Mixture (Step 5) to the **GolgiTag** in 1.5ml Eppendorf tube (Step 3) and **ControlTag** in 1.5ml Eppendorf tube (Step 4) mixtures.
- 8 Mix and incubate each mixture at Room temperature for 00:30:00 .
- 9 Add each mixture dropwise using a P1000 sterile pipette into a single HEK293FT containing plate.
- 10 Incubate cells 37 °C for 24:00:00 .
- 11 Replace media with fresh Growth Media (discard old media) and incubate cells at 37 °C for further 24:00:00 .
- 12 Collect media that contains the lentivirus and pass through 0.45 $\mu$ m syringe filter. This is now the lentivirus infection media. This could be used immediately or snap frozen in liquid nitrogen and stored at -80 °C .

## Infecting (Transducing) and Selecting cells stably expressing GolgiTag

13

**Safety information**

This is done under sterile condition in a CATEGORY 2 biological safety cabinet

14 Mix  5 mL infection media produced in step 12 with  5 mL fresh Growth Media (in a 15ml Eppendorf tube).



15 Add Polybrene reagent to Mix produced in step 14 using a stock of 10mg/ml dissolved in MilliQ water (which had been sterile filtered) to a final concentration of  10 Mass Percent Polybrene .

16 Gently add media to 10cm plate of HEK293 cells already at ~60% confluency.

**Note**

Other types of mouse and human cells could be infected with this media too

17 Incubate at  37 °C for  24:00:00 .

1d



18 Change media to Growth Media and incubate  37 °C for another  24:00:00 .

1d



19 To select for cells stably expressing GolgiTag or ControlTag, replace media with  10 mL freshly prepared Selection Media .

20 After  24:00:00 , there will be observable death and a considerable number of dead floating cells will be observed.

1d

21 Change Selection Media each 24h for 3-5 days (depending on the efficiency of the transduction process). At this stage colonies of cells stably expressing the GolgiTag or empty vector should have reached 100% confluency.

22 Split cells into more dishes maintaining cells in Selection Media.

**Note**

Cells should be grown only in Selection Media, going forward. Cells can be frozen down and stored long term in liquid nitrogen. Make note in cell freezer catalogue to culture cells in Growth Media when thawed until cells are fully recovered, then grow in Selection Media.

## Cell preparation for GolgiTag and ControlTag immunoprecipitation (To be undertaken 1-2 days prior to immunoprecipitation)

- 23 Split cells stably expressing GolgiTag or ControlTag in 15cm Petri dish. Allow to grow near confluency (approximately 24-48h depending on the seeding density).

## Pre-clearing of anti-HA beads (To be done on the day of immunoprecipitation before cell lysis)

2m 30s

24

**Note**

Steps should be done separately for both GolgiTag and ControlTag.

- 25 Pipette  200  $\mu$ L anti-HA bead slurry into 1.5ml Eppendorf tube. 

- 26 Insert tube containing bead anti-HA bead slurry onto a DynaMag™ Spin Magnet for  00:00:30 and carefully remove supernatant using a P1000 pipette. 

- 27 Remove bead from DynaMag™ Spin Magnet, add  1 mL ice cold KPBS , gently pipette up and down to disperse any clumps, put on DynaMag™ Spin Magnet for  00:00:30 and carefully remove supernatant using a P1000 pipette. 

- 28 Repeat wash step described in step 27 a further 3 times: 

- 28.1 (Wash 1/3): Remove bead from DynaMag™ Spin Magnet, add  1 mL ice cold KPBS , gently pipette up and down to disperse any clumps, put on DynaMag™ Spin Magnet for 



00:00:30 and carefully remove supernatant using a P1000 pipette.

28.2 (Wash 2/3): Remove bead from DynaMag™ Spin Magnet, add 1 mL ice cold KPBS , gently pipette up and down to disperse any clumps, put on DynaMag™ Spin Magnet for 00:00:30 and carefully remove supernatant using a P1000 pipette.

30s

28.3 (Wash 3/3): Remove bead from DynaMag™ Spin Magnet, add 1 mL ice cold KPBS , gently pipette up and down to disperse any clumps, put on DynaMag™ Spin Magnet for 00:00:30 and carefully remove supernatant using a P1000 pipette.

30s

29 After the last wash step, incubate the bead slurry in 50 µL KPBS and keep beads On ice .



## Cell lysis for immunoprecipitation

30

### Note

Steps should be done separately for both GolgiTag and ControlTag.

31 Place cells On ice slighted tilted to enable medium to be easily aspirated.

32 Remove media and immediately wash plates with 10 mL cold (4°C) PBS .



33 Immediately, gently aspirate PBS from cells.

34 Gently wash one more time repeating steps 32 and 33:



34.1 Remove media and immediately wash plates with 10 mL cold (4°C) PBS .

34.2 Immediately, gently aspirate PBS from cells.



- 35 Carefully remove PBS and add  1 mL cold (4°C) KPBS . 
- 36 Whilst cells are maintained  On ice , scrape cells off dish with a disposable cell lifter.
- 37 Carefully transfer scrapped cells using an P1000 pipette into a 1.5ml Eppendorf tube.
- 38 Pellet cell by centrifugation at  1000 x g, 4°C, 00:02:00 . 
- 39 Carefully discard supernatant using an P1000 Pipette.
- 40 Resuspend pellet in  950 µL ice cold KPBS .
- 41 Pipette cells up and down with a P1000 pipette ensure cells are in suspension and remove  50 µL into a 1.5 ml Eppendorf tube marked as "**cell lysate INPUT**". Keep this  On ice
- Place the tube containing the remaining ~900µl of resuspended cells on ice and take sample in the cold room.
- Note

Remaining steps of this section and all of the next section to be undertaken in cold room.
- 42 Pipette the  900 µL cells into a clean 2ml Dounce homogeniser. 
- 43 Homogenise cells with exactly 25 strokes, making sure bubbles are avoided.
- To verify efficacy of homogenisation, take 10µl on a slide, briefly check under the microscope. There should be about at least 80% cell lysis, indicated by the visible presence of nuclei without surrounding cytoplasm.

**Note**

If cells are not sufficiently lysed undertake another 10 strokes of homogenisation and check for cell lysis a second time-repeat as necessary but avoid over homogenisation, to maintain integrity of organelles.

- 44 Centrifuge lysate at  1000 x g, 4°C, 00:02:00 . This is to remove any non-lysed cells, nuclei and debris such as cell membranes. 
- 45 Remove supernatant (which should contain organelles) into a new 1.5 ml Eppendorf tube. The lysate at this stage will be between 700-800µl.
- 46 Remove  50 µL from the homogenate into a separate 1.5ml Eppendorf labelled as **IP INPUT**. Keep  On ice .

**GolgiTag and ControlTag immunoprecipitation**

17m

47

**Note**

**This must be done in a cold room.**

- 48 Take out Eppendorf tube containing the washed anti-HA beads from step 29.
- 49 Add the cell homogenate from step 45 into the eppendorf tube containing the HA-beads.
- 50 Ensure that bead clumps are dispersed by gently pipetting up and down 5 times using a P1000 Pipette.  
- 51 Incubate the cell homogenate-HA-bead slurry for  00:05:00 on orbiter.  5m
- 52 Separate bead by putting tube on DynaMag™ Spin Magnet for  00:00:30 .  30s



53 Collect  50  $\mu$ L of supernatant into separate 1.5ml Eppendorf tube label as **FT Supernatant (FT = flow-through)**. Keep  On ice .

54 Wash bead with 1000  $\mu$ L of ice cold KPBS 3 times:

54.1 (Wash 1/3): Wash bead with  1000  $\mu$ L ice cold KPBS .

54.2 (Wash 2/3): Wash bead with  1000  $\mu$ L ice cold KPBS .

54.3 (Wash 3/3): Wash bead with  1000  $\mu$ L ice cold KPBS .

55 Incubate beads with  500  $\mu$ L ice cold KPBS and prepare a slurry by pipetting up and down 5 times with a P1000 pipette and transfer slurry into a clean 1.5 ml Eppendorf tube.

#### Note

This step is important as it ensure contamination from cell extract sticking to tube plastic is minimized.

56 Insert tube containing anti-HA bead slurry onto a DynaMag™ Spin Magnet for  00:00:30 and carefully remove supernatant using a P1000 pipette.

30s

57 Add  80  $\mu$ L ice-cold lysis buffer to the HA beads containing the immunoprecipitated Golgi.

#### Note

Depending on application different solubilization buffers could be deployed at this stage.

58 Incubate for  00:10:00  On ice .

10m



59 Place tube on DynaMag™ Spin Magnet for 00:00:30 to separate beads.

30s

60 Collect supernatant into new Eppendorf tube.

#### Note

This is important to avoid contamination.

61 To ensure total removal of beads, add a further 20 µL ice-cold lysis buffer to HA-Beads and place the beads again on DynaMag™ Spin Magnet for 00:00:30 and collect supernatant into Eppendorf tube containing the solubilized Golgi in step 60.

30s

62 Determine the protein concentration of the IP using a BCA Protein Assay Kit. Use 5 µL solubilised Golgi (do not dilute) for quantification purpose. Quantification should be done in duplicate.



63 Depending on future application samples can be snap frozen in liquid nitrogen and stored at -80 °C or diluted into 4xLDS loading buffer supplemented with fresh 5% (by vol) 2-mercaptoethanol prior to analysis on SDS-polyacrylamide gel electrophoresis and immunoblot analysis. For Immunoblot purpose, 2-4 µg of solubilised Golgi is enough to detect Golgi markers (GM130, ACBD3 and GOLGIN97).



## Preparation of INPUT sample

20m

64 Pellet cells by centrifugation at 1000 x g, 4°C, 00:02:00 .



65 Discard supernatant and resuspend cell pellet in 80 µL ice-cold lysis buffer .

66 Pipette up and down 10-20 times with a P200 pipette.



67 Incubate lysate On ice for 00:20:00 .

20m



68 Centrifuge lysate at  13000 x g, 4°C, 00:10:00 .



69 Collect clarified cell lysate supernatant in fresh tube and determine protein concentration using a BCA Protein Assay Kit. Dilute an aliquot of the supernatant in lysis buffer at a ratio 1:10, before protein quantification.



70 Depending on future application the clarified cell lysate can be snap frozen in liquid nitrogen and stored at  -80 °C or diluted 4X into 4xLDS loading buffer supplemented with fresh 5% (by vol) 2-mercaptoethanol prior to analysis on SDS-polyacrylamide gel electrophoresis and immunoblot analysis.

## Preparation of Flow-through and IP INPUT sample

20m

71 Add  100 µL ice-cold lysis buffer to tubes marked **FT** and **IP INPUT**.

72 Mix thoroughly by pipetting up and down using P1000.



73 Incubate  On ice for  00:20:00 .

20m



74 Centrifuge lysate at  13000 x g, 4°C, 00:10:00 .



75 Collect clarified cell lysate supernatant in fresh tube and determine protein concentration using the BCA method. Dilute an aliquot of the supernatant in lysis buffer at a ratio 1:10, before protein quantification.



76 Depending on future application the clarified cell lysate can be snap frozen in liquid nitrogen and stored at  -80 °C or diluted 4X into 4xLDS loading buffer supplemented with fresh 5% (by vol) 2-mercaptoethanol prior to analysis on SDS-polyacrylamide gel electrophoresis and immunoblot analysis.

## ASAP Protocol Particulars: context and tips

77



[https://www.youtube.com/embed/-9sLvpuW3gg?si=VciFWaOQwS\\_hqH0D](https://www.youtube.com/embed/-9sLvpuW3gg?si=VciFWaOQwS_hqH0D)