

May 17, 2025

iPSC-derived organoid generation

 In 2 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.n92ldnpx7v5b/v1>

Irene Fernández-Carasa^{1,2}, Veronica Testa^{1,2}, Antonella Consiglio^{1,2}, Miquel Vila³

¹Department of Pathology and Experimental Therapeutics, Bellvitge University Hospital- IDIBELL, 08908 Hospitalet de Llobregat, Spain;

²Institute of Biomedicine of the University of Barcelona (IBUB), Carrer Baldri Reixac 15-21, Barcelona 08028, Spain.;

³VHIR-CIBERNED-ASAP

Vilalab Public



Miquel Vila

VHIR-CIBERNED-ASAP

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldnpx7v5b/v1>



Protocol Citation: Irene Fernández-Carasa, Veronica Testa, Antonella Consiglio, Miquel Vila 2025. iPSC-derived organoid generation. protocols.io <https://dx.doi.org/10.17504/protocols.io.n92ldnyp7v5b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 16, 2025

Last Modified: May 17, 2025

Protocol Integer ID: 218435

Keywords: derived organoid generation, organoid generation, ipsc

Funders Acknowledgements:

ASAP

Grant ID: ASAP-020505

Abstract

iPSC-derived organoid generation

- 1 iPSCs are grown until 90% of confluence and 10 μ M of RI is added to culture plate 4/6 hours before starting the protocol. Then, cells are dissociated in 1:1 accutase/DPBS, detached and centrifuged at 800xrpm for 5 minutes.
- 2 Cell pellet is resuspended in 1 mL of fresh mTeSRTM1, passed through a 40 μ m cell strainer (Clear Line[®]) and counted. Dissociated iPSCs are plated in ultra-low attachment 6-well plastic plate (Corning[®]) at a density of 4×10^6 cells per well, with mTeSRTM1 medium supplemented with 10 μ M of RI.
- 3 From then on, cells are kept under rotation on CO₂ resistant orbital shaker (Leica) at 120xrpm and maintained under normal incubation conditions.
- 4 After 2 days, medium is replaced with SRM supplemented with LDN and SB, a dual SMAD inhibition to promote neuronal commitment.
- 5 On day 4, and every 3 days, half of SRM medium is replaced with fresh factors and the addition of SAG, Purmorphamine, CHIR and AA.
- 6 On day 12, SRM medium is replaced by NBN2B27-vitA, which consisted of NeurobasalTM supplemented with 1x N2, 2x B27TM-vitA, 1% of P/S and 1% of Glut. N2B27-vitA is changed every other day by the addition of 100 nM of LDN, 1 μ M of SAG, 2 μ M of Purmorphamine, CHIR and AA, 50 ng/mL of FGF-2, and 50 ng/mL of Endothelial Growth Factor (EGF; R&D Systems), to promote conversion towards a DAn fate.
- 7 From day 22 onwards, spheres are matured to vmO in B27-vitA supplemented with BDNF, GDNF, TGF- β 3, cAMP and AA. Half of the medium is changed every 2-3 days and vmO are cultured until day 45, before being used for experiments.